Advincula Arch™ Handle

Reusable Uterine Manipulator Handle

UMH750

Directions for Use (English)

**NON-STERILE
STERILIZE BEFORE USE**

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.


CE0086

DEVICE DESCRIPTION

The Advincula Arch™ Uterine Manipulator Handle is a reusable device used with a sterile, single-use RUMI® Tip.

The Advincula Arch Uterine Manipulator Handle (Handle) along with the RUMI Uterine Manipulator Tip (Tip) is used for positioning the uterus during laparoscopic surgery and in dye delivery where chromopertubation is necessary.

The Grip, Shaft and Snap Drum are constructed of medical grade autoclavable materials. The Snap Drum serves for attachment of the Tip.

WARNINGS

- The uterus must be sounded for both depth and direction prior to application of the device.
- **DO NOT** use the Handle as a uterine sound.
- **NEVER** attempt uterine manipulation without a clear view of the uterus.
- As with all uterine manipulating devices, a careful clinical evaluation should be performed prior to use.
- Certain clinical conditions may present a uterus, which is more prone to perforation and bleeding.
- Dye injection should be performed SLOWLY. Because of the efficient air/liquid seal created at the internal cervical os by the balloon, rapid injection of fluids (dye) may create intrauterine pressures, which could cause uterine damage and/or result in fallopian tube spasms.
- The Handle has only been tested for use with patients who have been anesthetized. Due to the need to dilate the cervix, the Handle is not recommended for use in non-anesthetized patients.

CAUTION

U.S. Federal law restricts this device to sale by or on the order of a physician.

INTENDED USE / INDICATIONS FOR USE

The Handle is designed for use in operative endoscopy (laparoscopy), where a uterus is present and where positioning of the uterus, fallopian tubes and ovaries or vagina is desirable. These types of surgeries include laparoscopic tubal ligation, diagnostic laparoscopy and/or operative laparoscopy. The Tips also provide for dye delivery in those procedures requiring chromopertubation.

CONTRAINDICATIONS

The Handle should not be used in patients who are pregnant or who are suspected of being pregnant, planned gamete intrafallopian transfer procedures, in patients who have an IUD in place, in patients with suspected pelvic infections and in cases where the surgeon deems it inadvisable or finds it difficult to insert the silicone tip into the cervix or uterus.

PRECAUTIONS

- Sterilize before use.
- Dilate the cervix to Hegar/Hank 8 (French 24) to ease patient insertion.
- Refer to the RUMI® Uterine Manipulator Tip Directions for Use for further precautions.
- Inspect Handle prior to use for proper operation.

ADVERSE REACTIONS

The following adverse reactions have been suspected or reported to be associated with all uterine manipulators. The order of listing does not indicate frequency or severity: cramping, infection, uterine and fallopian tube spasm with associative temporary physiological blockage of patient fallopian tubes and uterine perforation.

DIRECTIONS FOR USE

Note: Refer to the RUMI Uterine Manipulator Tips Directions for Use to ensure proper selection of Tip size and Tip attachment / detachment instructions.

1. Select a Tip which is less than or equal to the sounded depth of the uterus and attach it to the Handle.
2. Ensure the Tip Catheters are secure in the Catheter Channels of the Handle (see Figure 1).



Figure 1

RECOMMENDED CLEANING INSTRUCTIONS

Before cleaning the Handle, remove the Tip as per the RUMI® Uterine Manipulator Tips Directions for Use.

1. Prepare the enzymatic detergent according to the manufacturer's recommendations and submerge the Handle in the detergent. Soak for ten (10) minutes.
2. Scrub the submerged Handle with a soft sponge and agitate. Clean lumens and crevices, using a pipe cleaner/brush or flush with a syringe when appropriate.
3. Rinse in warm 100 °F - 120 °F (38 °C - 49 °C) tap water for one (1) minute. Thoroughly flush all lumens and other difficult to reach areas.
4. Ultrasonically clean the Handle for ten (10) minutes in a neutral pH detergent (Neutrad® or acceptable alternative). Prepare the detergent according to the manufacturer's recommendations.
5. Rinse the Handle with clean tap water for one (1) minute; repeat rinse twice.
6. Dry the Handle thoroughly with a clean, lint free cloth.
7. Repeat this cleaning procedure if the Handle appears to be soiled after initial cleaning and until there is no sign of soil on the device.

RECOMMENDED STEAM AUTOCLAVE STERILIZATION PARAMETERS

- The instrument must be thoroughly cleaned of all foreign matter prior to sterilization.
- Follow the manufacturer's instructions for operation and loading of steam autoclave.
- There must be direct steam exposure to all surfaces of the instrument being sterilized.

Pre-vacuum Cycle:

Preconditioning Pulses: 4

Exposure Temperature: 270 °F - 273 °F (132 °C - 134 °C)

Exposure time: 4 minutes

Gravity Displacement Cycles:

Exposure Temperature: 270 °F (132 °C)

Exposure Time: 15 minutes

Exposure Temperature: 250 °F (121 °C)

Exposure Time: 30 minutes

WARRANTY

The Advincula Arch™ Uterine Manipulator Handle is warranted by CooperSurgical for a period of 120 days from shipment. This warranty applies only to the original purchaser and only against defects in workmanship or materials, which under normal use render the instrument inoperable. The original purchaser shall prepay the shipping costs for returning a CooperSurgical instrument for warranty service. CooperSurgical will, at its option and without charge, either repair or replace any instrument, which CooperSurgical determines to be defective in material or workmanship. Instrument damage caused by misuse or accident, shall void this warranty.

CooperSurgical disclaims any liability for special, incidental, consequential, punitive or exemplary damages arising out of the use of any CooperSurgical instrument. CooperSurgical's liability in all events is limited to, and shall not exceed, the purchase price paid. Except as expressly provided in this limited warranty section, CooperSurgical makes no representations or warranties, expressed or implied, as to the instruments, including warranties of merchantability and fitness for a particular purpose, all of which are hereby expressly disclaimed and excluded. No warranty or affirmation of fact, express or implied, other than as set forth in this limited warranty, is made or authorized by CooperSurgical.

EXPLANATION OF SYMBOLS

 Reorder Number

 Serial Number

 Latex Free

 **ATTENTION:**
See instructions for use.

 **CE0086** Product conforms to the Medical Device Directive 93/42/EEC

 Authorized Representative in the European Community.

Neutrad® is a registered trademark of Decon Labs, Inc.

CooperSurgical is a registered trademark of CooperSurgical, Inc.

All other marks bearing the symbols "TM" and "®" are trademarks and registered trademarks of CooperSurgical, Inc.

© 2011 CooperSurgical, Inc. All rights reserved.

Made in the USA

95 Corporate Drive • Trumbull, CT 06611 USA
Phone: (800) 243-2974
Fax: (800) 262-0105
www.coopersurgical.com

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Neutrad® é uma marca registada da Decon Labs, Inc.
CooperSurgical é uma marca registada da CooperSurgical, Inc.
© 2011 CooperSurgical, Inc. Todos os direitos reservados.
Leisegang Feinmechanik GmbH