
Integrated
CerviPATH™ Management System

ASSEMBLY INSTRUCTIONS



Distributed by:
CooperSurgical LEISEGANG
Part #35834 • Revised 9/04

CooperSurgical LEISEGANG
95 Corporate Drive, Trumbull, CT 06611
800-243-2974

CerviPath:

1. In the main CerviPath window, select “Options” > “User Preferences” from the menu bar.
2. In the “Main Settings” tab, change the Login name and input the new password and remaining information.
3. Click the “Save” button.
4. Click the “Close” button.
5. Shut down and re-start the CerviPath program.
6. In the “Login” window, click on the “Login” drop-down menu and choose the new “Login” name. Enter the new password.

STEP 20

CD Formatting Instructions:

When a CD backup is required, the user must first format the CD-RW. Before starting the CerviPath software, place a blank CD-RW disc in the DVD/CD tray. A window will automatically pop up. Click on “Autorun empty CD using InCD” (Fig 1).




Fig 1

In the Format screen click on “Start” (Fig 2).




Fig 2

The formatting could take up to several minutes; the status bar displays the percentage completed (Fig 3).




Fig 3

When finished, a window will pop up stating “Format complete” (Fig 4). The DVD/CD tray will open for you to remove the CD-RW disc, until needed.




Fig 4

STEP 1


Assemble the Colposcope by inserting the scope stand [BOM (Balance-O-Matic) or Swing] into the five-wheel, weighted roller base. Attach the weighted base to the stand by inserting and tightening the hex bolt and washer through the hole located at the center of the weighted base from the underside.



STEP 1

STEP 2

Insert the Colposcope head into the top of the stand and secure into place by tightening the gross height adjustment knob. Attach the Scope head power supply plug to the transformer. Attach Colposcope power cord to the transformer.



STEP 2

STEP 3

Remove Pole Mount Bracket box from tray assembly box and disassemble Bracket.



STEP 3

STEP 4

Remove Black Plastic Split bushing from Laptop System package, peel off backing of double-sided tape and place onto each side of Pole Bracket.



STEP 4

STEP 5

Attach Pole Bracket onto outside pole of BOM, 10 inches from top of pole and at 90° to top (see picture for orientation). Tighten screws using Allen keys provided. (On Swing Scope attach bracket just below top of pole.)



STEP 5

STEP 6

Place both plastic washers from Allen key set onto bracket post (6a) and attach Tray assembly onto pole bracket (6b). Tighten set screws with small Allen key from set.



STEP 6a



STEP 6b

STEP 7

Loop Notebook Cable Lock through bracket under Tray and pull tight.

STEP 8

Place NoteBook on Tray, attach cable lock to Notebook by placing teeth of lock into security slot of Notebook and push in button on back of lock to secure (use key to unlock).

STEP 9

Screw Video Tubus into camera, insert camera into head with **connections facing downward** and tighten thumb screw to secure. (Make sure camera is inserted all the way into head and is level.)

STEP 10

Attach power cord to NoteBook.

STEP 11

Attach footpedal to the lower USB port on the back of the NoteBook.

STEP 12

Attach camera cable to the **Upper USB Port** on the back of the NoteBook.

STEP 13

Feed camera cable thru arm as pictured and connect to camera.

STEP 14

Remove arm slot cover, feed camera cable, power cord and Foot Pedal cable from Notebook into slot and replace cover.



STEP 7



STEP 8



STEP 9



STEP 10



STEP 11



STEP 14



STEP 12



STEP 13

STEP 15

Using Expandable Cable Wrap, wrap Camera cable from camera and power cord from transformer together.

STEP 16

Using Expandable Cable Wrap, wrap Notebook power cord, Foot Pedal cable and Transformer power cord together.

STEP 17

CerviPath Software and printer drivers have been pre-loaded. Connect your CerviPath Hardware Key to the USB port on the side of the Notebook.

STEP 18

Remove Printer from box and set up per enclosed HP instruction chart. Turn on Notebook/Printer. Align the IR port of the Notebook with the IR on the printer (within 6 FEET and no more than a 30 degree angle from the IR port). A connection message will appear on the lower right hand side of the Notebook screen "Another Computer is Nearby" indicating that an IR connection has been established. A connection symbol will be displayed on the lower right side of the Notebook screen when the two units are connecting.

STEP 19

Changing the User Names and Passwords:

Windows XP:

1. Click Start > Control Panel.
2. Click "User Accounts".
3. In the "Or Pick an Account to Change" section, select the account name to be changed. (This will usually be the "Admin" account)
4. Click "Change My Name". Enter the new User Name and click the "Change Name" button.
5. The password for this account can be created and modified by selecting "Create a Password". Follow the directions and click the "Create a Password" button.
6. Click "Yes", "Make Private" for added file security.
7. Close all windows when finished.



STEP 15



STEP 16



STEP 18

NOTE: “Directions for Use” are included in the CerviPATH™ Management Software. Your Desktop will display 3 icons related to your CerviPATH™ Management System.

1. To use the system, click on the “CerviPATH™” icon.
2. For directions for Clinician’s use, click on the “Clinician Training Guide” icon.
3. For directions for Administrator’s use, click on the “Administrator Training Guide” icon.