# 10057 • Mystic® II M-Style® Mushroom® Cup

## Vacuum-Assist Delivery System


### Directions for Use (English)

CE 0086

**mityvac®**  
VACUUM-ASSIST DELIVERY SYSTEM

#### DEVICE DESCRIPTION

The Mityvac® Mystic® II M-Style Mushroom® Cup is designed to assist a practitioner in the delivery of an infant during childbirth. The Mityvac Mystic II is an integrated pump and cup system, packaged, and sterilized for immediate use.

#### WARNINGS

- **DO NOT exceed recommended vacuum levels.**
- **Never apply the vacuum cup to any portion of the infant's face.**
- **Vacuum-assisted delivery should only be performed or supervised by a trained and experienced operator.**
- **"There should be a willingness to abandon attempts at vacuum extraction if satisfactory progress is not made."** (ACOG Technical Bulletin #196, Aug. 1994)
- **Abandon vacuum-assisted delivery if:**
  - a. Vacuum cup becomes disengaged (pops off) three times
  - b. Vertex has not advanced substantially with each traction attempt
  - c. There is evidence of fetal scalp trauma
  - d. Cumulative traction time exceeds 10 minutes or total vacuum procedure time exceeds 15 to 30 minutes
- **If the extractor cup becomes disengaged, check for trauma to the fetal scalp before reapplying.**

**CAUTION:** U.S. Federal law restricts this device to sale by or on the order of a physician.

#### INDICATIONS FOR VACUUM-ASSISTED VAGINAL DELIVERY

**No indication is absolute.**

##### Fetal:

- Non-reassuring fetal status
- Failure to deliver spontaneously following an appropriately managed second stage of labor

##### Maternal:

- Need to avoid voluntary expulsive efforts
- Inadequate expulsive effort

#### CONTRAINDICATIONS

Do not initiate vacuum if any of the following conditions exist:

- Non-vertex positions (breech or transverse lie/position) or face or brow presentation
- Suspected cephalopelvic disproportion
- Previous scalp sampling
- Suspected macrosomia, or risk of shoulder dystocia
- Failed vacuum or forceps attempt
- Less than 34 weeks gestation
- Unengaged vertex
- Incompletely dilated cervix
- Need for active device rotation
- Suspected fetal bleeding abnormalities

#### ADVERSE EVENTS

**Fetal Injuries:** Head trauma, bruises, contusions, lacerations, scalp edema, skull fracture, cephalhematoma, subgaleal hematoma, subdural hemorrhage, parenchymal hemorrhage, intracranial hemorrhage, retinal hemorrhage.

**Maternal Injuries:** Vaginal, cervical, uterine, bladder, rectal tissue tears.

#### PREREQUISITES FOR VACUUM-ASSISTED VAGINAL DELIVERY

- Term infant
- Ruptured membranes
- Empty bladder and nondistended rectum
- Engaged vertex position
- Complete cervical dilation and effacement
- Willingness to abandon procedure and plan for alternate delivery route

**These instructions are intended as general guidelines.**

**Practitioners should refer to current institutional and recognized guidelines that address vacuum-assisted delivery.**

#### VACUUM-ASSISTED VAGINAL DELIVERY GUIDELINES

1. Check the integrity of the Mityvac Mystic II Obstetrical Vacuum-Assist Delivery System by pressing the cup to the palm of your gloved hand and applying vacuum. The gauge needle should remain steady. (NOTE: The gauge needle should rest in the "ZERO" box when vacuum is not applied.)


**CAUTION: DO NOT USE DEVICE IF GAUGE NEEDLE DOES NOT INDICATE ZERO IN THE ABSENCE OF A VACUUM.**

2. Carefully examine fetal presentation and position prior to positioning the cup over the flexion point.

**CAUTION: DO NOT USE PRODUCT IF FLEXION POINT IS NOT ACCESSIBLE.**


3. Wipe the scalp as clean as possible.
4. To insert the M-Style Mushroom Cup, fold the cup against the stem at a 90° angle (see **Figure 1**).

5. With the stem in the bent position, grasp the mushroom shaped portion of the cup with the fingers. The stem should be parallel to palm and wrist. Separate the labia with the other hand. Continuing with the stem in the bent position, press the cup portion downward and inward over the fourchette to make contact with the fetal scalp. Apply the center of the cup over the flexion point (see **Figure 2**).

**Figure 1 - M-Style Mushroom Cup**

The flexion point can be located by identifying the posterior fontanelle and then moving the finger anteriorly along the sagittal suture approximately **1.2 in. (3 cm)**. The sagittal suture should pass under the middle of the cup.

**Figure 2 - Infant head**

**CAUTION: NEVER APPLY THE CUP TO ANY PORTION OF THE INFANT'S FACE.**

6. Gently move a finger around the rim of the cup to remove intrusion of extraneous tissues and to ensure proper cup placement.

- With the cup placed over the flexion point, raise the vacuum level to approximately **3.9 in. Hg (10 cm Hg)** (yellow shaded area on gauge face) to initiate vacuum. Recheck the rim of the cup for interposed extraneous tissue.
- With the onset of contraction, rapidly raise the vacuum to **15 - 22.8 in. Hg (38 - 58 cm Hg)** (green shaded area on gauge face) and begin applying traction.<sup>1,2</sup> Always apply traction in harmony with contractions and along the pelvic axis.

**CAUTION: DO NOT EXCEED RECOMMENDED VACUUM LEVELS.**

- When the contraction is no longer effective, discontinue traction. You may reduce the vacuum to approximately **3.9 in. Hg (10 cm Hg)**

(yellow shaded area on gauge face) while awaiting the next contraction or maintain vacuum at current level. Reduce vacuum by slowly pushing or pulling on the vacuum release lever located directly beneath the gauge (see **Figure 3**).

Figure 3 - Vacuum release lever

- Recheck for interposed extraneous tissue prior to each tractive effort.
- If traction is misaligned or too forceful, the vacuum cup may disengage (pop off). In case of pop off, check fetal scalp for trauma before reapplying vacuum cup.

**CAUTION: ABANDON VACUUM-ASSISTED DELIVERY IF THE VACUUM CUP DISENGAGES (POPS OFF) THREE TIMES.**

- With each successive contraction, draw the head gently over the perineum. Once the head is delivered, release vacuum (see **Figure 3**) and remove the cup. Continue delivery in the usual manner. **In the unlikely event that vacuum does not release, with caution, use scissors and cut the white flexible joint between the stem and cup, being careful to avoid the fetal scalp.**

**CAUTION: DO NOT OPERATE VACUUM PUMP AT TRACTIVE LEVELS FOR MORE THAN TEN CUMULATIVE MINUTES OF TRACTION OR TOTAL PROCEDURE TIME OF 15 TO 30 MINUTES.**

- Carefully examine infant's head and observe vital signs at regular intervals to ensure infant's well-being.
- Dispose of the Mityvac Mystic II delivery system in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.
- Document the use of the Mityvac Mystic II delivery system and notify nursery staff per hospital protocol.

#### PREREQUISITES FOR VACUUM-ASSISTED CESAREAN DELIVERY

- Proper indication for operative delivery
- Proper indication for use of the vacuum-assisted delivery device

#### VACUUM-ASSISTED CESAREAN DELIVERY GUIDELINES

- Check the integrity of the Mityvac Mystic II Obstetrical Vacuum-Assist Delivery System by pressing the cup to the palm of your gloved hand and applying vacuum. The gauge

needle should remain steady. (NOTE: The gauge needle should rest in the "ZERO" box when vacuum is not applied.)

**CAUTION: DO NOT USE DEVICE IF GAUGE NEEDLE DOES NOT INDICATE ZERO IN THE ABSENCE OF A VACUUM.**

- Enter the uterus routinely and assess fetal head position.
- If the fetal head is high and accessible beneath the uterine incision, wipe the scalp as clean as possible and place vacuum cup over the flexion point. If the head is low and inaccessible, place gloved fingers under the head and flex upward to bring the scalp under the uterine incision. Wipe scalp as clean as possible and apply the vacuum cup to the flexion point.

**CAUTION: NEVER APPLY THE CUP TO ANY PORTION OF THE INFANT'S FACE.**

- Raise the vacuum to **15 - 22.8 in. Hg (38 - 58 cm Hg)** (green shaded area on gauge face) and begin applying traction to gently bring the head upward through the incision.<sup>1,2</sup>
- Once the head is delivered, release vacuum by slowly pushing or pulling on the vacuum release lever (see **Figure 3**) and remove the cup. Continue delivery in the usual manner. In the unlikely event that vacuum does not release, with caution, use scissors and cut the white flexible joint between the stem and cup, being careful to avoid the fetal scalp.
- Carefully examine infant's head and observe vital signs at regular intervals to ensure infant's well-being.
- Dispose of the Mityvac Mystic II delivery system in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.
- Document the use of the Mityvac Mystic II and notify nursery staff per hospital protocol.

<sup>1</sup> Bofill, James A., MD, et al., "The Mississippi Operative Vaginal Delivery Trial: Lessons Learned", *Contemporary OB/GYN*, October 1998.

<sup>2</sup> Schifrin, Barry S., MD, et al., "Cutting Your Legal Risks with Vacuum Assisted Delivery", *OBG Management*, June 1994.

#### EXPLANATION OF SYMBOLS

 Reorder number	 Keep dry
 Batch code	 <b>ATTENTION:</b> See instructions for use.
 Use by-date	 Product conforms to Medical Device Directive 93/42/EEC
 Latex free	 Sterilized Using Ethylene Oxide
 Do not reuse	 Authorized Representative in the European Community.
 Manufacturer	

Mityvac®, Mystic® II, M-Style® and Mushroom® are registered trademarks of CooperSurgical, Inc., its subsidiaries and/or affiliates. CooperSurgical is a registered trademark of CooperSurgical, Inc. © 2013 CooperSurgical, Inc.

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CooperSurgical è un marchio depositato di CooperSurgical, Inc.