
Advincula Arch™ Adapter

for use with the Advincula Arch

53469

Directions for Use (English)

**NON-STERILE
STERILIZE BEFORE USE**

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Advincula Arch™ Adapter is a reusable device that is used to attach the Advincula Arch Uterine Manipulator Handle to the Uterine Positioning System™ exclusively.

PRECAUTIONS AND WARNINGS

- Failure to properly follow instructions for the Advincula Arch Adapter (Adapter), including instructions supplied with the Uterine Positioning System (UPS, 371550-03) and the Advincula Arch Uterine Manipulator Handle (Manipulator Handle, UMH750), may lead to injury and result in improper functioning of the system.
- The Adapter is only intended for use with the Advincula Arch Manipulator Handle and the UPS. Use of this device with any other Manipulator is not recommended, and could lead to injury and improper functioning of the device.
- Surgical procedures requiring vaginal instrumentation are not sterile. Conventional operating room procedures for maintaining sterility must be observed when the UPS is in use.
- The Adapter is supplied non-sterile and must be cleaned and sterilized prior to use.
- Verify that the correct Adapter is being used with the Manipulator Handle.
- Refer to the UPS Directions for Use for further instructions.
- ALWAYS position the patient and the OR table prior to attaching the UPS to the Manipulator Handle.
- ALWAYS have the patient under general endotracheal anesthesia when the UPS is attached to the Manipulator Handle.

CONTRAINDICATIONS

Refer to the Manipulator Handle (UMH750) Directions for Use for applicable contraindications of the device.

CAUTION

U.S. Federal law restricts this device to sale by or on the order of a physician.

DIRECTIONS FOR USE


Figure 1. Advincula Arch™ Adapter

Before Attaching the Adapter

The following steps **must** be completed prior to attaching the Adapter to the UPS™ arm:

1. Patient is positioned, prepped and the Manipulator Handle is inserted.
2. UPS arm and patient are properly draped.

Attaching the Adapter to the UPS Arm

1. Align the reinforced punch-out hole in the UPS Sterile Drape with the Mating Hole in the UPS arm.
2. Insert the Adapter's Engagement Pin (see Figure 1) into the UPS arm's Mating Hole until secure. Use the Alignment Pin for proper orientation of the Adapter.
3. Ensure that the Adapter is securely in place on the arm.

Before Attaching the Manipulator Handle

The following steps **must** be completed prior to attaching the Manipulator Handle to the Adapter:

1. Patient is positioned for surgery.
2. Draped UPS arm and Adapter are placed in approximate position for attachment near the proximal end of the Manipulator Handle.

Attaching the Manipulator Handle

CAUTION: DO NOT attach the Manipulator Handle to the UPS arm until after the patient is in final position and sedated.

1. Open the Latch on the Adapter (see Figure 2).

CAUTION: Securely hold the UPS handle before the foot pedal is depressed to avoid inadvertent movement of the Manipulator Handle.

Note: For optimal motion, grasp only the handle of the UPS to move the device.

2. Grasp the handle of the UPS, press the foot pedal, and move the arm to align the Adapter with the Manipulator Handle. Take care to minimize movement of the Manipulator Handle within the patient.
3. Position the Mating Pin on the Adapter into the Mating Hole of the Manipulator Handle.


Figure 2. Attaching the Advincula Arch Handle

4. Close and lock the Latch, and ensure that the Manipulator Handle is securely attached to the Adapter with no free movement.

Removing the Manipulator Handle from the Adapter

CAUTION: When securing the Latch, use care to avoid contacting the patient or other accessories in the surgical field.

To detach the Manipulator Handle, unlock and open the Latch on the Adapter while holding the Manipulator Handle and separate it from the Adapter.

Note: There is no need to depress the foot pedal and move the UPS arm to separate the Manipulator Handle from the Adapter.

Removing the Adapter from the UPS

To detach the Adapter from the UPS Arm after surgery is complete, press the Release Button on the Arm (see Figure 3) and pull the Adapter straight out without rotation.


Figure 3. UPS Arm

RECOMMENDED CLEANING INSTRUCTIONS

1. Prepare the enzymatic detergent according to the manufacturers recommendations and submerge the Adapter in the detergent. Soak for ten (10) minutes.
2. Scrub the submerged Adapter with a soft sponge or brush and agitate. Clean crevices, using a pipe cleaner/brush or flush with a syringe when appropriate. Actuate any moving parts to loosen trapped soil.
3. Rinse in warm (100 °F – 120 °F, 38 °C – 49 °C) tap water for one (1) minute. Thoroughly flush difficult to reach areas. Actuate while rinsing.
4. Place the Adapter into a bath containing warm (100 °F – 120 °F, 38 °C – 49 °C) water. Agitate the Adapter by hand for at least one (1) minute. Repeat this process two (2) additional times.
5. Ultrasonically clean the Adapter for twenty (20) minutes in a neutral pH detergent (Neutrad® or acceptable alternative). Prepare the detergent according to the manufacturer's recommendations.
6. Rinse the Adapter with clean, tap water, actuating moving parts while rinsing, for one (1) minute; repeat rinse twice.
7. Dry the Adapter thoroughly with a clean, lint free cloth.
8. Repeat this cleaning procedure if the Adapter appears to be soiled after initial cleaning and until there is no sign of soil on the device.

RECOMMENDED STEAM AUTOCLAVE STERILIZATION PARAMETERS

Pre-vacuum Cycle

Preconditioning Pulses: 4

Exposure temperature: 270 °F – 273 °F (132 °C – 134 °C)

Exposure time: 4 minutes

Gravity Displacement Cycles

Exposure temperature: 270 °F – 275 °F (132 °C – 135 °C)

Exposure time: 15 minutes

Exposure temperature: 250 °F – 253 °F (121 °C – 123 °C)

Exposure time: 40 minutes

Following steam sterilization, allow all components to cool to room temperature. Sudden changes in temperature may damage the components.

Note: The parameters listed on the previous page are only a recommendation. Sterility of the device is the responsibility of the hospital or site performing the sterilization. (Refer to the sterilizer's manufacturer's manual for detail about the sterilization unit.)

DISPOSAL

Dispose of in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.

EXPLANATION OF SYMBOLS

 Reorder Number

 Batch Code

 Latex Free

 **ATTENTION:** See instructions for use.

Rx ONLY **CAUTION:** U.S. Federal law restricts this device to sale by or on the order of a physician.

 Product conforms to the Medical Device Directive 93/42/EEC

 Authorized Representative in the European Community.

Neutrad[®] is a registered trademark of Decon Labs, Inc.

CooperSurgical is a registered trademark of CooperSurgical, Inc.

All other marks bearing the symbols "TM" and "®" are trademarks and registered trademarks of CooperSurgical, Inc.

© 2011 CooperSurgical, Inc. All rights reserved.

Manufactured for:

CooperSurgical

Made in the USA

95 Corporate Drive
Trumbull, CT 06611 USA
Phone: (800) 243-2974
Fax: (800) 262-0105
www.coopersurgical.com

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Neutrad[®] é uma marca registada da Decon Labs, Inc.
CooperSurgical é uma marca registada da CooperSurgical, Inc.
© 2011 CooperSurgical, Inc. Todos os direitos reservados.
Leisegang Feinmechanik GmbH