Endosee® Hysteroscope
with Disposable Diagnostic
(Dx) and (Cx) Cannula for
Cystoscopy Only

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About This Manual

This manual provides instructions on how to use the Endosee® Hysteroscope with the Disposable Diagnostic (Dx) Cannula safely and effectively for its intended use. It is important to read and observe the information provided in this manual prior to use for proper performance, correct operation, and to ensure patient and operator safety. The user is responsible for operating this device as indicated in this guide.

Endosee Hysteroscope Part Number:

ES8000 (USA, Canada, Mexico):

Contains IEC plug type A

ES8000-INT (Rest of World):

Contains IEC plug types A, C, G, and I

Endosee Dx Cannula Part Numbers:

ESDX5 - Cannula with endometrial suction curette

ESDY5 - Cannula only

Endosee Cx Cannula Part Number:

ESCX5 - Cystoscopy Cannula

CooperSurgical

95 Corporate Drive • Trumbull, CT, 06611 • USA

Phone: (800) 243-2974 • Fax: (800) 262-0105

www.coopersurgical.com

Table of Contents

Warranty	4
1 General Description and Intended Use	5
1.1 General Description.....	5
1.2 Indications for Use.....	6
1.3 Contraindications.....	6
2 Warnings and Precautions	7
2.1 General.....	7
2.2 Inspection, Use and Disposal	8
2.3 Battery and Power Supply	10
2.4 Environmental.....	11
3 Symbol Legend	12
4 Preparation for Use	13
4.1 Unpacking and Inspection of Cannula	13
4.2 Power Requirements.....	13
4.3 Environmental Requirements.....	13
5 Description of Components	14
5.1 Product Description	14
5.2 Reusable Handheld Monitor	14
5.3 Charging/Docking Station	15
5.4 Sterile, Single-use Disposable Diagnostic (Dx) Cannula	16
5.5 Device Classification, Technical and Safety Specification.....	17
5.5.1 Technical Specifications.....	17
6 Basics of Operation and Procedure	18
6.1 Cautions for Handling and Use	18
6.2 Fluid Delivery System.....	18
6.3 Cannula Attachment and Detachment.....	19
6.4 Power Button Function.....	20
6.5 Brightness Button	20
6.6 Video/Single Frame Button	20
6.7 Touch Screen Functions	21
6.8 Patient Examination Procedure	24
6.9 Removal of the System	25
6.10 Disposal of Biohazardous Material	25
6.11 Battery Charging	26
7 External Interface Information	26
7.1 Charging/Docking Station	26
8 Maintenance	28
8.1 Cleaning/Disinfecting Instructions.....	28
8.2 Intermediate Level Combined Cleaning/Disinfecting Instructions.....	28
8.3 Low Level Combined Cleaning/Disinfecting Instructions	30
8.4 Disposable Cannula Handling/Care	31
8.5 Maintenance	31
9 Storage & Shipping	32
10 Technical Assistance	32
11 Troubleshooting Guide	32

Warranty

CooperSurgical, Inc. warrants that the device(s) supplied is free from defects in materials or workmanship. This warranty is valid only if the product is supplied to the end user by a CooperSurgical, Inc. approved agent or distributor and has been maintained in accordance with procedures documented in the User Manual. If failure occurs from manufacturing defects within 12 months of purchase, CooperSurgical, Inc. will repair or replace the defective item.

1 General Description and Intended Use

1.1 General Description

The Endosee Hysteroscope with the Disposable Diagnostic (Dx) Cannula is intended for use in viewing the cervical canal, uterine cavity, or female urinary tract, including the bladder, for the purpose of performing diagnostic procedures in an outpatient or in an office setting. Hysteroscopy can be used for the evaluation of the uterine cavity for pathology such as submucous myomas, endometrial polyps, malignancy and other conditions.

The Endosee Hysteroscope is a lightweight, handheld, battery operated portable system. It is used with a single-use Disposable Diagnostic (Dx) Cannula with a camera and light source at the distal end to illuminate the area for visualization and image and video capture. The video signal is electronically transferred to the main body of the hysteroscope via an electrical connector. An LCD touch screen display monitor on the hysteroscope is used for viewing.

Note: The Endosee Hysteroscope cannot be used with an external monitor.

System Components


Figure 1: Ready-to-use assembly

1.2 Indications for Use

The Endosee Hysteroscope and the Disposable Diagnostic (Dx) Cannula are used to permit viewing of the adult cervical canal and uterine cavity for the purpose of performing diagnostic procedures.

Generally recognized indications for diagnostic hysteroscopy include:

- Abnormal bleeding
- Infertility and pregnancy wastage
- Evaluation of abnormal hysterosalpingogram
- Intrauterine foreign body
- Amenorrhea
- Pelvic pain

The Endosee Hysteroscope and disposable diagnostic (DX) cannula can also be used to permit viewing of the adult urinary bladder through a minimally invasive approach by utilizing natural orifices to access the diagnostic site.

1.3 Contraindications

- Acute Pelvic Inflammatory Disease

Hysteroscopy may be contraindicated by the following conditions depending on their severity or extent:

- Inability to distend the uterus
- Cervical stenosis
- Cervical/Vaginal infection
- Uterine bleeding or menses
- Known pregnancy
- Invasive carcinoma of the cervix
- Recent uterine perforation
- Medical contraindication or intolerance to anesthesia

Cystoscopy is contraindicated in febrile patients with Urinary Tract Infections (UTIS) and those with severe Coagulopathy.

2 Warnings and Precautions

This section contains important safety information related to use of the Endosee Hysteroscope and the cannula. Other important safety information is repeated throughout this manual in sections that relate specifically to the precautionary information. Read all text surrounding all warning and precautionary information prior to performing any procedures with this equipment.

2.1 General

- Federal Law (USA) restricts the device to sale by or on the order of a physician.
- The Endosee Hysteroscope with the cannula is intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms.
- The Endosee Hysteroscope with the cannula is for use only by physicians with adequate training in hysteroscopy and/or cystoscopy.
- When a liquid distension medium is used, strict fluid intake and output surveillance should be maintained. Intrauterine instillation exceeding 1 liter should be followed with great care to reduce the possibility of fluid overload.
- Potential complications of continuous flow hysteroscopy:
 - Hyponatremia
 - Hypothermia
 - Uterine perforation resulting in possible injury to bowel, bladder, major blood vessels and ureter
 - Pulmonary edema
 - Cerebral edema
 - Infection
 - Bleeding
 - Pain

The use of normal saline as a distending medium, and limiting the infused volume to less than 1000 mL is recommended to decrease the risk of the above complications.

Intrauterine distension can usually be accomplished with pressures in the range of 35 - 75 mmHg.

Unless the systemic blood pressure is excessive, it is seldom necessary to use pressures greater than 75-80 mmHg.

- Vaginal ultrasonography before hysteroscopy may identify clinical conditions that will alter patient management.
- Since most urological endoscopy procedures utilize irrigation fluids, it is important that the correct fluid be chosen and that the volumes of inflow and outflow be strictly observed. Appropriate measures should be undertaken to prevent the possibility of fluid overload of the patient.

2.2 Inspection, Use and Disposal

- Inspect the integrity of the unit and condition before powering-on the fully assembled Endosee Hysteroscope. Do not use the system if indications of external damage are observed.
- The Endosee Hysteroscope contains no operator serviceable components within the enclosure. To avoid electrical shock or damage, the unit must not be disassembled. Doing so will void the warranty.
- Do not operate the Endosee Hysteroscope if any shipping damage or other defects to the equipment are noted during inspection. Immediately notify CooperSurgical, Inc. Customer Service if any defect is found.
- Do not use the cannula if the package has been damaged or the expiration date on the label has passed. The cannula is single use only. Do not reuse or resterilize the cannula.
- The cannula must be disposed of as biohazardous waste according to the safety guidelines of user facility/institution.
- Follow all cleaning procedures provided for the level of contamination of the device.

- The Endosee Hysteroscope can only be used with the LCD touch screen display on the Handheld Monitor. A separate video monitor cannot be used.
- Do not allow any fluids or substances to be spilled into or around the removable cannula connector on the Handheld Monitor, as damage to the electronics may occur.
- Do not misalign the connectors between the removable cannula and the Handheld Monitor during assembly, as damage may occur.
- When detaching the cannula from the Handheld Monitor, be sure that fluid does not contact the connector in the Handheld Monitor. Do not rotate or twist the cannula while removing.
- Do not alter the shape of the distal end of the cannula.
- No portion of the Endosee Hysteroscope should be opened for cleaning or disinfecting.
- Do not pour water or any cleaning solution directly onto Handheld Monitor during cleaning or disinfection. Doing so may damage electronic components and reduce product life.
- Do not immerse any component of the Handheld Monitor into solutions during the cleaning or disinfecting processes, as damage to the electronic components may occur.

2.3 Battery and Power Supply

Comply with the following to prevent damage to the battery:

- Use only the power supply provided with the device. Use of other chargers may damage the device and will void the warranty.
- Charge the battery at room temperature.
- Do not charge the battery near a heat source.
- Do not disassemble or alter the battery in any way.
- Battery is not user-replaceable.
- The Handheld Monitor should be fully charged before use. A fully charged battery will provide at least 2 hours of continuous operating time under normal operating conditions.
- Fully charge the Handheld Monitor after each procedure.
- Refer to battery power indicator on the touch screen monitor to determine remaining battery power.
- Do not attempt to use any power source or battery not directly provided by CooperSurgical, Inc.

Comply with the following to reduce the risk of burns, fire, electric shock, or injury to persons:

- Do not operate the power supply if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or if it has fallen into water.
- Do not charge the battery where oxygen is being administered.
- Do not unplug power supply by pulling on the cord. To unplug, pull on the plug, not the cord.

2.4 Environmental

- Failure to ensure that all of these environmental requirements are met may result in improper performance or damage to the electronic components within the Endosee Hysteroscope.
- Due to the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare environments (for example cellular phones, etc.) it is possible that high levels of such interference due to close proximity or strength of a source may result in a disruption of the performance of this device.
- The Endosee Hysteroscope is not designed for use in environments in which strong interference may impact the performance of the equipment. If this occurs, the site of use should be surveyed to determine the source of this disruption, and the following actions may be taken to eliminate the source:
 - Remove, reorient or relocate the interfering equipment
 - Increase the separation between the interfering equipment and the Endosee Hysteroscope
 - Incrementally turn off equipment in the vicinity to identify the interfering device
- Do not attempt to transport or ship the Endosee Hysteroscope without using proper packaging to protect the product.

3 Symbol Legend

The following symbols are used in this manual, and/or on the product labeling.

Reorder number

Batch code

Serial number

Sterilized using ethylene oxide.

U.S. Federal law restricts this device to sale by or on the order of a licensed practitioner.

Do not resterilize.

WARNING:

This symbol is used to indicate information that may relate to serious adverse reactions and potential safety hazards, limitations in use imposed by them, and steps that should be taken if they occur.

CAUTION:

This symbol is used to indicate information regarding any special care to be exercised by practitioner and/or patient for the safe and effective use of the equipment.

Consult instructions for use.

Do not re-use.

Nonsterile

Not made with natural rubber latex.

Do not use if package is damaged.

BF Applied Part

Power On/Off Switch

Battery

Manufacturer

MEDICAL EQUIPMENT
WITH RESPECT TO
ELECTRICAL SHOCK,
FIRE AND MECHANI-
CAL HAZARD ONLY IN
ACCORDANCE WITH
AAMI/ANSI ES 60601-1
(2005), IEC 60601-2-18
AND CAN/CSA C22.2
NO. 60601.1:08

4 Preparation for Use

4.1 Unpacking and Inspection of Cannula

Follow below to unpack, inspect, and install the cannula and the Endosee Hysteroscope. Each single-use, disposable cannula comes sealed in a pouch and is Ethylene Oxide (ETO) sterilized (Figure 2). Check that the expiration date has not passed before use.

Figure 2: Cannula in a pouch

- **Packing Materials:** Save the packing materials for potential future transportation and storage of the Endosee Hysteroscope.
- **Inspection:** Inspect all components for damage during shipment, or discrepancies upon arrival.

4.2 Power Requirements

The Endosee Hysteroscope operates on DC voltage, from an internal battery source. The Handheld Monitor must be fully charged before the initial use.

- Place the Handheld Monitor on Charging/Docking Station, and connect the two.
- Attach power supply to the Charging/Docking Station.
- Plug the power supply to a wall outlet.
- Continue to charge until the battery power indicator on the touch screen monitor shows 100%, indicating that the battery is fully charged. The Handheld Monitor is ready for use.

4.3 Environmental Requirements

All of the environmental requirements listed in the Warnings and Precautions section must be met in the location where the Endosee Hysteroscope is intended to be operated for patient use. If ANY of the seen environmental parameters are present, the Endosee Hysteroscope should not be operated until modifications to the environmental conditions are met.

5 Description of Components

5.1 Product Description

The Endosee Hysteroscope is comprised of the following components:

5.2 Reusable Handheld Monitor

The Handheld Monitor contains electronics including a power on/off button, a video processor, a display unit (LCD display), a rechargeable battery, battery management electronics, microcontrollers, storage for images and video, and firmware. The device uses a rechargeable battery. The battery is charged when the Handheld Monitor is connected to the Charging/Docking Station and the Charging/Docking Station is plugged into a wall outlet. A power supply is provided with the device. Figure 3 shows the reusable Handheld Monitor. The detailed battery requirements are listed in section 4.2.

1. LCD Touch Screen Monitor
2. Power On/Off Button
3. Brightness Button
4. Video/Single Frame Button

Figure 3: Reusable Handheld Monitor

5.3 Charging/Docking Station

The Charging/Docking Station has dual functions. It allows the battery to be charged using the power supply and allows video recordings and single frame photographs to be downloaded from the internal memory of the Handheld Monitor to an external computer.

WARNING: Use only the power supply provided with the device. Use of other chargers may damage the device and will void the warranty.

WARNING: Charge the battery at room temperature.


WARNING: Do not charge the battery near a heat source.

Figure 4: The Handheld Monitor on the Charging/Docking Station

5.4 Sterile, Single-use Disposable Diagnostic (Dx) Cannula

The cannula contains a miniature CMOS camera and light-emitting diode (LED) illumination module at the tip. The cannula connects to the Handheld Monitor through an electrical connector for: DC power, image data transfer and removal and control of the CMOS camera and illumination module. The cannula has an inner channel for the infusion of fluid. Fluid infusion is achieved by attaching either a pre-filled syringe or a pressurized bag of distending medium to the inflow port through an extension tube or a tubing set. The outflow port is at the distal tip of the cannula and empties around the cannula and out the vaginal vault. The cannula has a pre-curved tip section at an angle of $25^\circ \pm 5^\circ$. The cannula tip section should not be altered. The cannula maximum temperature will not exceed 41°C. Figure 5 shows the cannula.

1. Handheld Monitor Connector
2. Inflow Port
3. Tip Section
4. CMOS Camera & LED Illumination Module

Figure 5: Disposable Diagnostic (Dx) Cannula

5.5 Device Classification, Technical and Safety Specification

5.5.1. Technical Specifications

HANDHELD MONITOR

Length (Handheld Monitor): 117.5 mm (4.6 in)

Diameter (Handheld Monitor): 29.8 mm (1.2 in)

Weight (cannula, battery, display): 90 g (0.2 lb)

STERILE, SINGLE-USE CANNULA

Length: 276 mm (10.9 in)

Outside Dimensions: 4.3 mm (0.2 in)

Pre-fixed tip angle: 25° ± 5°

LCD DISPLAY PANEL

Diagonal Size of Display Area: 88.9 mm (3.5 in)

Thickness: 13.2 mm (0.5 in)

Horizontal: 93.0 mm (3.7 in)

Height: 72.0 mm (2.8 in)

POWER ADAPTER REQUIREMENTS

Input Power Range: 100-240V AC/0.5A

Frequency: 50/60Hz

Output: 5.0V DC/2.0A

Charging Connector: AC Adapter incorporated into Charging/Docking Station

BATTERY SPECIFICATIONS

Battery Type: AA Lithium-ion rechargeable battery

Capacity: 3.7V rating, 2500mAh

Maximum Voltage when Battery is 4.2V fully charged:

Voltage when system power is turned off: ≤ 3.7V

Full Charge Time: Approx. 4 hours

Operating Time: >2 hours

SITE REQUIREMENTS — OPERATING CONDITIONS

Operating Temperature: Room Temperature

6 Basics of Operation and Procedure

6.1 Cautions for Handling and Use

The following cautions should be strictly observed in the handling and use of the Endosee Hysteroscope.

WARNING: *The Endosee Hysteroscope is for use only by physicians with adequate training in hysteroscopy/cystoscopy.*

6.2 Fluid Delivery System

Fluid infusion is achieved by attaching either a pre-filled syringe or a pressurized bag of distending medium to the inflow port of the cannula through an extension tube or a tubing set.

WARNING: *When a liquid medium is used, strict fluid intake and output surveillance should be maintained. Intrauterine instillation exceeding 1 liter should be followed with great care to reduce the possibility of fluid overload.*

WARNING: *Potential complications of Continuous Flow Hysteroscopy include:*

- *Hyponatremia*
- *Hypothermia*
- *Uterine perforation resulting in possible injury to bowel, bladder, major blood vessels and ureter*
- *Pulmonary edema*
- *Cerebral edema*
- *Infection*
- *Bleeding*
- *Pain*

The use of normal saline as a distending medium, and limiting the infused volume to less than 1000 mL is recommended to decrease the risk of the above complications. Intrauterine distension can usually be accomplished with pressures in the range of 35-75 mmHg. Unless the systemic blood pressure is excessive, it is seldom necessary to use pressures greater than 75-80 mmHg.

CAUTION: Do not allow any fluids or substances to be spilled in or around the cannula connector on the Handheld Monitor, as damage to the electronics may occur.

6.3 Cannula Attachment and Detachment

In preparation for use, attach the cannula to the electrical connector on the Handheld Monitor per the following:

1. Peel apart the proximal end of the package and remove the cannula.
2. Rotate the cannula to a position where the “red dot” marking on the cannula aligns with the matching dot on the Handheld Monitor. Press it fully into place until a subtle “click” of the connector detent locking mechanism is felt.

CAUTION: Do not misalign the connectors between the cannula and the Handheld Monitor during assembly, as damage may occur.

Figure 6: Cannula Attachment

6.4 Power Button Function

The function of the Handheld Monitor Power Button on the Hysteroscope is described below.

WARNING: *Inspect the integrity of the unit and condition before powering-on the fully assembled unit. Do not use the system if visual signs of external damage are observed.*

Figure 7: Power Button

Touch Function	Action
Power On	Press-and-Release
Power Off	Press-Hold (1s) and Release

6.5 Brightness Button

The LED Brightness control is a cycle-through push button. There are four (4) steps of LED brightness adjustments.

Figure 8: Brightness Button

6.6 Video/Single Frame Button

The Video/Single Frame Button allows user to save still images and video images. The image capture operations are described below.

Figure 9: Video/Single Frame Button


Video/Single frame Touch Function	Action
Capture Still Image	Press-and-Release
Capture Video Clip	Press-and-Hold: See Page 22
Stop Video Capture	While system is capturing video, press button one more time

6.7 Touch Screen Functions

Touch Screen Description/Touch Screen Action


HOME MENU

Home menu has a battery power usage indicator and (4) four soft buttons each allowing a discrete function. Pressing each soft button enters a sub menu for a specific function. The instructions for touch screen user menus are described in this section.

Battery Power Usage Indicator

Battery usage indicator provides the battery power usage status. When battery power usage reaches 30%, user should charge the Handheld Monitor.

New Patient Menu


- Press the appropriate numbers to enter patient ID. The number is displayed in the window with a white background. User may enter only four (4) digits.
- Press to return to **Home Menu** and no Patient ID is entered.
- Press **C** to cancel previously entered number.
- Press **OK** to confirm Patient ID and start the procedure. A confirmation message "Patient ID# is entered" along with the newly created Patient ID will be displayed on the LCD screen.
- After three (3) seconds, the confirmation message is erased and the real time image will be displayed on the LCD Screen.

Touch Screen Description/Touch Screen Action




Preview Menu

- Pressing the **Preview** soft button in the **Home Menu** starts to display real time image. A green circle appears in the upper right corner of the display.
- Press the **Video/Single frame Button** on the Handheld Monitor (Refer to 6.6 for detailed instructions).
- Press **Right Arrow** soft button entering **Playback Menu** or Press  to return to **Home Menu**.

Playback Menu


- Press the Patient ID; the selected patient ID is highlighted.
- Press **Up Arrow** or **Down Arrow** to select another Patient ID.
- Press highlighted Patient ID a second time; the captured image will be displayed.
- User can also double click the Patient ID to directly review the image.
- Press **DEL** to delete a Patient ID. Before a Patient ID can be deleted, a confirmation screen is displayed. In the confirmation menu, press **Yes** to delete and **No** to cancel Patient ID deletion.
- **Reviewing Saved Images/Video.** After selecting the patient ID in the **Playback** menu, you may review the saved images/video.
- Press **Left Arrow** to select previous image; press **Right Arrow** to select next image.
- Press **Garbage Can** soft button to delete an image. Before an image can be deleted, a confirmation screen is displayed. In confirmation menu, press **Yes** to delete and **No** to cancel image deletion.
- Press  to return to **Home Menu**

Touch Screen Description/Touch Screen Action

Setup Menu

- Press **Date & Time** soft button to enter the date and time.
- Press **NTSC** or **PAL** soft button to select display format. The format chosen is highlighted in blue.
- Press **Format** button to format the SD card in the Handheld Monitor.
- Press **Language** button to change the system language.
- Press  to return to **Home Menu**.

Date and Time

- Press one of the following white soft buttons; the white window edge will be highlighted with red lines.
- Press **Down Arrow** or **Up Arrow** to increase or decrease digits.
- Press **OK** soft button to save settings.
- Press  to cancel the changes and return to **Setup Menu**.


Language Settings

- Language settings must be set at initial boot up of the Endosee Hysteroscope.
- Select the language by pressing on the button for the specific language.

Touch Screen Description/Touch Screen Action

- Confirm your selection by pressing **YES** or cancel by pressing **NO**.

- The language settings can also be changed at any other time by pressing the **Language** button within the Setup Menu screen.

6.8 Patient Examination Procedure

The following steps describe the process for using the Endosee Hysteroscope.

CAUTION: A fully charged battery will provide at least 2 hours of continuous operating time under normal operating conditions. The Handheld Monitor should be fully charged before initial use. It is recommended to fully charge the Handheld Monitor after each procedure.

CAUTION: Refer to battery power indicator on the touch screen monitor to determine remaining battery power.

- Ensure that the Handheld Monitor is thoroughly cleaned and disinfected using the procedures, as described in section 8.
- Ensure that the Handheld Monitor is charged, as described in section 6.7.
- Remove the sterile cannula from the sterilized pouch.
- Attach the cannula to the Handheld Monitor using the procedures, as described in section 6.3.
- Attach the fluid tubing to the cannula, as described in section 6.2.

- Press and Release POWER ON/OFF button.
- The light source at the distal tip of the cannula will now be illuminated. If the tip of the cannula does not illuminate, system hardware troubleshooting procedures should be performed, as described in section 11.
- Flush the fluid through the cannula until all air bubbles have cleared from the tubing.
- Using the Handheld Monitor touch screen, enter the patient ID or touch the preview soft button as described in section 6.7.
- The Handheld Monitor is now ready for use.
- The cannula comes with a pre-curved tip at an angle of $25^\circ \pm 5^\circ$. The curve of the tip is similar to that of most conventional cervical dilators.
- The cannula is inserted through the cervical canal or uretha, with fluid flowing, while viewing the image as is standard procedure for introduction of a hysteroscope, or cystoscope.

6.9 Removal of the System

- After the procedure is completed, remove the hysteroscope, and Press-hold (1s), and release the Handheld Monitor POWER ON/OFF.
- Hold the Handheld Monitor vertically so the cannula is below the Handheld Monitor. Gently pull the cannula straight out from the Handheld Monitor. The cannula will detach from the Handheld Monitor.

CAUTION: When detaching the cannula from the Handheld Monitor, be sure that fluid does not contact the connector in the Handheld Monitor. Do not rotate or twist the cannula while removing.

6.10 Disposal of Biohazardous Material

WARNING: The used cannula must be disposed of as biohazardous waste according to the safety guidelines of user facility/institution.

6.11 Battery Charging

The Endosee Hysteroscope has an internal battery that powers the unit during use. A fully charged battery will provide at least 2 hours of continuous operating time under normal operating conditions.

To recharge the battery, follow the procedure in section 4.2

7 External Interface Information

7.1 Charging/Docking Station

The Charging/Docking Station provides the capability for battery charging and image uploading.

Figure 10: Handheld Monitor on the Charging/Docking Station

DOWNLOADING STORED IMAGES

Download to a PC

1. Plug Charging/Docking Station charge cable to Handheld Monitor.
2. Connect the USB cable to the computer and Charging/Docking Station.
3. Turn on Handheld Monitor power.
4. Boot up computer.
5. Click Start, and then click My Computer.
6. The Handheld Monitor shows up as a removable disk in "My Computer."
7. Left click on that removable disk to find a folder named with a patient ID.
8. Click folder named with patient ID.

- 9.** Click folder with exam date.
- 10.** The files in exam date folder are recorded pictures, and their names end in “jpg.”
- 11.** Copy those files by selecting them and copying them. The keyboard shortcut (control-A) marks all the photos in a folder and the keyboard shortcut (control-C) copies them to the computer’s clipboard.
- 12.** Navigate to the folder where you want to save patient pictures and paste (control-V) your photos into place.
- 13.** After you copy files or folders to the hard drive, it is useful to view the hard drive to confirm that the files are copied.
- 14.** Do not erase photos from your memory card until you are positive you’ve saved them on the hard drive!

Download to a Mac

- 1.** Plug Charging/Docking Station charge cable to Handheld Monitor.
- 2.** Connect the USB cable to the computer and Charging/Docking Station.
- 3.** Turn on Handheld Monitor power.
- 4.** Boot up computer.
- 5.** Click the Finder on the application dock.
- 6.** The Handheld Monitor shows up as a removable disk under “Devices.”
- 7.** Left click on that removable disk to find a folder named with a patient ID.
- 8.** The files in that folder are recorded pictures, and their names end in “jpg.”
- 9.** Copy those files by selecting them and copying them. The keyboard shortcut (command-A) marks all the photos in a folder and the keyboard shortcut (command-C) copies them to the computer’s clipboard.
- 10.** Navigate to the folder where you want to save patient pictures and paste (command-V) your photos into place.

11. After you copy files or folders to the hard drive, it is useful to view the hard drive to confirm that the files are copied.
12. Do not erase photos from your memory card until you are positive you've saved them on the hard drive!

8 Maintenance

8.1 Cleaning/Disinfecting Instructions

After each patient use, follow the combined cleaning/disinfecting instructions provided below for the Handheld Monitor. If the Handheld Monitor was contaminated with blood or body fluids during use, you must follow the intermediate level instructions. If the Handheld Monitor was not contaminated with blood or body fluids, follow the low level instructions.

WARNING: *Follow all cleaning procedures provided for the level of contamination of the device.*

CAUTION: *No portion of the Endosee Hysteroscope should be opened for cleaning or disinfecting.*

8.2 Intermediate Level Combined Cleaning/Disinfecting Instructions

- After each patient use, power off and unplug all electric connections.
- Place Handheld Monitor on a clean, dry, disinfected or sterile basin or surface.

CAUTION: *Throughout the cleaning/disinfection process, orient the connector in a manner that prevents liquid from entering the electrical connection area of the connector (e.g., with the opening to the electrical connections facing upward). Do not wet the opening or inside of connector housing the electrical connection components.*

CAUTION: When cleaning/disinfecting the LCD screen of the Handheld Monitor, pay careful attention to the seams around the edge of the display and the back of the screen where it attaches to the grip. When cleaning the grip section, pay careful attention to the seams around the edge of the button panels and labels.

- Follow the instructions below to clean and disinfect the handheld monitor. These instructions were validated using the following wipes: CaviWipes®, PDI SANI-CLOTH® AF3, and Clinell® universal sanitizing wipes. If you choose a wipe not listed here, you must validate it before use.
- Using a fresh, moist disinfectant wipe, thoroughly wipe the front and back of the LCD screen of the Handheld Monitor. Discard the wipe. Repeat with a fresh wipe for a minimum of thirty (30) seconds. Discard the wipe.
- Using a fresh, moist wipe, thoroughly wipe the grip section of the Handheld Monitor. Discard the wipe. Repeat with a fresh wipe for a minimum of thirty (30) seconds. Discard the wipe.
- Using a fresh, moist wipe, thoroughly wipe the seams of the Handheld Monitor for one (1) minute.
- Using a fresh, moist wipe, thoroughly wipe all surfaces of the Handheld Monitor for a minimum of ten (10) seconds. Discard the wipe.
- Allow the wiped and wetted Handheld Monitor to stand for a minimum of three (3) minutes. (For SANI-CLOTH, use additional wipes if needed to ensure the device is visibly wet for three (3) continuous minutes.)
- At the end of the three (3) minutes, use a sterile swab saturated with 70% Isopropyl alcohol (IPA) to thoroughly swab the seams around the edge of the button panels for a minimum of fifteen (15) seconds.

- Use a sterile 70% IPA wipe, or a sterile, lint-free wipe, saturated with sterile 70% IPA, to thoroughly wipe the entire Handheld Monitor for a minimum of thirty (30) seconds.
- Allow the device to air dry. Once dry, visually inspect the Handheld Monitor for any residual soil. Repeat cleaning/disinfection instructions if visible soil is present.
- Place the cleaned/disinfected device in an appropriate dry storage area.

8.3 Low Level Combined Cleaning/Disinfecting Instructions

- After each patient use, power off and unplug all electric connections.
- Place Handheld Monitor on a clean, dry, disinfected or sterile basin or surface.

CAUTION: Throughout the cleaning/disinfection process, orient the connector in a manner that prevents liquid from entering the electrical connection area of the connector (e.g., with the opening to the electrical connections facing upward). Do not wet the opening or inside of the connector housing the electrical connection components.

CAUTION: When cleaning/disinfecting the LCD screen of the Handheld Monitor, pay careful attention to the seams around the edge of the display and the back of the screen where it attaches to the grip. When cleaning the grip section, pay careful attention to the seams around the edge of the button panels and labels.

- Using a fresh, moist 70% Isopropyl Alcohol (IPA) wipe (minimum size of 5.6in x 7in), thoroughly wipe the front and back of the LCD screen of the Handheld Monitor for one (1) minute. Discard the wipe.

- Using a fresh, moist IPA wipe, thoroughly wipe the grip section of the Handheld Monitor for one (1) minute. Discard the wipe.
- Using a fresh, moist IPA wipe, thoroughly wipe the seams of the Handheld Monitor for one (1) minute. Discard the wipe.
- Using a fresh, moist IPA wipe, thoroughly wipe all surfaces of the Handheld Monitor for one (1) minute. Discard the wipe.
- Allow the device to air dry. Once dry, visually inspect the Handheld Monitor for any residual soil. Repeat cleaning/disinfection instructions if visible soil is present.
- Place the cleaned/disinfected device in an appropriate dry storage area.

8.4 Disposable Cannula Handling/Care

The cannula is a single-use, disposable item. As such, no formal cleaning or disinfection procedures apply to the cannula.

8.5 Maintenance

The Endosee Hysteroscope has no user - serviceable parts.

Note: The Endosee Hysteroscope requires no calibration or service in the field. All electronics are tested and validated for performance at the manufacturer. If the performance of the Endosee Hysteroscope is in doubt, contact CooperSurgical, Inc. Customer Support for assistance.

Maintenance Intervals

Basic preventive maintenance of the Endosee Hysteroscope is an important function that should be performed on schedule to ensure safe and effective operation.

If the integrity or performance of any portion of the Endosee Hysteroscope is in doubt, the system should NOT be used until the issue is resolved.

The following maintenance intervals are recommended:

Maintenance Type	Interval	Personnel
General Inspection	Before each Patient Use	Users
Safety Label Integrity	As Required	All

9 Storage & Shipping

The Endosee Hysteroscope is shipped in a protective box. Do not dispose of protective box or any other shipping materials. These items should be retained for future storage or transport of the equipment.

CAUTION: *Do not attempt to transport or ship the Endosee Hysteroscope without using proper packaging to protect the product.*
Store at room temperature.

10 Technical Assistance

For technical information, assistance or to order parts or manuals, contact CooperSurgical, Inc. Customer Support Department, or your local CooperSurgical, Inc. representative.

Information regarding system serial number is available from Customer Support.

CooperSurgical

95 Corporate Drive • Trumbull, CT 06611 • USA

Phone: (800) 243-2974 • Fax: (800) 262-0105

www.coopersurgical.com

11 Troubleshooting Guide

The information in this section is intended to provide simple steps that can be performed by a user for identifying the primary cause, and possible simple solutions that can be resolved on-site, to basic problems that may be encountered while operating the Endosee Hysteroscope.

Any issues determined to be beyond the scope of the basic user troubleshooting steps provided in this User Manual should be communicated to the CooperSurgical, Inc. Customer Support department.

WARNING: *The Endosee Hysteroscope contains no user-serviceable components within the enclosure. To avoid electrical shock or damage, the unit must not be disassembled. Doing so will void the warranty.*

WARNING: *Do not disassemble or alter the battery in any way.*

WARNING: *Battery is not user replaceable.*

Problem	Test	Action
(1) When switched ON, there is no display on screen	Is the cannula tip LED illuminated?	<ul style="list-style-type: none">• If NO, change to new cannula• If YES, contact CooperSurgical, Inc. Customer Service.
(2) Poor picture quality	N/A	<ul style="list-style-type: none">• Clean cannula tip with a sterile, clean, soft wipe.• Ensure cannula is fully inserted in the Handheld Monitor <p><i>If steps above do not correct picture quality, contact CooperSurgical, Inc. Customer Support.</i></p>
(3) Image on screen flickers or has lines across it	N/A	<ul style="list-style-type: none">• Ensure cannula is fully inserted in the Handheld Monitor <p><i>If this does not correct picture quality, contact CooperSurgical, Inc. Customer Support.</i></p>
(4) Loose component, poorly fitting connection	N/A	<ul style="list-style-type: none">• Do not use device. Contact CooperSurgical, Inc. Customer Support for product return.

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© 2018 CooperSurgical, Inc.
SANI-CLOTH® is a registered trademark of PDI, Inc.
CLINELL® is a registered trademark of Gama Healthcare
CaviWipes® is a registered trademark of Metrex Research Corporation.
CooperSurgical, Inc. • 95 Corporate Drive • Trumbull, CT 06611 USA