

C SAFE - Uterine Safety Scalpel

Performing a C-section, particularly an emergent procedure, can pose several challenges. Accidental and preventable fetal laceration should not be one of them. C SAFE is designed to make fetal lacerations a "never event." The product features:

- Distal tip that ensures blunt entry into the lower uterine segment
- Shielded blade that cuts up and away from the fetus to prevent accidental nicks or cuts

C SAFE helps reduce the risk of laceration posed by traditional scalpels and provides added safety during C-section for baby, mother and the surgical team.

TRADITIONAL SCALPEL



A traditional scalpel exposes the fetus to the risk of laceration.

C SAFE SAFETY SCALPEL



C SAFE is designed to shield the fetus from laceration.





Shielded surgical blade cuts up and away from the fetus to prevent the risk of laceration.



Specially designed distal tip enables blunt entry into the lower uterine segment.

Ergonomic, ridged handle ensures firm grasp under any intraoperative conditions.

Quick Guide to Using C SAFE Safety Scalpel







- 1. Prep the patient and expose uterine wall using standard abdominal incision techniques.
- 2. Make initial entry into uterine wall with C SAFE distal tip using an oblique or slightly angled position.
- 3. Immediately rock the C SAFE handle horizontally to point distal tip in direction of intended incision and slightly up toward the inner uterine wall.

Ordering Information

Product#	Contents
BX101	C SAFE Safety Scalpel – Box of 10
900-517	OB/Mobius Medium-Large Retractor, 17 cm l.D. ring, for incisions to 15 cm $-$ Box of 5
900-520	OB/Mobius Large Retractor, 22 cm I.D. ring, for incisions to 17 cm — Box of 5

The Ideal C-Section Safety Kit:

Pair C SAFE with the OB/Mobius® C-Section Retractor to achieve maximum exposure and an unobstructed surgical field.



OB/Mobius®

We Can Help Bring C SAFE Into Your Hospital

- Identify a key champion or decision-maker within your hospital.
- We will work with key decision-makers in your facility to provide Value Analysis Committee or other technical documentation required for approval of the device in your hospital system.
- Set up a C SAFE product evaluation in your hospital through one of CooperSurgical's highly-trained OB Specialists.

CoperSurgical

CooperSurgical is the leading manufacturer of innovative medical devices and procedure oriented solutions that advance the standard of care. Our highly reliable and clinically relevant products, as well as nearly 150 clinically trained utilization sales specialists facilitate the delivery of enhanced outcomes for your patients, regardless of the clinical setting.

See us on YouTube at www.YouTube.com/CooperSurgical.

800.243.2974 | 203.601.5200 | www.coopersurgical.com 95 Corporate Drive, Trumbull, CT 06611

©2015 CooperSurgical, Inc.

82560 Rev. 1/15