RUMI® Uterine Manipulator Tips

for use with the

RUMI® and Advincula Arch™ Handles

Instructions for Use (English)

STERILE • SINGLE-USE

R_x Only

CE 0086

UMY514 • UML516 • UMW676 • UMB678 • UMG670 • UMO672

DEVICE DESCRIPTION


The RUMI® Uterine Manipulator Tip (Tip) is a sterile, single-use device, which is attached to a reusable RUMI, RUMI Arch or Advincula Arch™ Uterine Manipulator Handle. The two components together create The RUMI System®, which provides uterine manipulation.

There are six different size Tips available, which are pictured above and listed in the table on the following page.

Each single use RUMI Uterine Manipulator Tip is manufactured with a USP Class VI (medical grade) silicone body with a balloon and a stainless steel internal support rod. When inflated properly, the balloon secures the Tip in the uterus.

Two silicone Catheters extend from the Tip. Both Catheters have a Luer Lock fitting for syringe attachment. The Balloon Catheter has a hose clamp to retain balloon pressure. The other Catheter is for infusing dye into the uterus.

Not made with natural rubber latex.

To accommodate different uterine sizes, the tips are available in the following shaft diameters and shaft lengths:

Part Number	Color	Diameter	Length	
UMY514	Yellow	5.1 mm	3.75 cm	
UML516	Lavender	5.1 mm	6 cm	
UMW676	White	6.7 mm	6 cm	
UMB678	Blue	6.7 mm	8 cm	
UMG670	Green	6.7 mm	10 cm	
UMO672	Orange	6.7 mm	12 cm	

For further Information on how to use The RUMI® System, please refer to the appropriate Handle Directions for Use.

WARNINGS

- **NEVER** use a tip that is longer than the sounded depth of the uterus.
- **NEVER** attempt to use the device *in utero* without first performing the balloon test inflation procedure. **Inflate the balloon with sterile saline only. DO NOT** use air for inflation of the intrauterine balloon.
- **NEVER** use the device with the intrauterine balloon deflated.
- Contents supplied sterile. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

CAUTIONS

- U.S. Federal law restricts this device to sale by or on the order of a physician.
- **THE TIP IS SINGLE-USE ONLY. NEVER REUSE A RUMI TIP.**

INDICATIONS FOR USE

The Tip is intended to be used only with the RUMI, RUMI II, RUMI Arch® and Advincula Arch™ Handles. Please refer to the appropriate Handle Directions for Use for further indications.

CONTRAINDICATIONS

DO NOT use if the balloon will not inflate during test inflation, or if sterility has been compromised. For further contraindications, please refer to the appropriate Handle Directions for Use.

PRECAUTIONS

Dilate the cervix using standard technique, as required, to accommodate the tip diameter selected. Lubricate the tip with sterile solution to ease patient insertion if necessary.


INSTRUCTIONS FOR USE

Tip Selection


1. Sound the uterus for depth and direction.
2. Select the Tip which is less than or equal to the sounded depth of the uterus.

Tip Attachment

1. Using aseptic technique, grasp the base of the tip and handle shaft as shown (photo 1).
2. Insert the tip support rod into the hole in the snap drum.
3. With a clockwise twisting motion of the tip, push the tip base completely over the snap drum nose (photos 2 & 3).
4. Visually inspect to ensure that the base is seated completely and that the catheters exit to the sides. A slight tug on the tip will verify that it is on securely (photo 4).


1. Insert rod into Snap Drum

2. Twist clockwise while pushing on Tip Base

3. Snap into place

4. Pull to verify

Balloon Test Inflation Procedure

Inflate and fully deflate the balloon with 5 cc (ml) of saline prior to patient insertion.

Insertion

1. Use a sterile instrument to grasp the lip of the cervix.
2. Sufficiently dilate the cervix to Hegar/Hank 8 (French 24) to ease Tip insertion.
3. Following the direction the uterus was sounded, SLOWLY insert the Tip into the uterus until the Tip Base meets flush against the cervix.
4. While maintaining a slight pressure of the Tip base against the cervix, SLOWLY inject 3 cc (ml) to 5 cc (ml) of sterile saline into the balloon, close the hose clamp and then remove the syringe.
5. Gently pull on the device to test for proper engagement of the Tip in the uterus.
6. If dye injection is desired, attach a dye syringe to the dye catheter. Closely follow the media manufacturer's instructions.

Tip Disengagement

1. Open hose clamp and withdraw the saline from the balloon.
2. Using care, gently withdraw the tip from the cervix and remove from the patient.
3. Remove the tip from the Handle (using a counter clockwise motion) and dispose of in accordance with all applicable Federal, State and local Medical / Hazardous waste practices.

EXPLANATION OF SYMBOLS

REF Reorder number

LOT Batch code

 Use-by date

 Do not re-use

 Caution

 Consult instructions for use

 Do not use if package is damaged

 Keep dry

 Manufacturer

 Not made with natural rubber latex.

STERILE EO Sterilized using ethylene oxide

R_xOnly U.S. Federal law restricts this device to sale by or on the order of a physician.

EC REP Authorized Representative in the European Community.

Marks bearing the symbol TM and © are trademarks and registered trademarks of CooperSurgical, Inc.

CooperSurgical is a registered trademark of CooperSurgical, Inc.

© 2018 CooperSurgical, Inc. All rights reserved.

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CooperSurgical é uma marca registada da CooperSurgical, Inc.
© 2018 CooperSurgical, Inc. . Todos os direitos reservados.