Vaginal Arcing and Omnidflex Diaphragm Fitting Kit

Instructions for Use (English)

Product No: MXWSFS Arcing Diaphragms
MXWFFS Omnidflex Diaphragms
Sizes: 65 mm thru 85 mm

CE0086

The diaphragms are made of Medical Grade silicone, are latex free, and are autoclavable. Thoroughly rinse food grade powder off fitting diaphragm before initial use. Cleanse the diaphragm before initial use by washing it with a mild, non-perfumed detergent and warm water, rinsing and drying it carefully.

DEVICE DESCRIPTION

- **The Silicone Arcing Diaphragm** is a molded, flesh-colored, latex-free, silicone diaphragm which is compressible on one plane only. It is folded for insertion by bringing the two notched sites on the inner lip together.
- **The Silicone Omnidflex Diaphragm** is a molded, flesh-colored, latex free, silicone non-arcing diaphragm which can be bent for insertion, at any point.

Arcing and Omnidflex Diaphragms, in conjunction with an approved spermicidal gel or cream are indicated for the prevention of pregnancy in women who elect to use diaphragms as a method of contraception.

A properly fit Arcing and Omnidflex Diaphragms remains in place not only by virtue of tension between the spring and vaginal wall, but also because of suction created by the close contact of the lip of the Arcing and Omnidflex Diaphragms with the vaginal mucosa. Because this double action helps assure the diaphragm will remain in place, it may be possible to fit the patient with a smaller size should she fall between two sizes. This helps avoid the discomfort which occurs because of too much pressure on the urethra and bladder neck, when fitted with too large a size.

How Supplied

The Arcing and Omnidflex Diaphragms are available in sizes 60 mm thru 95 mm, in 5 mm increments.

The Fitting Kits include a set of diaphragms from 65 mm to 85 mm, in 5 mm increments.

INDICATIONS FOR USE

The Arcing and Omnidflex Diaphragms, in conjunction with an approved spermicidal gel or cream are indicated for the prevention of pregnancy in women who elect to use diaphragms as a method of contraception.

WARNINGS

- **IMPORTANT** - For contraceptive effectiveness, the diaphragm should remain in place for six to eight hours (depending on the spermicide that is used) after intercourse and it should be removed thereafter. Continuous wearing of a contraceptive diaphragm for more than twenty-four hours is not recommended. Removal of the diaphragm before six to eight hours may increase risk of becoming pregnant.
- Wearing the diaphragm for any period of time may encourage the growth of certain bacteria in the vaginal tract. These bacteria can lead to symptoms of Toxic Shock Syndrome (TSS). An association has been reported between diaphragm use and Toxic Shock Syndrome (TSS), a serious condition which can be fatal.

SYMPTOMS OF TOXIC SHOCK SYNDROME (TSS)

1. Sudden high fever (102 °F/39 °C or over).
2. Nausea, vomiting and/or diarrhea.
3. Fainting, feeling faint or dizziness when standing up.
4. Rash over much of body (may resemble sunburn).
5. Weakness, tenderness and/or aching or pain of muscles and/or joints.
6. Sore throat, red eyes, and general feeling of weakness.

Instruct your patients that if they experience a sudden high fever and one or more of the above symptoms, the diaphragm should be removed and to contact their healthcare professional immediately.

Diaphragm users should be instructed to consult their physician or other healthcare provider:

1. If they or their partner feel, or are made uncomfortable by the presence of the diaphragm.
2. If they experience any discomfort or pain while the diaphragm is in place. This may be due to incorrect diaphragm insertion, an abnormal pelvic condition, constipation or incorrect diaphragm size.
3. If the diaphragm slips out of place when walking, coughing, sneezing or straining.
4. If the diaphragm no longer fits snugly above the pubic bone.
5. If at times other than menstruation there is blood on the diaphragm when it is removed.
6. If there are any holes, tears or other deterioration of the diaphragm.
7. If unable to remove the diaphragm.
8. The size and shape of the vagina change and this may require a new size diaphragm. As a matter of routine, each time a pelvic examination is performed, refitting should be done. Even if the diaphragm size does not change, it is advisable to replace the diaphragm every two years or sooner.
9. Diaphragms may increase the risk of urinary tract infections especially if not properly fitted. Patients should be instructed to consult their physician if they experience any of the signs or symptoms of this type of infection, which include pain on urination, blood in the urine, elevated temperature, frequent urination, or a sensation of obstruction while urinating.
10. Persons sensitive to silicone or spermicides used with the diaphragm should discontinue use of the spermicide/diaphragm method of contraception.

CONTRAINDICATIONS

- A known hypersensitivity to silicone
- A history of toxic shock syndrome (TSS)

INSTRUCT PATIENTS

1. Not to douche or remove diaphragm for at least six hours after coitus. Diaphragm should be removed as soon after this as possible. Removal of diaphragm before six hours may increase the risk of pregnancy.
2. Not to wear diaphragm for more than 24 consecutive hours. Continuous prolonged wearing of the diaphragm may encourage the growth of certain bacteria in the vagina. It has been suggested that under certain conditions, this over growth of bacteria can lead to Toxic Shock Syndrome (TSS). Primary symptoms of TSS are sudden high fever (102 °F/39 °C or higher), nausea, vomiting, fainting or feeling faint, weakness, muscle or joint pain, sore throat and red eyes.
3. To be resized for a diaphragm if there is a gain or loss of ten or more pounds, if the patient has had a baby, an abortion or any kind of pelvic/vaginal surgery. It is also advisable to have a refitting each time patient has a routine physical examination.

DIAPHRAGM FITTING KIT

This device is intended for the sole purpose of properly fitting a diaphragm. It is intended to be inserted only for the time necessary only for the diaphragm to be sized.

DIAPHRAGM FITTING KIT INSTRUCTIONS FOR USE

To measure a patient for a diaphragm it is necessary to measure the distance from the inside of the pubic symphysis to the posterior fornix.

To do this:

1. Bring your index finger and middle finger together and insert into vagina until you reach the posterior fornix. Now lift your hand so that your index finger makes contact with the pubic arch. Mark the point directly below the inferior margin of the pubic bone with the tip of your thumb. Withdraw your fingers in this position.

2. Place one end of the fitting diaphragm on the tip of your middle finger. The opposite end should lie in front of the thumb tip. This gives you the approximate diameter for the diaphragm needed.

3. Select a fitting diaphragm closest to this diameter and insert into vagina. The proper size will fit snugly in the posterior fornix and the anterior rim snugly under the pubic arch. Feel the cervix to make sure it is completely covered by the dome of the diaphragm.

Suggestion: Try fitting the patient with one size smaller and one size larger diaphragm before your final decision as to the correct size.

Remember: The vagina lengthens and widens during arousal so it is important the diaphragm fits snugly allowing just enough room for finger to pass easily between the diaphragm and the vaginal wall.

4. Ask the patient if she can feel the diaphragm; it should not be felt once in position.
5. Ask patient to do a Valsalva maneuver, cough or move about to check for comfort and to be sure diaphragm is not displaced during normal activities.
6. If there is undue pressure or patient can feel the diaphragm, a smaller size should be tried. If the diaphragm has been displaced a larger size or a diaphragm with a different type spring should be tested.

WARNINGS AND PRECAUTIONS

The fitting Diaphragms are not intended for contraception.

CLEANING/DISINFECTION INSTRUCTIONS

After each patient use, follow the combined cleaning/disinfection instructions below for the diaphragm fitting kit.

1. CLEANING

Prepare a cleaning solution by mixing a mild soap (such as Dawn® or equivalent) with tap water using the soap manufacturer's recommended concentration. Prepare this solution in a container large enough to fully submerge the device.

2. SOAK AND SCRUB

- a. Soak the device in the container of prepared soap solution for a minimum of 5 minutes.
- b. Following the 5-minute soak period, scrub the device for a minimum of 15 seconds with a soft-bristled brush, such as a tooth brush and/or pipe brush. Scrub device below water line to prevent aerosolization of contaminants.
- c. Following scrub, inspect device for visible soil residue.

3. RINSE

- a. Remove the device from the soap solution and thoroughly rinse under flowing tap water for a minimum of 30 seconds.
- b. Allow the device to dry.

HIGH LEVEL DISINFECTION

1. Prepare a container of Revital-Ox™ RESET™ STERIS, Cidex® OPA or equivalent disinfectant in a container large enough to fully submerge the device.
2. Transfer the device to the container to soak. Ensure that the device is completely immersed in the disinfecting solution. Soak the device for at least 9 minutes in RESET or 12 minutes in Cidex OPA.
3. Rinse the device thoroughly in sterile water.
4. Place the cleaned/disinfected device in an appropriate dry storage area.

STERILIZATION

STERILIZATION PROCESS	EXPOSURE TEMPERATURE	EXPOSURE TIME	DRY TIME
Pre-Vacuum	273° F/134° C	3 minutes	30 minutes
Pre-Vacuum	270° F/132° C	4 minutes	30 minutes
Gravity	250° F/121° C	30 minutes	30 minutes

Pregnancy Probability

Contraception Method	6-month Pregnancy Probability	12-month Pregnancy Probability
Surgical Sterilization Injectable Hormones IUDs	Less than 1%	Less than 1%
Hormone Pills Vaginal Ring	1-2%	1-2%
Male Condom	7%	11% ¹
Contraceptive Diaphragm	8%	17%
Cervical Cap	11%	17%
Female Condom	13%	21% ¹

¹ These 1-year probabilities are projected since most barrier studies today are conducted as 6-month studies.

EXPLANATION OF SYMBOLS

REF	Reorder Number
LOT	Batch Code
	Consult instructions for use
	Caution
Rx Only	U.S. Federal law restricts this device to sale by or on the order of a physician.
	Manufacturer
EC REP	Authorized Representative in the European Community.
	Not made with natural rubber latex.
	Do not use if package is damaged

Dawn® is a registered trademark of Procter and Gamble Company.

Cidex® is a registered trademark of Johnson & Johnson

CooperSurgical® is a registered trademark of CooperSurgical, Inc.

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© 2019 CooperSurgical, Inc.
CooperSurgical® un marchio depositato di CooperSurgical, Inc.