# Cooper Surgical

## 47-9005

(English)

Instructions for Use  
**Adult Esophageal Balloon Catheter Set**

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**C €0086**

# 47-9005 • Adult Esophageal Balloon Catheter


## Instructions for Use (English)

### DEVICE DESCRIPTION

The set contains an 33.9 in. (86 cm) closed-end catheter with a balloon of 3.7 in. (9.5 cm) length. The catheter has depth markings to aid in positioning the balloon in the lower third of the thoracic cavity. (See Figure 1.)

A stylet with a “Y” connector is provided to allow pressure readings during placement to aid in positioning.

A pressure extension tube with male/female luer and a 3-way stopcock is provided for connection to pressure measurement instruments. Not made with natural rubber latex.

### WARNINGS

- If the catheter meets obstruction, **DO NOT FORCE THE CATHETER**. Remove it and insert it through the other naris.
- Avoid placement of the catheter in the trachea to avoid risk of airway obstruction or trauma to the respiratory system.
- Particular care should be taken if the patient has an endotracheal tube in place, as devices of this nature tend to direct the catheter into the trachea.
- Do not insert the stylet while the catheter is in the patient as there is a chance that the stylet could exit through a catheter shaft balloon hole, puncture the balloon and cause injury.
- Contents supplied sterile. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

**CAUTION:** U.S. Federal law restricts this device to sale by or on the order of a physician.

### INTENDED USE/INDICATIONS

The balloon catheter is inserted by nasal intubation into the esophagus for monitoring and measuring intrathoracic pressure changes related to respiration. It is intended for use in pulmonary function studies, sleep apnea diagnostic studies and for patients on mechanical ventilation.

### CONTRAINDICATIONS

Use of the esophageal catheter is contraindicated in patients with diseases such as esophageal ulcerations, tumors, diverticulitis, bleeding varices or in patients with sinusitis, epistaxis or recent nasopharyngeal surgery.

### PRECAUTIONS

The stylet and “Y” connector are bonded into a single unit and should be removed from the catheter before taking final pressure measurements (pressure measurements for placement can be done with the stylet in place).

### ADVERSE EVENTS

None known.


### PATIENT COUNSELING INFORMATION

Have the patient be NPO for 4 hours prior to catheter insertion. The patient should be instructed to wait for numbness to vanish before eating.

### INSTRUCTIONS FOR USE

#### Patient/Catheter Preparation

1. Have available for use an emesis basin, tissues, a protective drape, topical anesthetic, water-soluble lubricant and a 1 to 5 cc (ml) glass syringe.
2. Select a naris with the best airflow for catheter insertion.
3. If necessary apply a suitable topical anesthetic (e.g. 2-4% Lidocaine Spray) to the patient's nasal passage and throat.
4. Remove the sterile radiopaque catheter with the stylet from its protective package.
5. Remove the yellow protective sleeve from the catheter and discard.
6. Apply water-soluble lubricant to the distal tip of the catheter.
7. With the patient's head in a neutral position or flexed slightly forward, slowly insert the catheter through the naris and hypopharynx using a gentle advancing motion.
8. Avoid placement of the catheter in the trachea. Tracheal placement can be identified by patient choking or airway obstruction causing an increase in airway resistance and pressure.
9. To estimate depth in which to place catheter calculate the product of the patient's height x 0.288.
10. Advance the catheter to the calculated depth mark. (If the catheter meets obstruction, **DO NOT FORCE THE CATHETER**. Remove it and insert it through the other naris.) At this depth, the balloon will be entering the thoracic cavity.
11. Attach the extension tubing to the “Y” connector of the stylet, a syringe and an isolated physiologic transducer to the 3-way stopcock. (See Figure 2.)

**Catheter Placement**

12. Turn the 3-way stopcock open to the syringe and extension tube. Evacuate all the air from the balloon by pulling back on the syringe plunger and then allowing the plunger to return to a non-vacuum position. Use of a glass syringe avoids creating a vacuum in the balloon catheter.
13. Turn the 3-way stopcock off to the extension tube, remove the syringe and fill the syringe with 1 cc (ml) of air.
14. Attach the syringe to the 3-way stopcock; open the ay stopcock from the syringe to the extension tube. Introduce 1 cc (ml) of air into the balloon. The balloon will now be semi inflated. An incorrect amount of air in the balloon will adversely affect pressure wave performance. (see **Trouble Shooting Guide**).
15. After 1 cc (ml) of air is injected, turn the stopcock off to the syringe, and open from the extension tube to the transducer to read pressure from the catheter.

If no or a damped pressure signal is seen, the catheter may need to be advanced further into the thoracic cavity or may be kinked on itself and needs to be withdrawn.

In the absence of diaphragmatic paralysis, the pressure recorded should be negative on inspiration. A positive inspiratory pressure recording may indicate gastric placement and the catheter should be pulled back.

Pressures taken through the “Y” connector are for balloon placement only. The stylet, luer cap and “Y” connector are bonded into a single unit and must be removed from the catheter before taking pressure measurements for clinical purposes.

16. Once the catheter has been positioned properly, disconnect the extension tube from the stylet “Y” connector and remove the stylet assembly from the catheter. Excessive curvature of the catheter may cause the stylet to bind in the catheter making removal difficult. **If the stylet binds in the catheter during removal, instruct the patient to raise their head to straighten the catheter.**

**Pressure Data Acquisition**

17. After removing the stylet assembly, reattach the extension tube to the luer of the catheter and repeat steps 12 through 15.
18. Radiographic placement of the balloon is recommended to verify proper balloon placement.
19. When the catheter is properly positioned it can be secured with tape to prevent extubation or movement.
20. Take pressure measurements.  
**Note:** For tolerance specifications and degree of accuracy for the pressure readings, refer to the pressure measuring device manufacturer’s Instructions for Use (IFU).
21. Upon completion of the pressure measurements, deflate the catheter prior to removal.
22. Dispose of in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.

TROUBLE SHOOTING GUIDE	
Problem	Cause
Pressure waveform flat on top	Not enough air in balloon.
Pressure waveform flat on bottom	Too much air in balloon.
Pressure waveform all position	Too much air in balloon. Balloon may be in stomach. Transducer not zeroed to atmosphere. Patient position may need adjustment.
Pressure waveform is dampened	Too much air in balloon. Air/fluid interface in pressure line. Do not use fluid.
Pressure waveform is flat with CPAP	Added pressure from CPAP requires additional air in balloon - 1 cc (ml) more for a CPAP of 5-10 cm H <sub>2</sub> O.

**EXPLANATION OF SYMBOLS**

 Reorder Number

 Batch Code

 Use by Date

 Sterilized using Irradiation

 U.S. Federal law restricts this device to sale by or on the order of a physician.

 Do not re-use

 Caution

 Consult instructions for use

 Do not use if package is damaged

 Authorized Representative in the European Community.

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© 2018 CooperSurgical, Inc.