44-8000 • NEO-prep™ Neonatal Skin Prep Tabs • Small

44-8001 • NEO-prep™ Neonatal Skin Prep Tabs • Large

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DEVICE DESCRIPTION

Hypoallergenic skin barrier tabs pre-cut to small and large sizes for use during nasal cannulae procedures. Non-sterile. Not made with natural rubber latex.

WARNINGS

- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.
- NEO-prep should not be used if infant has skin irritation.
- Discontinue use if NEO-prep interferes with the administration of oxygen or other therapies.

CAUTION

U.S. Federal law restricts this device to sale by or on the order of a physician.

INTENDED USE / INDICATIONS

To act as a skin barrier between nasal cannulae and the skin of the upper lip and nasal septum to minimize skin irritation or damage.

CONTRAINDICATIONS

Do not use on irritated or compromised skin.

PRECAUTIONS

For proper adhesion insure that the skin is clean, dry and free of grease or oil prior to application. Not intended to be used to secure nasal cannulae. Use only as a protective skin barrier.

DIRECTIONS FOR USE

1. Remove the NEO-prep™ tab from the backing strip.
2. Position the NEO-prep tab on the upper lip with the narrow portion aligned on the septum of the nose. See figure.

Note: The NEO-prep tab may be trimmed according to the patient's anatomy.

EXPLANATION OF SYMBOLS

REF Reorder Number

Do Not Reuse

Consult instructions for use

Rx Only U.S. Federal law restricts this device to sale by or on the order of a physician.

LOT Batch Code

CE Product conforms to the Medical Device Directive 93/42/EEC

Use By

Caution

Manufacturer

EC REP Authorized Representative in the European Community.

Not made with natural rubber latex.

NEO-prep™ is a trademark of CooperSurgical, Inc.
CooperSurgical is a registered trademark of CooperSurgical, Inc.

© 2014 CooperSurgical, Inc..

95 Corporate Drive • Trumbull, CT 06611 USA
Phone: (800) 243-2974 • Fax: (800) 262-0105
www.coopersurgical.com
003223LS-DFU • Rev. A • 2/14

International
Phone: (203) 601-9818
Fax: (203) 601-4747

EC REP EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands

Made in the USA

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Hypoallergenic skin barrier tabs pre-cut to small and large sizes for use during nasal cannulae procedures. Non-sterile. Besteht nicht aus Latex aus Naturkautschuk.

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NEO-prep™ è un marchio commerciale di CooperSurgical, Inc.
CooperSurgical è un marchio registrato di CooperSurgical, Inc.