ABBI™ • ABBI™-KIT

Air Bubble Based Infuser

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C€0086

DEVICE DESCRIPTION

The ABBI™ (Air Bubble Based Infuser) device is a sterile, single-use, two-chamber syringe that connects to an intrauterine infusion catheter for uterine cavity and/or fallopian tube evaluation. The ABBI device comes with a male luer for connection to an intrauterine infusion catheter, a saline input port, a 10 mL air chamber and a 10 mL saline chamber, a selection dial (to switch between saline only or a saline and air bubble mixture) and a two-position comfort grip thumb plunger.


Figure 1

CAUTION

U.S. Federal law restricts this device to sale by or on the order of a physician.

WARNINGS/PRECAUTIONS

- For single use only. Do not reuse, reprocess, or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure that, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient. After use, dispose of the device and packaging in accordance with hospital, administrative, and/or local government policy.
- Do not use oil-based contrast media
- To minimize the risk of air embolism, during Sono HSG, do not infuse more than 60 cc of air (six (6) filled air-saline syringe volumes).
- Make sure to prime device before SIS or Sono HSG procedures.
- ABBI has only been validated for use with CooperSurgical, Inc.'s intrauterine infusion catheters part numbers 61-4005 and 61-4005S. Use with any other intrauterine infusion catheter has not been validated.

INDICATIONS FOR USE

The ABBI (Air Bubble Based Infuser) is indicated for the following:

- Instillation of a consistent alternating pattern of saline and air as a continuous stream of contrast media into the uterus and fallopian tubes to be used in conjunction with an intrauterine catheter for performance of sono-hysterosalpingogram (SonoHSG).
- Instillation of saline as a continuous stream of contrast media into the uterus to be used in conjunction with an intrauterine catheter for performance of saline infusion sonography (SIS).

CONTRAINDICATIONS

- Suspected pregnancy, profuse bleeding, possible pelvic inflammatory disease (PID) or vaginal infection, cervical stenosis, sexually transmitted disease, and/or uterine malignancies.
- Sono HSG is contraindicated in women who are pregnant or who have been pregnant in the previous 6 weeks (including miscarriage). These women may be at increased risk for air embolism.

INSTRUCTIONS FOR USE

Note: It is normal to hear a slight whistling sound when depressing the thumb plunger to expel the contents of the air and saline chambers of the ABBI device. This whistling sound may be heard during use, but it poses no harm to the patient.

Note: ABBI provides the option of attaching a sterile saline bag or a sterile saline-filled syringe (both not included) to fill the saline chamber using the saline input port.

A. General Examination

- Perform a pelvic examination.
- Perform a preliminary transvaginal ultrasound to exclude any uterine or adnexal pathology and to localize the ovaries and the interstitial portion of the fallopian tubes.

B. ABBI Preparation

- Connect the sterile saline bag/syringe to the ABBI saline port.
- Pull back the thumb plunger until the saline chamber is filled completely with sterile saline. While the saline chamber is filling, ABBI is also drawing air into the air chamber through an internal 0.2 micron air filter.
- Prime the device before use by holding the ABBI with the luer up so that the air bubbles rise to the top of the syringe. Then depress the device until the air is evacuated from the saline chamber and a small amount of contrast media is released from the device.

Figure 2
Saline Input Port

C. Position the Intrauterine Infusion Catheter

See intrauterine infusion catheter Instructions for Use for proper catheter placement.

D. Injection of Saline Only for the SIS Procedure (Note: SIS is an optional procedure.)

It may be necessary to prime ABBI so that air isn't introduced into the uterine cavity during the SIS procedure.

- Attach the ABBI male Luer to the intrauterine catheter hub. When using a balloon type intrauterine infusion catheter, take special care to make sure that ABBI is connected to the infusion channel of the catheter, and not the balloon filling channel, to avoid mistakenly overinflating or rupturing the balloon. See balloon catheter's Instructions for Use for more information.

Figure 3
Male luer for
Connection
to Intrauterine
Infusion
Catheter

2. Make sure the ABBI dial is set to saline. Perform an SIS procedure according to your practice protocol using ABBI as your syringe.
3. Repeat step 2 if necessary. To refill ABBI with sterile saline, pull back the thumb plunger when ABBI is connected to a sterile saline bag or syringe filled with sterile saline.


Figure 4 Saline

E. Injection of Saline/Air Bubble Media for the Sono HSG/HyCoSy Procedure

1. Make sure the ABBI dial is set to Saline/Air.
2. When performing an air bubble/saline media based Sono HSG/HyCoSy procedure, slowly and consistently inject the saline/air bubble media into the uterine cavity while evaluating the flow of the media under ultrasound imaging using a transvaginal ultrasound probe in the transverse orientation. The air bubbles are highly echogenic and assist in confirming tubal patency by allowing observation of the media flow into and through the cornual fundal junctions, into and through the fallopian tubes, and out of the fimbria, or by confirming tubal occlusion by allowing observation of the media pool around either of the tubal ostia or within either of the fallopian tubes.
3. Repeat step 2 if necessary by refilling ABBI with sterile saline and pull back the thumb plunger.
4. When procedure is completed, remove ABBI from the intrauterine catheter.

Figure 5 Saline/Air


Figure 6

F. Remove Intrauterine Infusion Catheter per the catheter's Instructions for Use.

G. Disposal

Dispose of all products in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

H. Patient Monitoring

At the end of the procedure, it is recommended that the patient be monitored for about 15 minutes in the supine position in order to prevent the onset of vasovagal reactions or scapular pain, resulting from the possible irritation of the phrenic nerve.

EXPLANATION OF SYMBOLS

Reorder number

Batch code

Use-by date

Caution

Do not re-sterilize

Consult Instructions for Use

Do not use if package is damaged.

Do not re-use

Sterilized using ethylene oxide

R_x Only

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

Not made with natural rubber latex.

Authorized Representative in the European Community

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© 2018 CooperSurgical, Inc.
CooperSurgical, Inc. | 95 Corporate Drive | Trumbull, CT 06611 USA