RING Folding Pessaries

For uterine prolapse (Stage I and Stage II)

AVAILABLE MODELS

RING without SUPPORT	RING with SUPPORT
REF MXPER (Pessary Only)	REF MXPRS (Pessary Only)
MXPER (Pessary, TRIMO-SAN™, Jel Jector)	MXPRS (Pessary, TRIMO-SAN™, Jel Jector)
RING with KNOB without Support	RING with KNOB with SUPPORT
REF MXPKR (Pessary Only)	REF MXPRSK (Pessary Only)
MXPKR (Pessary, TRIMO-SAN™, Jel Jector)	MXPRSK (Pessary, TRIMO-SAN™, Jel Jector)

Single-Patient Use Rx Only

Wash pessary with mild soap and thoroughly rinse prior to initial use.

DESCRIPTION

The Milex® RING Pessary is a medical device made from silicone that is inserted into the vagina to function as a supportive structure of the uterus, bladder and/or rectum.

WARNINGS

- Sizes 11 through 13 contain a wire coil - remove before X-ray or MRI.
- Do not leave pessaries in place for long periods of time, as serious complications can occur which may require surgical intervention.
- Do not use these pessaries on a patient with a known silicone allergy.

CooperSurgical

38588-CR • Rev. A • 04/21

• Chemicals in various vaginal preparations can interact with the pessary material, resulting in discoloration or deterioration of the pessary. TRIMO-SAN™ does not interact with the pessary material. Use only TRIMO-SAN. Other materials have not been tested for compatibility.

INDICATIONS FOR USE

- The RING without SUPPORT pessary is indicated for support in Stage I and Stage II prolapse.
- The RING with SUPPORT pessary is indicated for support in Stage I and Stage II prolapse complicated by a mild cystocele.
- The RING with KNOB pessary, without SUPPORT is indicated for stress urinary incontinence complicated by Stage I or Stage II uterine prolapse.
- The RING with KNOB with SUPPORT pessary is indicated for stress urinary incontinence complicated by Stage I or Stage II uterine prolapse and mild cystocele.

CONTRAINDICATIONS

- The presence of pelvic infections or lacerations
- A noncompliant patient
- Endometriosis
- Pregnant patient

RECOMMENDED CLEANING INSTRUCTIONS*

- Prepare a cleaning solution by mixing a mild soap with tap water using the soap manufacturer's recommended concentration. Prepare this solution in a container large enough to fully submerge the device.
- Soak and Scrub
 - Soak the device in the container of prepared soap solution for a minimum of 5 minutes.
 - Following the 5-minute soak period, scrub the device for a minimum of 15 seconds with a soft-bristled brush, such as a tooth brush and/or pipe brush. Scrub device below water line to prevent aerosolization of contaminants.
 - Following scrub, inspect device for visible soil residue.
- Rinse
 - Remove the device from the soap solution and thoroughly rinse under flowing tap water for a minimum of 30 seconds.
 - Allow the device to dry.

TRIMO-SAN is available at the prescription counter of your pharmacy. A prescription is required.

REPORT ANY OF THE FOLLOWING SYMPTOMS TO YOUR PHYSICIAN

- Any difficulty in urinating
- Any discomfort
- Any changes in the color or consistency of vaginal discharge
- Any increase in the amount of vaginal discharge or vaginal bleeding
- Any foul odor associated with vaginal discharge
- Vaginal itching
- If the pessary falls out

For medical emergencies and for all medically related advice, consult your healthcare professional.

FOR THE PHYSICIAN / HEALTHCARE PROFESSIONAL

Review these instructions with the patient to establish use regimen.

INSTRUCTIONS

- Wear dry gloves. When necessary, lubricate only the entering end of the pessary with TRIMO-SAN. Hold as illustrated in Figure 1.
- The pessary is folded along the axis of the bigger outer holes (or the inner notches for the RING without SUPPORT) by bringing


Figure 1


Figure 2

Figure 3

the small round holes together. The KNOB will be at the top of the arch (see Figure 2). The arch formed points downward as shown with the RING pessary.

- Direct the pessary past the cervix into the posterior fornix. Allow the pessary to open again into the ring shape after passing the introitus.
- The index finger is inserted deep into the vagina to turn the pessary approximately 90° so that the KNOB is resting behind the symphysis pubis (see Figure 3). The RING pessaries in this position cannot be folded and pushed out.

- Ask the patient to sit, and bear down slightly. If there is no leakage, and the patient is comfortable with the pessary in position, have her empty her bladder. A properly fitted pessary takes up slack in redundant tissue, holding the uterus higher in the vagina.
- If the patient can void without difficulty, and the pessary remains in position upon re-examination, and the patient is comfortable with the pessary in place, this is a good indication that the correct size may have been selected. Patient experience may vary.

Note: If the patient is unable to urinate with the pessary in position, remove it and fit her with the next smaller size. Repeat as necessary.

- Examine the patient while she is in the standing position to ensure the pessary has not shifted position. The patient should not feel the pessary once it is in position. The pessary should not be too loose as it may turn or be expelled and it should not be too tight as it may cause discomfort.
- The healthcare professional should be able to sweep one finger between the pessary and the vaginal wall. If there is not enough space to do this, the next smaller size should be tried. If excessive space exists, the pessary will not be effective and may rotate or even be expelled.

- It is sometimes necessary to re-fit the patient with a different size or type of pessary after a period of time. Do not assume that a replacement will always be the same size as the previous one. Check the fitting to ensure continued patient comfort and relief of symptoms. The useful life of a pessary is limited. Examine frequently for signs of deterioration (such as cracks or breaks in silicone outer surface). A pessary should be replaced if damaged.

- Ulcerations and erosions frequently occur in cases of complete prolapse due to irritation of the exteriorized cervix or vaginal wall. Whenever possible, reducing the mass and treating the irritation are primary steps

before using a pessary. Prolapse reduction may resolve cervical vaginal irritation. Verification of cervical cytology (pap) and or biopsy as clinically indicated.

- During each visit, the vagina should be carefully inspected for evidence of pressure or allergic reaction. The patient should be questioned concerning douching, discharge, disturbance of bowel function or urination. It may be necessary to fit another size or an entirely different type of pessary.
- At the physician's discretion, the patient can be instructed in the proper removal, cleaning and reinsertion techniques for her own pessary. This process can be performed nightly or even weekly by the patient under ideal circumstances.

TO REMOVE

- Use one finger to depress the perineum. Turn the pessary until the notches face the introitus. Fold the pessary and gently ease it out.

RECOMMENDED FOLLOW-UP

- Have patient return within 24 hours for first examination.
- Have patient return for second examination within 3 days.
- Have patient return for examination every few months.

Note: The above schedule of follow-up examinations may be altered to fit the needs of the individual patient at the discretion of the healthcare provider.

EXPLANATION OF SYMBOLS

REF Reorder number

LOT Batch code

Date of manufacture

Use-by date

Consult instructions for use

Not made with natural rubber latex

Rx Only CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician

Manufacturer

TRIMO-SAN™ is a trademark of CooperSurgical, Inc. Milex® and CooperSurgical are registered trademarks of CooperSurgical, Inc.

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© 2021 CooperSurgical, Inc.