# Sacrocolpopexy/Sacrocervicopexy Tips

## Instructions For Use (English)

CE 2797

Sacrocolpopexy Tip

Sacrocervicopexy Tip

### DEVICE DESCRIPTION / INDICATIONS FOR USE

The Sacrocolpopexy Tips (Sacro Tips) are used for laparoscopic vaginal prolapse repair procedures where the cervix has been removed.

The Sacrocervicopexy Tip (Sacro Tip) is used for laparoscopic vaginal prolapse repair procedures where the cervix is intact.

The Sacro Tips are intended to be connected to any of the RUMI® or Advincula Arch™ Uterine Manipulator Handles (Manipulator Handle) and utilized for any laparoscopic or open procedure where elevating vaginal tissue into the operative field would be beneficial. The single-use Sacro Tips are made from hard plastic that is intended to provide a solid base, or “thimble-like” surface to suture against.

Not made with natural rubber latex.

### WARNING

Do not use if the pouch seal is damaged or broken.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

### CAUTION


U.S. Federal law restricts this device to sale by or on the order of a physician.

### CONTRAINDICATIONS

The Sacro Tips should not be used when a patient has a non-patent vagina or vaginal septum.

### INSTRUCTIONS FOR USE

Prior to use, inspect the Sacro Tip.

1. Grasp the Sacro Tip with one hand and the shaft of the Manipulator Handle with the other. Line up and insert the Rod (on the end of the Sacro Tip) into the hole on the end of the Snap Drum.

2. With a clockwise twisting motion of the Tip, push the Tip completely over the Snap Drum. The Tip will snap into place around the Snap Drum.

3. Visually inspect to ensure that the Tip is seated completely over and around the Snap Drum.

4. A slight tug on the Tip will verify that it is on securely.

Dispose of in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.

## EXPLANATION OF SYMBOLS

**[REF]** Reorder number

**[LOT]** Batch code

Use-by date

Consult instructions for use

Do not re-use

Do not resterilize

Do not use if package is damaged

Caution

Not made with natural rubber latex

**[STERILE] [EO]** Sterilized using ethylene oxide

**[EC] [REP]** Authorized Representative in the European Community

**RxOnly** U.S. Federal Law restricts this device to sale by or on the order of a physician

Manufacturer

RUMI® is a registered trademark of CooperSurgical, Inc.  
CooperSurgical is a registered trademark of CooperSurgical, Inc.  
Advincula Arch™ is a trademark of CooperSurgical, Inc.

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© 2021 CooperSurgical, Inc.
RUMI® è un marchio registrato di CooperSurgical, Inc.
CooperSurgical è un marchio registrato di CooperSurgical, Inc.
Advincula Arch™ è un marchio di CooperSurgical, Inc.