Cervical Sizer

for use with the RUMI® Koh-Efficient® or Koh Colpotomizer™ System

KC-SIZER

Instructions for Use (English)

NON-STERILE • STERILIZE BEFORE USE

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Cervical Sizer (Sizer) is a medical device used to determine the diameter of a cervix and for correctly sizing a colpotomy cup.

Not made with natural rubber latex.

INDICATIONS FOR USE

The Sizer is designed for use in conjunction with the Koh-Efficient® or Koh Colpotomizer™ System, where the cervix is present.

CONTRAINDICATIONS

The absence of a cervix.

PRECAUTIONS

- Do not use if the measuring platforms on the distal end of the device are not reading properly or are damaged.
- Sterilize before use.

INSTRUCTIONS FOR USE

Calibration test prior to insertion: To confirm that the Sizer is within its calibrated parameter, set the Scale on the Sizer to the 4 cm position and hold it up to a 4 cm Koh Cup™ or metric scale. If the device is not calibrated correctly, do not use.

1. Under direct visualization, place the Sizer around the cervix at the fornix in the 4.0 position (Figure 1).
2. Close until you approximate the cervix (Figure 2). **Caution: Compression of the cervix may cause a false small measurement.**
3. Take note of the cup size indicated on the Scale of the Sizer (Figure 3).
4. Once the correct cup size is determined, remove the Sizer.

REPROCESSING INSTRUCTIONS

1. INTRODUCTION

This section is intended to provide detailed instructions for effectively processing the Sizer. All reusable instruments should be thoroughly cleaned and sterilized to prepare them for use.

CooperSurgical, Inc. has validated the processes provided in these instructions to be capable of effective instrument processing. Equipment, operators, cleaning agents and procedures all contribute to the efficacy of processing. Healthcare facilities should ensure that selected processing steps are safe and effective within their systems.

Alternative methods of processing these instruments outside the instructions described in this document may be suitable for reprocessing; however, these have not been evaluated by CooperSurgical, Inc. Operators and healthcare facilities which choose to perform processes outside the instructions described in this document must validate these processes before use. In the event that national or regional government requirements conflict with the recommendations provided here, these shall override the recommendations of CooperSurgical, Inc., but must be validated before use.

Reusable instruments must be rinsed thoroughly after use and after cleaning to ensure removal of gross soil and residual cleaners and solvents. Gross soil is damaging to instrument surfaces and inhibits thorough cleaning and subsequent sterilization. Residual cleaners and solvents can impede further processing and may affect the surface finish of the instrument over time.

2. CLEANING



There are two methods for cleaning that have been validated by CooperSurgical, Inc. Any healthcare facility should be able to perform the manual cleaning process.

For those facilities that use automated washer-disinfectors, an automated method has also been validated. Personnel are cautioned to use personal protective apparel due to the unknown and potentially hazardous nature of biological fluids and soils present. Specific detergents or cleaning solutions vary in requirements for concentrations and temperatures required for optimal cleaning performance. Use temperatures in this document as guidance, to be superseded by the detergent/cleaner manufacturer's directions for use. Any processes outside the recommendations in this document must be validated before use.

The quality of water used to prepare solutions and rinse reusable instruments can affect the efficacy of the process and the instrument reuse life. CooperSurgical, Inc. recommends the use of freshly prepared purified water or sterile water for dilutions of solution and rinsing of instruments. Non-purified water can add mineral deposits and recontaminate instruments with microorganisms. Mineral deposits can impede sterilization and affect the condition of the instrument, resulting in staining, corrosion, and/or other damage.

Equipment and Materials Required for Manual Cleaning

- Personal protective apparel/equipment as recommended by cleaning solution supplier.
- Cleaning bath, sink, or other vessel large enough to accommodate full immersion of instruments.
- Freshly prepared cleaning solution intended for manual cleaning (use enzymatic, neutral, or alkaline detergents).
- Soft-bristled brushes and sponges.
- Clean, low-lint or lint-free cloth.

Equipment and Materials Required for Automated Cleaning

- Personal protective apparel/equipment as recommended by cleaning solution supplier.
- Legally marketed (FDA-cleared and/or CE-marked) medical device washer or washer-disinfector.
- Cleaning solution intended for automated cleaning (use enzymatic, neutral, or alkaline detergents).
- Racks and accessories to hold and support instruments during processing.

Point of Use

During and immediately after the clinical procedure, instruments should be treated to reduce and remove gross soil materials. Instruments should be wiped with low-lint or lint-free absorptive materials, such as low-lint or lint-free wipes to remove solid tissues and the majority of bodily fluids. Immediately after use, all disposable single-use accessories should be removed and discarded into appropriate biohazardous waste containers. Instruments should be placed into a transport container that will protect them from damage during transport and should be kept moist.

Transport to Processing Area

Transport the instruments to the processing area as soon as practical. Minimize holding time before removing organic debris. Ensure that lighter instruments are not mixed with heavier instruments to avoid damage to the instruments due to mechanical abrasion. Seal or enclose the transport container to prevent fluid loss and the potential for cross-contamination of other areas.

Preparation

- Remove all single-use accessories and discard in an appropriate waste container designed for biohazardous materials if this has not been performed at point of use.
- There are no reusable accessories for the Sizer.

Cleaning

Cleaning should take place as soon after use as practical. Soils should not be allowed to dry onto the instrument surfaces. The manual process and the automated process were validated independently and are considered equivalently effective when performed according to the instructions in this document. Do not use metal-bristled brushes on the instruments.

Manual Cleaning Processes

1. Prepare a neutral or alkaline enzymatic detergent specifically designed for manual cleaning of medical instruments. Prepare a fresh solution of detergent at the detergent manufacturer's recommended concentration using water at the detergent manufacturer's recommended temperature.
2. Ensure that the instruments are fully immersed in the prepared solution. Soak the instrument for a minimum of ten (10) minutes at room temperature (68-77°F/20-25°C) or longer if indicated on the detergent manufacturer's instructions for use.
3. After the soak time, scrub the handle with a soft sponge for large surfaces. Using a soft-bristled brush, scrub the regions. Scrub regions with crevices. Rinse the instrument in warm tap water (100-120°F/ 38-49°C) for one (1) minute. Ensure that all crevices are flushed.
4. Place the instrument into a prepared ultra-sonicator containing a neutral or alkaline detergent specifically designed for manual cleaning of medical instruments and labeled for use with an ultra-sonicator. Ultra-sonicate the instrument in the detergent solution for ten (10) minutes.
5. Rinse the instrument with purified or sterile water for one (1) minute or until all visible signs of detergent residue are removed, whichever is longer. Ensure that all crevices, lumens and mated surfaces are flushed. If there are any signs of residual soil on the instrument, repeat the cleaning process.
6. Dry the instrument with a low-lint or lint-free towel or wipe.

Automated Cleaning Process

1. Ensure that all gross soil has been wiped or rinsed from the surfaces of the instruments.
2. Use only a legally marketed (FDA-cleared and/or CE-marked) medical device washer or medical device washer-disinfector. Industrial use and household use dishwashers are not acceptable for the cleaning of medical devices.
3. Place the instruments loosely within a rack designed for the medical device washer or medical device washer-disinfector, using a hold-down screen if indicated by the washer/washer-disinfector manufacturer.
4. Load the rack with the instruments into the washer/washer-disinfector and close the door.
5. Select a cycle intended for instruments that has the following parameters at a minimum:

PARAMETER	TOLERANCE
Detergent	Neutral or Alkaline medical device detergent
Pre-wash rinse	Minimum of one
	Minimum of 00:15 (mm:ss)
	Cold tap water or better quality
Wash	Minimum of one repetition
	Minimum of 02:00 (mm:ss)
	Warm tap or purified water (minimum of 140°F/60°C)
	Minimum concentration per detergent manufacturer's IFU
Post-wash rinse	Minimum of one
	Minimum of 00:15 (mm:ss)
	Warm tap or purified water (minimum of 110°F/43.3°C)
Thermal rinse*	Minimum of one
	Minimum of 01:00 (mm:ss)
	Purified water (minimum of 180°F/82.2°C)
Drying**	Minimum of default cycle setting

* Only heat-based medical device washer-disinfectors will include a thermal rinse. Medical device washers that do not include a thermal disinfection phase should be set to have a minimum of a 01:00 (mm:ss) rinse following the post-wash rinse using purified water.

** Drying time will be variable dependent on the size and composition of the load, altitude, environmental conditions, and air temperature and source characteristics.

NOTE: automated wash processes that include additional phases or longer phase durations are anticipated to deliver equivalent or better cleaning efficacy as the minimal validated parameters. Additional phases and longer phase durations may be added if they are within the healthcare facility's normal procedures.

6. Start the cycle and allow it to run through conclusion.
7. Medical washer-disinfector default cycles often use high heat during the final stages of the processing cycle. Use caution and heat-protective gloves when opening and removing the wash rack at the end of the cycle.
8. Examine the instruments for signs of residual soil. If there are any signs of residual soil on the instrument, repeat the cleaning process.
9. If the instruments are not thoroughly dry following the automated process, dry the instruments with low-lint or lint-free towels or wipes.

3. INSPECTION

All instruments must be inspected for signs of damage, wear or residual soil each time they are reprocessed. If residual soil is observed during inspection, the cleaning process must be repeated for the instrument. Remaining soils can impede effective sterilization by shielding remaining microorganisms and can lead to tissue reactions in subsequent patients.

Fully articulate the moving components of the handle using the trigger and grip. The motion should be smooth. Jerking and catching within the range of motion may be signs of damage or retained material between the mated surfaces. Ensure that the trigger moves easily when depressed and that it automatically drops down when released. The un-actuated position should be locked into the articulating dial end of the shaft. When the trigger is engaged in the articulating end of the shaft, the instrument should not shift or rotate. Damage to the trigger locking mechanism may allow slippage during a procedure. Corrosion may weaken the materials leading to a break under normal load. Signs of damage or poor motion quality may indicate a need to replace the instrument or have it repaired.

Surface Finish

Inspect the surface finish on the instrument; the surface should be smooth. The shafts of the instrument have a polished finish that should be bright silver-gray in color. The finger-grip region is gold-plated with a polished finish that should be a bright yellow-gold color. The interlocking sizing plate has a matte finish that is gray in color with etched markings. Additionally there are etched markings near the hinge joint in the mid-shaft region. Pay attention to the etched markings on the surface as signs of corrosion will typically appear in these areas first.

Hinge Joint

Fully open the instrument and inspect the surface of the hinge joint for signs of trapped soil or corrosion. Look at the surfaces for signs of excessive wear between the mating surfaces. Holding the instrument up in front of a light source may allow easier visualization in the joint. Open and close the instrument, checking the feel. The motion should be smooth and there should not be any catching, jerking, or roughness until the interlocking sizing plate and opposing sizing lock meet. Rough motion may indicate damage or trapped soil within the hinge joint.

4. WRAPPING

Prior to terminal sterilization, the instruments must be packaged to preserve sterility after processing during storage. Legally marketed sterilization packaging (e.g. CSR wrap) that is compatible with and labeled for use with the chosen sterilization method should be used. CooperSurgical, Inc. has not evaluated the sterilization efficacy of the processes recommended in this document using containment systems (e.g. procedure/case trays, cassettes, rigid reusable sterilization containers), therefore, users who choose to use a containment system in the sterilization of these instruments should verify that the process is effective under these conditions. Containment systems can impede the ingress and egress of sterilant (moist heat). CooperSurgical, Inc. recommends individually wrapping the instruments prior to sterilization. The instruments should be wrapped using one of the common closure techniques, using properly sized material and secured with tape or adhesive designed and labeled for use with the sterilization process chosen. Instruments should be labeled with identifying information and sterilization date to facilitate selection for use after processing.

5. STERILIZATION - MOIST HEAT (STEAM)

Individually packaged instruments using legally marketed sterile barriers with labeling consistent with the following parameters may be sterilized prior to use and between subsequent uses. Legally marketed (FDA-cleared and/or CE-marked) moist heat (steam) sterilizers should be used. The following parameters have been validated for use with the instruments within the scope of this document:

STERILIZATION PROCESS	EXPOSURE TEMPERATURE	EXPOSURE TIME	DRY TIME
Gravity Displacement	250°F/121°C	30 minutes	30 minutes
Pre-vacuum	270°F/132°C	4 minutes	30 minutes
Pre-vacuum	273°F/134°C	3 minutes	30 minutes

6. STORAGE

Individually packaged instruments should be cooled to a safe handling temperature following moist heat sterilization prior to final storage. Individually packaged instruments should be transported to a clean, temperature and humidity controlled, sterile supply storage area. If instrument packaging is compromised during storage (becomes wet, torn, or seals are open), the instrument must be reprocessed through cleaning, inspection, and sterilization prior to use. Do not stack instruments or allow other product or instruments to be stacked on top of them. Protection from dust and moisture is recommended. Ensure that traceability and identity is maintained throughout storage and up to subsequent use.

7. VALIDATION OF PROCESSES

CooperSurgical, Inc. performed validations for the cleaning of the instruments. The following detergents were used in the validations demonstrating efficacy: Klenzyme™, Manu-Klenz®, ENZOL®, Neutrad®, Prolystica® Ultra-Concentrate Neutral Detergent and Prolystica® Ultra-Concentrate Alkaline Detergent.

CooperSurgical, Inc. does not endorse the use of the identified products in lieu of other similar products designed for use with medical devices, however, this information may provide a basis for comparing the formulations to select an appropriate alternate in facilities or locations where these detergents are not readily available.

The parameters listed in the automated cleaning section as the minimum parameters for automated cleaning are identical to the parameters used in the validation, with the exception of the dry time. Post-cleaning dry time was eliminated for the cleaning process validation to facilitate detection of soil residuals, if present, in order to reduce the potential for a false negative result.

WARRANTY

The Cervical Sizer is warranted by CooperSurgical, Inc. for a period of 1 year from date of shipment. This warranty applies only to the original purchaser and only against defects in workmanship or materials, which under normal use render the instrument inoperable.

The original purchaser shall prepay the shipping costs for returning a CooperSurgical, Inc. instrument for warranty service. CooperSurgical, Inc. will, at its option and without charge, either repair or replace any instrument, which CooperSurgical, Inc. determines to be defective in material or workmanship. Instrument damage caused by misuse or accident, shall void this warranty.

CooperSurgical, Inc. disclaims any liability for special, incidental, consequential, punitive or exemplary damages arising out of the use of any CooperSurgical, Inc. instrument. CooperSurgical, Inc.'s liability in all events is limited to, and shall not exceed, the purchase price paid. Except as expressly provided in this limited warranty section, CooperSurgical, Inc. makes no representations or warranties, expressed or implied, as to the instruments, including warranties of merchantability and fitness for a particular purpose, all of which are hereby expressly disclaimed and excluded. No warranty or affirmation of fact, express or implied, other than as set forth in this limited warranty, is made or authorized by CooperSurgical, Inc.

EXPLANATION OF SYMBOLS

Reorder number

Batch code

R_xOnly

U.S. Federal law restricts this device to sale by or on the order of a physician.

Non-Sterile

Consult instructions for use

Authorized Representative in the European Community.

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RUMI® and Koh-Efficient® are registered trademarks of CooperSurgical, Inc.
Koh Colpotomizer™ and Koh Cup™ are trademarks of CooperSurgical, Inc.
Klenzyme®, Manu-Klenz® and Prolystica® are trademarks of STERIS Corporation.
ENZOL® is a registered trademark of Johnson and Johnson.
Neutrad® is a registered trademark of Decon Labs, Inc.

© 2018 CooperSurgical, Inc. All rights reserved.

CooperSurgical

 95 Corporate Drive
Trumbull, CT 06611 USA
Phone: (800) 243-2974
Fax: (800) 262-0105

KCSIZER-IFU-B • Rev. A • 04/18

International:
Phone: +1 (203) 601-9818
Fax: +1 (203) 601-4747
www.coopersurgical.com

 EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands
Made in Germany