Cytobrush® Plus GT Scored

R_x Only ②

- cs Návod na použití
- da Brugsanvisning
- de Gebrauchsanweisung
- el Οδηγίες Χρήσης
- en Instructions for Use
- es Instrucciones de uso
- et Kasutusjuhend
- fr Mode d'emploi
- fi Käyttöohjeet
- hu Használati utasítás
- it Istruzioni per l'uso
- lt Vartojimo instrukcija
- lv Lietošanas pamācība
- nl Gebruiksaanwijzing
- pl Instrukcja użytkowania
- pt Instruções de utilização
- sk Návod na použitie
- sl Navodila za uporabo
- sv Bruksanvisning

CooperSurgical®

Cytobrush® Plus GT Scored

Nesterilní

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Cytobrush® Plus GT Scored

Non-sterile

Cytobrush® Plus GT Scored is a single use device.
Discard after use.

INDICATIONS FOR USE

Cytobrush® Plus GT Scored is used to collect cell samples from the endocervix for liquid cytology and microbiology. The device has a prescored shaft and should **NOT** be used for pap smear due to risk of unintentional shaft breakage.

CONTRAINDICATIONS

Do not use Cytobrush® Plus GT Scored after the first ten (10) weeks of gestation in pregnant patients.
Do not use Cytobrush® Plus GT Scored for endometrial sampling.

INTENDED USERS

The intended users of the Cytobrush® Plus GT Scored are healthcare professionals and other trained medical personnel in a professional environment such as a hospital, clinic, or doctor's office.

PATIENT POPULATION

The device can be used by any patient women between the ages of 21 years old up to 65 years.

WARNINGS

Never allow Cytobrush® Plus GT Scored to enter too far into the endocervix. A few bristles of the brush should at all times be seen in the lower part of the cervix.

CAUTION: For single use only. Do not reuse, reprocess, or re-sterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

ADVERSE REACTIONS

Slight bleeding may occur during sampling and is acceptable. Check that bleeding has ceased before end of gynaecological examination.

PLEASE NOTE: If a serious incident is suspected from using the Cytobrush® Plus GT Scored, report the details of the incident to CooperSurgical via phone number +1 203-601-5200 Ext 3100 or by email at ProductSurveillance@coopersurgical.com and to the local Health Authority in your country. A serious incident may have caused or contributed to a death, a delay in a procedure which resulted in death or serious injury, or a malfunction that could have caused an adverse event.

SAMPLING INSTRUCTIONS

1. Gently insert Cytobrush® Plus GT Scored into the endocervical canal. A few bristles of the brush should be seen outside the cervix.
2. Rotate the brush slowly one (1) turn.
3. Withdraw Cytobrush® Plus GT Scored.

PREPARATION OF SAMPLE

Consult your laboratory for detailed instructions regarding specific tests.

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CooperSurgical Distribution B.V.
CooperSurgical, Inc.