Fischer Cone Biopsy Excisor®

**Sterile
Single-Use Only**

Instructions for Use (English) Page 1

Models - 900-150 • 900-151 • 900-152 • 900-154 • 900-155 • 900-157 • 900-158

Sterile, Single-Use Only

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
Sterile, Single-Use Only


DEVICE DESCRIPTION

The Fischer Cone Biopsy Excisor® removes a cone-shaped specimen as an alternate for Cold Knife Conization, Lazer Conization or Top Hat/Cowboy Hat sampling used with the LEEP procedure. The model numbers, sizes, wattages and dimensions are shown below:

900-157* (Small Wide Angle): Recommended wattage setting 40-45

900-158* (Small Wide Angle Extended): Recommended wattage setting 40-45

Dimensions

Model Number	Size	A mm / in	B mm / in	C (wire length) mm / in
900-150	Small	9 / 0.349	12 / 0.472	15 / 0.590
900-151	Medium	12 / 0.468	15 / 0.590	19 / 0.756
900-152	Large	15 / 0.584	18 / 0.708	23 / 0.921
900-154	Medium Extended	12 / 0.468	20 / 0.787	23 / 0.918
900-155	Large Shallow	15 / 0.584	8 / 0.315	17 / 0.669
900-157	Small Wide Angle	9 / 0.349	12 / 0.472	17.9 / 0.705
900-158	Small Wide Angle Extended	9 / 0.349	20 / 0.787	25.8 / 1.020
900-156: one of each small, medium, large, medium extended and large shallow				

WARNINGS

- Contents supplied sterile. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

CAUTION

- U.S. Federal law restricts this device to sale by or on the order of a physician.
- This device SHOULD NOT be used without proper training and preceptorship.
- Sterile unless package has been opened or damaged. Ethylene Oxide gas sterilized.
- **THE INSULATION ON THE ELECTRODE MAY BREAKDOWN AT OUTPUT VOLTAGES ABOVE 6 KV. DO NOT ALLOW THE ELECTROSURGICAL UNIT OUTPUT TO EXCEED 6 KV.**

PRECAUTIONS

- LEEP is a powerful surgical procedure with occasional intra-operative and post-operative complications. You MUST obtain informed signed consent prior to the procedure.
- Severe bleeding may rarely develop after deep excisional procedures and may be difficult to impossible to control by direct pressure, electrofulguration and application of AstraGyn® gel or Monsel's paste. You must have the appropriate instruments for suturing blood vessels in the office.
- Adequate anesthesia and complete visualization of the lesional area should be obtained before initiating LEEP.
- Use only specialty-designed insulated instruments for performing LEEP procedures, to avoid burns to patient or physician.
- For maximum safety, the vaginal walls may be protected from activated excisor and inadvertent cuts and burns by a latex condom placed on any of the LEEP specula. The use of condoms or LEEP vaginal retractor (F400H) is of great help to retract vaginal walls allowing for good visualization of the cervix and vaginal fornix. Check first for patient allergy to latex before using condoms.
- Prior to any cone electrosurgical procedures, make a "dry pass" with the excisor under colposcopic guidance. This will determine that the path of the activated excisor will not be obstructed or come in contact with noninsulated or anesthetized tissue.
- Electrosurgical excision with cone excisor must be relatively slow but continuous.
- To increase good visualization during LEEP, it is necessary to have an adequate fume evacuation vacuum-suction system. Make sure it functions properly before initiating LEEP.
- Make sure that the electrosurgical unit (ESU) is maintained in the "Standby" mode until electrosurgery begins.
- Physicians with limited experience with LEEP will obtain a stable hand by placing the shaft of the cone excisor on the posterior blade of the LEEP speculum.
- The Fischer Cone Biopsy Excisors® are for single-use ONLY. After use, discard of in accordance with all applicable Federal, State and local Medical/Hazardous waste products.

INDICATIONS FOR USE

- The Fischer Cone Biopsy device is indicated for Large Loop Excision of The Transformation Zone (LLETZ) in the diagnosis and treatment of some Cervical Intraepithelial Neoplasias (CIN) and dysplasias.
- Diagnostic Usages (Electroexcision):
 - Intraepithelial neoplasia of lower genital tract
 - Condylomata acuminata
 - Lesions such as polyps, molluscum, contagiosum, nevi, seborrheic keratosis, achricordon, etc.
- Therapeutic Usages (Electrofulguration):
 - Intraepithelial neoplasia of lower genital tract
 - Condylomata acuminata
 - Symptomatic conservative treatment-resistant ectropion
 - Lesions such as polyps, molluscum, contagiosum, nevi, seborrheic keratosis, achricordon, etc.

CONTRAINdications

- Absolute contraindication include hypovolemic shock, cardiac decompensation (Class IV) and coagulopathy/blood dyscrasias; acute/chronic severe, active infection caused by Herpes, Chlamydia or Gonococcus and other aerobic and anaerobic organisms; intraepitheal disease extending to or close to lower uterine segment; severe allergic reaction to local anesthetic and iodine solutions; and a refusal to sign an informed consent form
- Relative contraindications include clinically apparent invasive cancer (punch biopsy is sufficient for histologic diagnosis); pregnancy (except in expert hands when invasion is suspected); less than 3 months postpartum and hypoplastic DES-cervix

**CAUTION: THE INSULATION ON THE ELECTRODE MAY BREAKDOWN AT
OUTPUT VOLTAGES ABOVE 6 KV. DO NOT ALLOW THE
ELECTROSURGICAL UNIT OUTPUT TO EXCEED 6 KV.**

INSTRUCTIONS FOR USE

1. Use the appropriate electrosurgical precautions, coated insulated instruments, dispersing pad, smoke and suction tubing.
2. If needed, stabilize the cervix with a tenaculum. Apply the tenaculum further up in the anterior lip to avoid closeness to the cutting action of the electrode.
3. Clean the cervix and vagina from accumulated secretions. This is done to enhance the electrical transmission.
4. Apply an intracervical local anesthesia of 2-5 cc of Lidocaine with epinephrine 1/100,000.
5. Establish the patency of the endocervical canal using a Q-tip.

6. Choose the proper size Fischer Cone Biopsy Excisor® (the electrode) to accommodate for the anatomical features (size and shape) of the patient's uterine cervix. General guidelines for choosing the correct size include:
 - Small:** small cervix, nulliparous patient
 - Medium:** moderate size cervix, post-menopausal patient
 - Large:** large cervix, multiparous cervix, suspicion of deep endocervical lesions
7. Without applying electricity, position the electrode top in the endocervical canal. Choose the point of entrance. If the cervix has any indentations, choose that as the Entrance point. (The specimen may be fragmented if not.)
8. Set the electrosurgical unit (ESU) at 35-50 watts cutting depending upon the size of the electrode being used:

Small: 30-35 watts	Large: 45-50 watts
Medium: 40-45 watts	Large Shallow: 40-45 watts
Small Wide Angle:	40-45 watts
Small Wide Angle Extended:	40-45 watts
Medium Extended:	40-45 watts

The voltage/power should be kept as low as possible to achieve the desired effect.

- a) The user should refer to the electrical surgical generator labeling for additional contraindications to electrosurgery.
- b) The Fischer Cone Biopsy Excisors are compatible with most ESUs; the CooperSurgical LEEP System 1000, Wallach Quantum 2000 and Valley Lab Force FX Electrosurgical Generators.
9. With the electrode positioned, press "Cut" and gently advance the instrument until the stop arm softly abuts the cervix, then turn clockwise until completing 360° to remove the specimen. While keeping the electrode in place, you may stop and restart the electricity after repositioning your wrist to facilitate completion of the 360° circle. This maintains quality and integrity of a single specimen. Techniques can be learned to obtain excisional biopsies, partial cones, and customized cones.
10. Apply hemostasis techniques as needed (ball electrode, Monsel's solution, sutures).
11. If larger cones are to be removed, preventative hemostatic sutures may be appropriate.

NO LONG TERM FOLLOW-UP

No long term follow-up studies with this device have been performed as to recurrence rates and the effects of loop electrosurgical excision procedures on pregnancy outcome are not known.

ELECTROSURGICAL TISSUE EFFECT

Delivery of continuous sinusoid wave form currents through a small electrode at appropriate power levels can cause rapid heating to the intracellular fluids in the cells in close proximity to the electrode, turning these fluids into steam. The significant increase in volume (approximately 5 times) causes cellular structure to rupture creating the clinical effect of (CUT) with little or no hemostatic effect along the margin of the divided tissue. Delivery of short duration pulses of R.F. currents through a small electrode at appropriate power levels can cause heating of intracellular fluids at a more gradual pace. This allows evaporation of these fluids without rupturing the cellular structure, creating the clinical effect of desiccation or (COAG) without the division of tissue. By varying the pulse to an intermediate duration, it is possible to get the clinical effect to combine or "blend", the clinical characteristics of (CUT) and (COAG) yielding the effect referred to as (BLEND), where tissue is divided with a desirable amount of hemostasis along the margins of the divided tissue.

NOTE: Some electrosurgical generators have output load characteristics that cause the electrosurgical effects to vary considerably throughout the procedure. When these are used, it may be necessary to readjust the relative power during the procedure.

TECHNIQUE GUIDANCE

- The endocervix is commonly not included in the loop excision and the results of the Endocervical Curettage (ECC) do not appear to be predictive of either residual or invasive disease after loop excision procedures. If the ECC is positive for dysplasia, a standard cone biopsy should be considered.
- Loop excision procedures performed with small diameter wire loop electrodes produce multiple small pieces of cervical tissue and provide a less acceptable tissue specimen for histopathologic analysis. The effectiveness of the procedure and influence of electrode design are not completely understood.
- Larger lesions involving multiple quadrants of the cervix are more difficult to remove with either the small or large diameter loop electrodes.

SMALL DIAMETER LOOP VS LARGE DIAMETER LOOP

The histological quality of the specimens obtained using the small diameter loop was inferior to that obtained with the larger diameter loop. This is because more of the epithelium comes in direct contact with the small diameter electrode than with the large diameter electrode and that part of the electrode is not insulated. This results in more thermal damage to the tissue. In 7% of the samples, the amount of thermal damage precluded histological evaluation of the excised tissue. In addition, since multiple, often irregular, strips of epithelium are excised from the cervix, it is often difficult for the pathologist to orient the specimens for optimal histopathological examination.

THERMAL INJURY AND DEFECTS OF TISSUE TREATED WITH LOOP ELECTRODES

Possible injury to cervical tissue may include: (I) thermal coagulation injury of the cervix, up to one-third the thickness of normal epithelium of the cervix, (II) fragmentation of squamous epithelium of the cervix attributable to long exposure periods along the excision site that allows heat to dissipate laterally, and (III) partial coagulation of the endocervical epithelium because of lateral radiation of heat. Loop electrosurgical excision procedures may also produce thermal defects at the periphery of the excised tissue and may make histopathological interpretation difficult or impossible and may prevent accurate diagnosis and the evaluation of the need for further treatment.

RELATED PRODUCTS

- KH1000 LEEP System 1000 Workstation
- 6050P Return Pad - single-use
- 6040 Handpiece - single-use
- 6084 Tubing set - single-use
- 6083 Reducer - single-use
- 52560 Prefilter - single-use
- 6082 ULPA Filter

EXPLANATION OF SYMBOLS

REF	Reorder number
LOT	Batch code
	Use-by date
	Do not resterilize
	Do not re-use
	Do not use if package is damaged
	Consult instructions for use
	Not made with natural rubber latex
	Caution
	Manufacturer
STERILE EO	Sterilized using ethylene oxide
EC REP	Authorized Representative in the European Community
R_x Only	Caution: U.S. Federal law restricts this device to sale by or on the order of a physician

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© 2020 CooperSurgical, Inc.