**INCA<sup>®</sup> Instructions for Use** (English)  
**INCA<sup>®</sup> Notice d'utilisation** (français, French)  
**INCA<sup>®</sup> Istruzioni per l'uso** (Italiano, Italian)  
**Instrucciones de uso de INCA<sup>®</sup>** (español, Spanish)  
**INCA<sup>®</sup> Instruções de uso** (Português, Portuguese)  
**INCA<sup>®</sup> Οδηγίες χρήσης** (Ελληνικά, Greek)  
**INCA<sup>®</sup> 使用説明** (中文, Chinese)  
**INCA<sup>®</sup> Kullanım Kılavuzu** (Türkçe, Turkish)

**CE<sub>2797</sub>**

**STERILE EO**

CooperSurgical<sup>®</sup>

# INCA® Complete Set, Replacement Set and Accessory Pack (English)

## DEVICE DESCRIPTION

The **INCA Complete Set (44-27XX Series)** consists of a nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters, one in-line “T” pressure adapter, a knit cap with hook/loop band, two foam blocks and a hook/loop chin strap (Figures 1 and 4). One INCA Sizing Gauge is included in every box of five sets (Figure 2).

The **INCA Replacement Set (44-07XX Series)** includes the nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters and one in-line pressure “T” adapter (Figure 1).

The **INCA Accessory Pack (44-20XX Series)** includes five knit caps with hook/loop band, five chin straps with ten hook fasteners (not attached) and ten foam blocks.


Figure 1

## INTENDED USE

For administering continuous positive airway pressure (CPAP) in neonates.

## WARNINGS

Use of this device can cause nasal septal erosion, necrosis or permanent nasal damage.

Frequent observation of the neonate and prong position is essential. If nasal irritation is observed, discontinue use immediately.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Dispose of in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.

**PLEASE NOTE:** If a serious incident is suspected from using the INCA Complete Set, Replacement Set and Accessory Pack, report the details of the incident to CooperSurgical via phone number +1 203-601-5200 Ext 3100 or by email at [ProductSurveillance@coopersurgical.com](mailto:ProductSurveillance@coopersurgical.com) and to the local Health Authority in your country. A serious incident may have caused or contributed to a death, a delay in a procedure which resulted in death or serious injury, or a malfunction that could have caused an adverse event.

## INTENDED USERS

Health care professionals in an acute care setting environment.

## PATIENT POPULATION

Preterm neonate, neonate, infant, and pediatric population.

## CONTRAINDICATIONS

Nasal septal erosion or necrosis

## HOW SUPPLIED

Nasal Cannula (only) is packaged sterile for single use.

Do not reuse Nasal Cannula. The rest of the items in set are non-sterile

## INSTRUCTIONS FOR USE – INCA COMPLETE SET

- For appropriate nasal prong fitting, use the INCA Sizing Gauge P/N 025-001 (Figure 2). Wipe gauge with sterile 70% isopropyl alcohol prior to each use. Hold the gauge up to the nasal area to determine appropriate cannula size.

**Note:** Do not place sizing gauge in neonate’s nares.


Figure 2

SIZES LISTED ARE SUGGESTED GUIDELINES ONLY	
NEONATE WEIGHT	CANNULA SIZE
<700g	7.5F or 9F
701 – 1000g	10.5F
1001 – 1300g	12F
1301 – 2000g	15F

- Remove the nasal cannula from its sterile package. Insert the nasal cannula several notches into the ends of the Stay-Flex tubing securely. Expand the tubing to the desired length (Figure 3).


Figure 3

- Disconnect “Y” adapter (Figure 1) from breathing circuit and do not use. Attach inspiratory and expiratory limbs to INCA Stay-Flex tubing. Remove Pressure “T” Adapter Insert (Figure 3) and connect pressure monitor line to INCA “T” adapter. If the ventilator system does not have a separate pressure monitoring line, remove “T” adapter and replace with extra white adapter included in tray.

# INCA<sup>®</sup> Complete Set, Replacement Set and Accessory Pack (English)

- Place knit cap over neonate's occipital area to ensure stability (cap can be cuffed for size.) Leave band open (Figure 4).


Figure 4

- Begin the prescribed gas flow. Fit nasal cannula in neonate's nares. Hold prongs in place until neonate becomes stable and continuous flow begins through the circuit. This may take a minute or two (Figure 5).


Figure 5

**The cannulae prongs should fit the nares snugly without pinching the septum. If the prongs do not fit snugly, the excess movement may cause damage to the mucosa and possible septal erosion.**

**The nasal prongs should fill the nasal opening completely without stretching the skin or putting undue pressure on the nares (blanching around the rim of the nostrils suggests that the prongs are too large).**

**Do not allow the bridge of the cannula to press up against the septum. There should be a small space between the tip of the septum and the bridge between the prongs.**

**Note:** Applying a skin barrier material where the nasal cannula makes contact with the skin may minimize tissue damage. (Ref. NEO-prep skin barrier, small P/N 44-8000 or large, P/N 44-8001).

- Make initial bends in Stay-Flex tubing to conform to neonate's head size and position (i.e., prone, supine).
- Place both left and right lengths of Stay-Flex tubing snugly into slit of foam block (Figure 6).


Figure 6

- Position foam blocks between knit cap and cotton band, and secure by wrapping the band around the head and securing it with the hook/loop tab. This will hold tubing securely and comfortably in place. The foam also provides a cushion between the tubing and the neonate's head (Figure 7).


Figure 7

- Make final adjustments in Stay-Flex tubing to fine-tune prong position. Ideally 1/3 to 1/2 of the prong should rest in neonate's nares.
- Alter the neonate's positioning from the supine to the prone or sideline positions as necessary. Re-adjust the Stay Flex tubing if required.
- Use an oral gastric tube to eliminate gastric insufflation if required.
- To help optimize the effect of CPAP, a chin strap may help prevent mouth leaks. To apply the chin strap to the neonate, place cotton side of the chin strap under the chin (blue lines on exterior side). Ensure blue line side does not touch neonate's face. Apply adhesive side of hook tab to cotton side of chin strap. Fasten both ends of chin strap to knit cap band. Be careful not to fasten the chin strap too tightly.

## INSTRUCTIONS FOR USE – INCA REPLACEMENT SECTIONS

To remove the used tubing and cannula:

- Gently remove the cannula from the nares.
- Lift the tubing out of the head band assembly, leave foam blocks.
- Disconnect the used tubing from the circuit at the white adapter on one side, and at the pressure "T" adapter on the other side. If the ventilator system does not have a separate pressure monitoring line, remove "T" adapter and replace with extra white adapter included in tray.
- Install replacement tubing and cannula per previous instructions.


## INSTRUCTIONS FOR USE – INCA ACCESSORY PACK

The Accessory Pack provides clinicians with extra hats when cannula and hat sizes do not match neonate. The pack also includes extra chin straps, hook-and-loop fasteners and foam blocks.

See "INCA Complete Set" above for instructions for use.

# INCA® Complete Set, Replacement Set and Accessory Pack (English)

## GLOSSARY OF SYMBOLS

Ref #	Symbol	Title	Description	Standard Development Organization
5.1.6		Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1
5.1.5		Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1
5.1.3		Date of manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1
5.1.4		Use-by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1
5.2.3		Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.	ISO 15223-1
5.4.5, Annex B.2		Not made with natural rubber latex	Indicates that natural rubber latex was not used in the manufacturing of the product, its container, or its packaging.	ISO 15223-1
5.4.2		Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1
5.2.8		Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged opened.	ISO 15223-1
5.4.3		Consult instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1
5.1.2		Authorized representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1
5.1.1		Manufacturer	Indicates the medical device manufacturer, as defined in EU Directive 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1
n.a.	<b>R<sub>x</sub> Only</b>	Prescription device	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.	21 CFR 801.109
n.a.		Not made with phthalates	n.a.	n.a.
n.a.		Medical Device	Medical Device	n.a.
n.a.		Product conforms to the Medical Device Directive 93/42/EEC	Signifies European technical conformity.	n.a.

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