
SM610 • MiniSpace®

IUI Catheter

For intrauterine artificial insemination

Instructions for Use (English)

Length: 17.2 cm (effective) 19.9 cm (overall)
O.D.: 2.6 mm stepped down to
1.6 mm for distal 4.5 cm of length
I.D.: 0.7 mm for the 17.2 cm effective length

WARNING: NEVER ATTEMPT TO RESTERILIZE.

MOUSE EMBRYO ASSAY (MEA) 2-CELL SURVIVAL >80%

This device is not intended for use for intrafallopian tube procedures, In vitro fertilization (IVF) and Gamete intrafallopian transfer (GIFT).

R_x Only

CE2797

DEVICE DESCRIPTION

The MiniSpace® IUI Catheter is a single-use, sterile, disposable, flexible catheter designed for use in intrauterine insemination. The opening is at the tip of the catheter which has been rounded and beveled for the comfort of the patient.

The device has an overall length of 19.9 cm including a 2.7 cm long polyethylene Luer-lock hub which is permanently affixed to its proximal end.

The effective portion of the catheter (intended for insertion) comprises a 17.2 cm long clear polyethylene tube with a 1.6 mm O.D. (outer diameter) and a 0.7 mm I.D. (inner diameter). When filled, the MiniSpace contains 0.07 ml of fluid.

The proximal 12.7 cm length of the catheter is tightly surrounded by a clear polypropylene sheath with a 2.6 mm O.D. for all but the distal 1 mm length, which joins the inner tube to the outer sheath. At this 1 mm joining length, the O.D. tapers inward from 2.6 mm to 1.6 mm over the 1 mm spindle-shaped portion. This leaves 4.5 cm of the inner catheter tube exposed at the distal end.

The proximal inner tube and the outer sheath circumference are molded directly onto the distal end of the Luer-lock hub so that only the lumen of the inner tube can provide a channel for flow. This design results in a single-channel catheter with two segments each having different degrees of flexibility (i) a highly flexible 1.6 mm O.D., 4.5 cm long distal portion and (ii) a more rigid 2.6 mm O.D., 12.7 cm long proximal portion.

Graduation markings are at 5, 6, 7, 8, 9, and 10 cm distance from the extreme distal tip of the catheter.

INDICATIONS FOR USE

Intrauterine artificial insemination

CONTRAINDICATIONS

The MiniSpace shall not be used in:

- pregnant patients or in patients suspected of being pregnant
- the presence of acute or chronic cervical infection
- the presence of, or after recent Pelvic Inflammatory Disease (PID)

WARNINGS

Always prepare the spermatozoa to separate the motile sperm from the rest of the semen sample before performing Intrauterine Insemination (including Intracervical Artificial Insemination) using, for example, the Sperm Select® system.

The introduction of unwashed spermatozoa into the uterus with the MiniSpace or by any other means will almost always result in severe adverse

reactions, which may include anaphylactic shock. Please refer to published medical literature for methods of preparing spermatozoa for Intrauterine Insemination before performing this procedure.

Contents supplied sterile. Do not use if sterile barrier is damaged.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

PRECAUTIONS

Care should be taken prior to catheter insertion to ascertain the depth of the uterus and any uterine retroflexion or antelexion present in each patient. (This should be performed during a previous examination and not on the day of insemination when disturbance of the uterus should be minimized.)

In no case shall the device be forced against resistance. If passage of the MiniSpace through the cervical canal cannot be achieved, the device should not be forced and the possibility of pathologic cervical stenosis should be considered. In no case should air be injected into the uterine cavity.

ADVERSE REACTIONS

As with any instrument which is passed through the internal cervical os, mild cramping may be experienced by patients. In every case the INSTRUCTIONS FOR USE shall be followed: taking note of any CONTRAINDICATIONS, WARNINGS and PRECAUTIONS. ALWAYS USE WASHED SPERMATOZOA. REFER TO WARNINGS.

INSTRUCTIONS FOR USE

1. a. With the patient in the dorsal lithotomy position, expose the uterine cervix to view with a vaginal speculum.
b. Firmly attach a syringe (preferably 3 ml, not included) with the plunger fully depressed, to the Luer-lock hub at the proximal end of the MiniSpace.
(NOTE: If Sperm Select System is used to prepare spermatozoa, proceed directly to Step 2.)
c. Place the distal tip of the MiniSpace into the specimen of PREPARED SPERMATOZOA

so that the tip of the catheter is immersed in the specimen. Aspirate the desired quantity (usually 0.5 ml) of the specimen into the lumen of the catheter by retraction of the syringe plunger. The clear tubular portion of the MiniSpace has a capacity of approximately 0.07 ml of fluid.

2. ALTERNATIVELY:

- a. A syringe filled with washed spermatozoa may be attached to the Luer-lock hub.
- b. Fill the catheter to the end with washed spermatozoa prior to inserting.
- c. Inseminate 0.5 ml of **WASHED SPERMATOZOA**.

3. Having previously determined the depth and direction of the patient's uterus, gently introduce the distal tip of the MiniSpace into the external cervical os and advance it through the cervical canal into the uterine cavity until resistance is felt.

The distal tip of the catheter should then be at the uterine fundus. Now retract the MiniSpace approximately 1 cm. The markings on the catheter tube provide further guidance in estimating the length of catheter insertion and retraction.

NOTE: If any difficulty of insertion is encountered and bending of the catheter occurs, this can frequently be overcome by rotating the catheter gently while advancing it, or by grasping it close to the cervix with a ring forceps which is then used for advancement. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction of the cervix with a tenaculum may be required to perform full insertion.

In no instance should force be used against resistance to achieve insertion.

4. With the distal tip of the MiniSpace now positioned 1 cm from the uterine fundus, SLOWLY INJECT THE PREPARED SPERMATOZOA into the uterine cavity by depressing the syringe plunger fully. Should there be any doubt about the complete expulsion of the specimen in the uterus, detach the syringe from the Luer-lock catheter hub, load the syringe with approximately 0.1 ml of sterile sodium chloride 0.9% (normal saline) solution, reattach the syringe to the Luer-lock catheter hub and, SLOWLY inject the sterile sodium chloride 0.9% solution in the insemination catheter.

5. Following injection of the specimen, slowly withdraw the MiniSpace from the uterus and discard in accordance with Federal, State and local Medical / Hazardous waste practices. Remove the vaginal speculum and allow the patient to rest in the supine position for at least 10 minutes following the insemination.

EXPLANATION OF SYMBOLS

- | | |
|--|--|
|  Reorder number |  Sterilized using ethylene oxide |
|  Batch code | R_xOnly CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician. |
|  Use-by date |  Authorized Representative in the European Community. |
|  Caution |  Not made with natural rubber latex |
|  Do not use if the package is damaged or opened |  Manufacturer |
|  Consult instructions for use | |
|  Do not re-use | |
|  Do not re-sterilize | |

MiniSpace® and Sperm Select® are registered trademarks of CooperSurgical, Inc.

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© 2021 CooperSurgical, Inc.
MiniSpace® e Sperm Select® sono marchi registrati di CooperSurgical, Inc.