
# 10068 • MityOne® MitySoft® Bell Cup

## Vacuum-Assist Delivery System • Instructions for Use (English)

**CAUTION:** U.S. Federal law restricts this device to sale by or on the order

**DEVICE DESCRIPTION**

The MityOne® MitySoft® Bell Cup is designed to assist a practitioner in the delivery of an infant during childbirth. The MityOne is an integrated pump and cup system, packaged, and sterilized for immediate use.

**INDICATIONS FOR VACUUM-ASSISTED VAGINAL DELIVERY**

**No indication is absolute.**

Fetal:

- Nonreassuring fetal status
- Failure to deliver spontaneously following an appropriately managed second stage of labor

Maternal:

- Need to avoid voluntary expulsive efforts
- Inadequate expulsive efforts

**CONTRAINDICATIONS**

Do not initiate vacuum if any of the following conditions exist:

- Non-vertex positions (breech or transverse lie position) or face or brow presentation
- Suspected cephalopelvic disproportion
- Less than 34 weeks gestation
- Unengaged vertex
- Incompletely dilated cervix
- Previous scalp sampling
- Suspected macrosomia, or risk of shoulder dystocia
- Failed vacuum or forceps attempt
- Need for active device rotation
- Suspected fetal bleeding abnormalities

**ADVERSE EVENTS**

Fetal Injuries: Head trauma, bruises, contusions, lacerations, scalp edema, skull fracture, cephalhematoma, subgaleal hematoma, subdural hemorrhage, parenchymal hemorrhage, intracranial hemorrhage, retinal hemorrhage

Maternal Injuries: Vaginal, cervical, uterine, bladder, rectal tissue tears

**WARNINGS**

- This product is for single use and not to be reused, reprocessed or resterilized due to risk of contamination and/or cause patient/fetal infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient/fetus to another.
- DO NOT exceed recommended vacuum levels.
- Never apply the vacuum cup to any portion of the infant's face.
- Vacuum-assisted delivery should only be performed or supervised by a trained and experienced operator.
- "There should be a willingness to abandon attempts at vacuum extraction if satisfactory progress is not made." (ACOG Technical Bulletin #196, Aug. 1994)
- Abandon vacuum-assisted delivery if:
  - a. Vacuum cup becomes disengaged (pops off) three times
  - b. Vertex has not advanced substantially with each traction attempt
  - c. There is evidence of fetal scalp trauma
  - d. Cumulative traction time exceeds 10 minutes or total vacuum procedure time exceeds 15 to 30 minutes
- If the extractor cup becomes disengaged, check for trauma to the fetal scalp before reapplying.

**PREREQUISITES FOR VACUUM-ASSISTED VAGINAL DELIVERY**

- Term infant
- Ruptured membranes
- Empty bladder and nondistended rectum
- Engaged vertex position
- Complete cervical dilation and effacement
- Willingness to abandon procedure and plan for alternate delivery route

These instructions are intended as general guidelines. Practitioners should refer to current institutional and recognized guidelines that address vacuum-assisted delivery.

**VACUUM-ASSISTED VAGINAL DELIVERY GUIDELINES**


1. Remove the protective lid from the cup before use.
2. Check the integrity of the MityOne Obstetrical Vacuum-Assist Delivery System by pressing the cup to the palm of your gloved hand and applying vacuum. The gauge needle should remain steady. (NOTE: The gauge needle should rest in the "ZERO" box when vacuum is not applied.)

**CAUTION: DO NOT USE DEVICE IF GAUGE NEEDLE DOES NOT INDICATE ZERO IN THE ABSENCE OF A VACUUM.**

3. Carefully examine fetal presentation and position prior to positioning the cup over the flexion point.

**CAUTION: DO NOT USE PRODUCT IF FLEXION POINT IS NOT ACCESSIBLE.**


4. Wipe the scalp as clean as possible.
5. To insert the MityOne MitySoft Bell Cup, grasp the cup with the fingers and fold edges inward to facilitate insertion (see **Figure 1**).

**Figure 1 – MitySoft Bell Cup**

6. Separate the labia with the other hand. Press the cup portion downward and inward over the

fourchette to make contact with the fetal scalp. Apply the center of the cup over the flexion point (see **Figure 2**). The flexion point can be located by identifying the posterior fontanelle and then moving the finger anteriorly along the sagittal suture approximately 1.2 in. (3 cm).

**Figure 2 - Infant head**

The sagittal suture should pass under the middle of the cup.

**CAUTION: NEVER APPLY THE CUP TO ANY PORTION OF THE INFANT'S FACE.**

7. Gently move a finger around the rim of the cup to remove intrusion of extraneous tissues and to ensure proper cup placement.
8. With the cup placed over the flexion point, raise the vacuum level to approximately 3.9 in. Hg (10 cm Hg) (yellow shaded area on gauge face) to initiate vacuum. Recheck the rim of the cup for interposed extraneous tissue.

- With the onset of contraction, rapidly raise the vacuum to **15-22.8 in. Hg (38-58 cm Hg)** (green shaded area on gauge face) and begin applying traction.<sup>1,2</sup> Always apply traction in harmony with contractions and along the pelvic axis.

**CAUTION: DO NOT EXCEED RECOMMENDED VACUUM LEVELS.**

- When the contraction is no longer effective, discontinue traction. You may reduce the vacuum to approximately **3.9 in. Hg (10 cm Hg)** (yellow shaded area on gauge face) while awaiting the next contraction or maintain vacuum at current level. Reduce vacuum by slowly pushing or pulling on the vacuum release lever

Figure 3 – Vacuum release lever

- Recheck for interposed extraneous tissue prior to each tractive effort.
- If traction is misaligned or too forceful, the vacuum cup may disengage (pop off). In case of pop off, check fetal scalp for trauma before reapplying vacuum cup.

**CAUTION: ABANDON VACUUM-ASSISTED DELIVERY IF THE VACUUM CUP DISENGAGES (POPS OFF) THREE TIMES.**

- With each successive contraction, draw the head gently over the perineum. Once the head is delivered, release vacuum (see Figure 3) and remove the cup. Continue delivery in the usual manner. **In the unlikely event that vacuum does not release, use scissors and cut the stem, being careful to avoid the fetal scalp.**

**CAUTION: DO NOT OPERATE VACUUM PUMP AT TRACTIVE LEVELS FOR MORE THAN 10 CUMULATIVE MINUTES OF TRACTION OR TOTAL PROCEDURE TIME OF 15 TO 30 MINUTES.**

- Carefully examine infant's head and observe vital signs at regular intervals to ensure infant's well-being.
- Dispose of the MityOne Delivery System in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.
- Document the use of the MityOne Delivery System and notify nursery staff per hospital protocol.

**PREREQUISITES FOR VACUUM-ASSISTED CESAREAN DELIVERY**

- Proper indication for operative delivery
- Proper indication for use of the vacuum-assisted delivery device

**VACUUM-ASSISTED CESAREAN DELIVERY GUIDELINES**

- Remove the protective lid from the cup before use.
- Check the integrity of the MityOne Obstetrical Vacuum-Assist Delivery System by pressing the cup to the palm of your gloved hand and applying vacuum. The gauge needle should remain steady. (NOTE: the gauge needle should rest in the "ZERO" box when vacuum is not applied.)

**CAUTION: DO NOT USE DEVICE IF GAUGE NEEDLE DOES NOT INDICATE ZERO IN THE ABSENCE OF A VACUUM.**

- Enter the uterus routinely and assess fetal head position.
- If the fetal head is high and accessible beneath the uterine incision, wipe the scalp as clean as possible and place vacuum cup over the flexion point. If the head is low and inaccessible, place gloved fingers under the head and flex upward to bring the scalp under the uterine incision. Wipe scalp as clean as possible and apply the vacuum cup to the flexion point.

**CAUTION: NEVER APPLY CUP TO ANY PORTION OF THE INFANT'S FACE.**

- Raise the vacuum to **15-22.8 in. Hg (38-58 cm Hg)** (green shaded area on gauge face) and begin applying traction to gently bring the head upward through the incision.<sup>1,2</sup>
- Once the head is delivered, release vacuum by pushing or pulling on the vacuum release lever (see Figure 3) and remove the cup. Continue delivery in the usual manner. **In the unlikely event that vacuum does not release, use scissors and cut the stem, being careful to avoid the fetal scalp.**
- Carefully examine infant's head and observe vital signs at regular intervals to ensure infant's well-being.
- Dispose of the MityOne Delivery System in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.
- Document the use of the MityOne Delivery System and notify nursery staff per hospital protocol.

**EXPLANATION OF SYMBOLS**

	Reorder number		Keep dry
	Batch code		Consult instructions for use.
	Use-by date		Sterilized using ethylene oxide
	Not made with natural rubber latex		Authorized Representative in the European Community
	Do not re-use		<b>Caution:</b> U.S. Federal law restricts this device to sale by or on the order of a physician
	Manufacturer		Do not use if package is damaged
	Caution		

© 2020 CooperSurgical, Inc.

<sup>1</sup> Bofill, James A., MD, et al., "The Mississippi Operative Vaginal Delivery Trial: Lessons Learned", *Contemporary OB/GYN*, October 1998.  
<sup>2</sup> Schifrin, Barry S., MD, et al., "Cutting Your Legal Risks with Vacuum Assisted Delivery", *OBG Management*, June 1994.

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