Mobius®

Elastic Retractor Systems

Instructions for Use

Elastische Retraktorsysteme

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


OR/Mobius® Elastic Retractor System

Product Image	Part Number	Description	Incision Range	Inner Diameter	Working Area	
	900-508	Mobius, Mini	2 - 6 cm	8 cm	2 cm Incision  3 cm ² Working Area	6 cm Incision  28 cm ² Working Area
	900-500	Mobius, Small	4 - 9 cm	11 cm	4 cm Incision  12.5 cm ² Working Area	9 cm Incision  64 cm ² Working Area
	900-515	Mobius, Medium	8 - 14 cm	15 cm	8 cm Incision  50 cm ² Working Area	14 cm Incision  154 cm ² Working Area

OB/Mobius® Elastic Retractor System


Product Image	Part Number	Description	Incision Range	Inner Diameter	Working Area	
	900-517	Mobius, Medium, Large	10 - 15 cm	17 cm	10 cm Incision  78 cm ² Working Area	15 cm Incision  176 cm ² Working Area
	900-520	Mobius, Large	12 - 17 cm	22 cm	12 cm Incision  113 cm ² Working Area	17 cm Incision  227 cm ² Working Area

DEVICE DESCRIPTION

The Mobius® Elastic Retractor provides 360° of circumferential retraction while lining the wound to protect against contamination. The Mobius Elastic Retractor is comprised of three components:

- an open-ended, flexible yet non-compliant (will not stretch) polyurethane sleeve ①
- an internal, flexible O-Ring (blue) ② and
- a Quad-Ring (yellow) ③

The O-Ring and Quad Ring components are made entirely of Ethyl-Vinyl-Acetate (EVA). The polyurethane sleeve is the single component contacting the patient. See the diagram below:

The yellow Quad-Ring is injection-molded from Ethyl-Vinyl-Acetate (EVA) into the form known as the “Mobius” configuration. This configuration creates a torsional force that both assists the physician in the operation of the device and helps to maintain desired wound dilation once installed.

Contents: One (1) Mobius Elastic Retractor per individual sterile pouch.
Each order contains five (5) individual units.

Storage: Store at room temperature.

WARNING

The Mobius Elastic Retractor is a Single Use Device. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of the device and packaging in accordance with hospital, administrative and/or local government policy.

Do not use if the packaging is open or damaged.

When using Mobius during a cesarean section, do not place the baby’s face against the sleeve of the retractor as it could result in suffocation.


CAUTION

U.S. Federal law restricts this device to sale by or on the order of a physician.

PRECAUTION

This document is designed to provide instructions for the use of the Mobius Elastic Retractor and is not intended as a guide to surgical technique.

PLEASE READ ALL INFORMATION CAREFULLY.

INDICATIONS FOR USE

- The OR/Mobius Elastic Retractor is indicated for use to assist in laparotomy procedures. It is intended to provide incision retraction and to protect against wound contamination during open surgery.
- The OB/Mobius Elastic Retractor is indicated for use to assist in non-urgent cesarean deliveries that are routine procedures. It is intended to provide incision retraction and to protect against wound contamination during a cesarean section. It is indicated for use as a surgical retractor for both vertical and transverse incisions.





CONTRAINDICATIONS

The OB/Mobius Elastic Retractor is contraindicated for use in women with significant adhesions.

INSTRUCTIONS FOR USE

OR/Mobius[®] Elastic Retractor System

Used in a variety of general, Colorectal Bariatric and OB-GYN procedures

Product Image	Part Number	Description	Incision Range	Inner Diameter	Working Area	
	900-508	Mobius, Mini	2 - 6 cm	8 cm	2 cm Incision  3 cm ² Working Area	6 cm Incision  28 cm ² Working Area
	900-500	Mobius, Small	4 - 9 cm	11 cm	4 cm Incision  12.5 cm ² Working Area	9 cm Incision  64 cm ² Working Area
	900-515	Mobius, Medium	8 - 14 cm	15 cm	8 cm Incision  50 cm ² Working Area	14 cm Incision  154 cm ² Working Area

Step 1: Prepare the surgical site according to standard hospital procedures.

Step 2: Using a skin marking pencil, mark a transverse incision line at the desired location.

Step 3: The fat is then cleared digitally from the midline superiorly and inferiorly to expose the fascia in the vertical axis.

Step 4: The anterior rectus fascia is then incised in a vertical direction through the full length of the previously cleared area.









Step 5: Grasp the blue O-Ring, collapse it into a column and insert through the incision into the abdominal cavity.

Step 6: Gently grasp the Yellow Quad-Ring of the Mobius Retractor at the 11 o'clock and 1 o'clock position and roll outward. Ensure an equal distribution of the sleeve around the entire circumference of the yellow ring and repeat until the retractor sits snug against the patient's skin.

Step 7: Proceed with surgery. Note: If the Mobius Elastic Retractor becomes loose during surgery, simply roll the yellow ring 1-2 more times until it is once again snug against the patient's skin.

REMOVAL:


Step 8: To remove the Mobius Elastic Retractor after use, insert a finger through the Mobius Elastic Retractor, grasp the blue O-Ring and pull it up and out of the incision.

INSTRUCTIONS FOR USE

OB/Mobius® Elastic Retractor System

Used in C-sections and open Bariatric and Colorectal procedures

Product Image	Part Number	Description	Incision Range	Inner Diameter	Working Area	
	900-517	Mobius, Medium, Large	10 - 15 cm	17 cm	 <p>10 cm Incision 78 cm² Working Area</p>	 <p>15 cm Incision 176 cm² Working Area</p>
	900-520	Mobius, Large	12 - 17 cm	22 cm	 <p>12 cm Incision 113 cm² Working Area</p>	 <p>17 cm Incision 227 cm² Working Area</p>

TECHNIQUE 1: ASSISTED ROLL-DOWN

- Step 1:** Prepare the surgical site according to standard hospital surgical procedure.
- Step 2:** Make the desired transverse (or vertical) incision at the desired location.
- Step 3:** The peritoneum is opened; introduce a hand through the incision and elevate the abdominal wall to ensure that no significant adhesions or obstructions are present. The inner, blue ring must be placed against the peritoneum a full 360° (see Figure 1A).
- Step 4:** The circulating Nurse opens the sterile pouch and places the OB/Mobius Retractor on the sterile field. When the attending Surgeon calls for the retractor, the surgical Nurse grasps the outer, yellow ring and raises the device to allow the sleeve to fully extend with the blue ring below the yellow ring. The Nurse hands the device to the Surgeon who grasps the outer, yellow ring at a convenient single point on the circumference with the non-dominant hand (see Figure 2A). With the other hand the Surgeon grasps the blue ring at two points opposite of one another and collapses the blue ring into a column. Note: The sleeve is drawn into a narrow band when the blue ring is collapsed (see Figure 3A).
- Step 5:** While maintaining the column shape of the blue ring, the Surgeon inserts the upper portion of the blue ring into the upper portion of the incision (see Figure 4A) making sure that it is in the peritoneal cavity and there is no obstruction to the placement of the OB/Mobius Retractor. The blue ring is then inserted interiorly far enough to allow the ring to clear the incision. A second hand is used to place the lower portion of the blue ring into the inferior portion of the incision (see Figure 5A). At this point the blue ring will automatically return to its original shape.
- Step 6:** The physician and the assistant grasp opposite sides of the outer, yellow ring gently pulling the yellow ring upward to remove any slack in the sleeve (see Figure 6A). This slight tension upon the incision will create a “window” in the incision at the center of the device. The Surgeon then places a hand through the center portion of the OB/Mobius Retractor to make sure that the blue ring is not trapping any organs or omentum (see Figure 7A). Any structures detected must be cleared before proceeding.

continued

Step 7: Roll-down procedure: The Surgeon grasps the outer, yellow ring between the thumb and forefinger of each hand. The thumbs are placed on the top of the ring with both hands near the distal end of the incision approximately at the “11 o’-clock and 1 o’clock” positions (see Figure 8A). The wrists are rotated outward and both hands are moved down the ring. As the hands move down, the assistant grasps the proximate end of the yellow ring (see Figure 9A) to prevent it from unrolling. The Surgeon then moves his hand downward along the yellow ring (see Figure 10A) until the first layer is completely folded (see Figure 11A). The OB/Mobius Retractor is then placed in slight tension to remove any looseness in the sleeve and a second roll is accomplished in a similar fashion. If the OB/Mobius Retractor is still noted to be slack, a third roll should be attempted. If the thickness of the abdominal wall doesn’t permit this, two rolls may be all that is necessary. Some patients may require four rolls. At this point in the procedure, the physician is ready to proceed with the surgery (see Figure 12A).

REMOVAL:

Step 8: Once surgery is completed, removal is accomplished by the Surgeon inserting a single digit through the wound and hooking the OB/Mobius Retractor inner blue ring. The inner, blue ring is then drawn up and through the wound (see Figure 13A).

Figure 1A

Figure 2A

Figure 3A

Figure 4A

Figure 5A

Figure 6A


Figure 7A

Figure 8A


Figure 9A


Figure 10A

Figure 11A


Figure 12A

Figure 13A

TECHNIQUE 2: UNASSISTED ROLL-DOWN

Step 1: Prepare the surgical site according to standard hospital surgical procedure.

Step 2: Make the desired transverse (or vertical) incision at the desired location.

Step 3: The peritoneum is opened; introduce a hand through the incision and elevate the abdominal wall to ensure that no significant adhesions or obstructions are present. The inner, blue ring must be placed against the peritoneum a full 360° (see Figure 1U).

Step 4: The circulating Nurse opens the sterile pouch and places the OB/ Mobius Retractor on the sterile field. When the attending Surgeon calls for the OB/ Mobius Retractor, the surgical Nurse grasps the outer, yellow ring and raises the device to allow the sleeve to fully extend with the blue ring below the yellow ring. The Nurse hands the device to the Surgeon who grasps the outer, yellow ring at a convenient single point on the circumference with the non-dominant hand (see Figure 2U). With the other hand the Surgeon grasps the blue ring at two points opposite of one another and collapses the blue ring into a column. Note: The sleeve is drawn into a narrow band when the blue ring is collapsed (see Figure 3U).

Step 5: While maintaining the column shape of the blue ring, the Surgeon inserts the upper portion of the blue ring into the upper portion of the incision (see Figure 4U) making sure that it is in the peritoneal cavity and there is no obstruction to the placement of the OB/ Mobius Retractor. The blue ring is then inserted inferiorly far enough to allow the ring to clear the incision. A second hand is used to place the lower portion of the blue ring into the inferior portion of the incision (see Figure 5U). At this point the blue ring will automatically return to its original shape.

Step 6: The physician grasps the opposite sides of the outer yellow ring and gently pulls the yellow ring upward removing any slack in the sleeve (see Figure 6U). This slight tension upon the incision will create a “window” in the incision at the center of the device. The Surgeon then places a hand through the center portion of the OB/ Mobius Retractor to make sure that the blue ring is not trapping any organs or omentum (see Figure 7U). Any structures detected must be cleared before proceeding.

Step 7: Roll-down procedure: The Surgeon grasps the outer, yellow ring between the thumb and forefinger of each hand. The thumbs are placed on the top of the ring with both hands near the distal end of the incision approximately at the “11 o’-clock and 1 o’clock” positions (see Figure 8U). Both wrists are then turned slightly inward creating a depression in the yellow ring (see Figure 9U). Next both wrists are rotated upward and outward (see Figure 10U). After the ring has inverted at the top, the Surgeon grasps the partially inverted bottom portion (the 6 o’clock position) of the ring (see Figure 11U) and twists it downward causing the outer yellow ring to completely invert and at the same time shortening the soft, flexible sleeve. The first roll is complete. The OB/Mobius Retractor is then placed in slight tension to remove any looseness in the sleeve and a second roll is accomplished in a similar fashion. If the OB/Mobius Retractor is still noted to be slack, a third roll should be attempted. If the thickness of the abdominal wall doesn’t permit this, two rolls may be all that is necessary. Some patients may require four rolls. At this point in the procedure, the physician is ready to proceed with the surgery (see Figure 12U)

REMOVAL:

Step 8: Once surgery is completed, removal is accomplished by the Surgeon inserting a single digit through the wound and hooking the OB/ Mobius Retractor inner blue ring. The inner, blue ring is then drawn up and through the wound (see Figure 13U).

Figure 1U

Figure 2U

Figure 3U

Figure 4U


Figure 5U

Figure 6U


Figure 7U

Figure 8U


Figure 9U


Figure 10U


Figure 11U


Figure 12U



Figure 13U

SPECIFICATIONS

Part Number	Product Description	Height (in/mm)
900-508	Mobius, Mini	6.75 (171.45)
900-500	Mobius, Small	6.00 (152.40)
900-515	Mobius, Medium	6.50 (165.1)
900-517	Mobius, Medium-Large	9.90 (251.46)
900-520	Mobius, Large	10.65 (270.51)

Part Number	Product Description	O-Ring c/s (in/mm)	Quad-Ring c/s (in/mm)
900-508	Mobius, Mini	0.222 (5.639)	0.252 (6.401)
900-500	Mobius, Small	0.312 (7.925)	0.400 (10.160)
900-515	Mobius, Medium	0.318 (8.077)	0.400 (10.160)
900-517	Mobius, Medium-Large	0.500 (12.700)	0.580 (14.732)
900-520	Mobius, Large	0.504 (12.802)	0.750 (19.050)

EXPLANATION OF SYMBOLS

REF Reorder number

LOT Lot number

 Use-by date

 Consult instructions for use

 Caution

 Do not re-use

 Do not re-sterilize

 Not made with natural rubber latex

 Do not use if package is damaged

STERILE EO Sterilized using ethylene oxide


R_xOnly Caution: U.S. Federal law restricts this device to sale by or on the order of a physician

EC REP Authorized Representative in the European Community


 Manufacturer

CooperSurgical is a registered trademark of CooperSurgical, Inc.

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CooperSurgical è un marchio depositato di CooperSurgical, Inc.