NEO-fit™ Neonatal Endotracheal Tube Holder

42-2540 • Instructions for Use • English


CooperSurgical®

DEVICE DESCRIPTION

Includes one plastic tube grip with adjustable hook and loop strap and adhesive foam foot pads, two foam "lollipop" securing tapes. Adjustable hook and loop strap allows rapid release and adjustment of ET tube. This product is non-sterile.

Figure 1

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

INTENDED USE

For stabilizing a 2.5 mm to 4.0 mm ID Neonatal ET tube.

INTENDED USERS

The NEO-fit™ is intended to be used by medical professionals only (including physicians, respiratory therapists, nurses and emergency medical staff); they may be used in both emergency situations or under controlled conditions where quick and simple stabilization and support are necessary.

WARNINGS

- Separation of the plastic base foot from the adhesive foam pad may result in potential ET tube extubation. If signs of delamination are observed, replace with a new NEO-fit.
- Usage of ointments or oils may adversely affect the adhesion properties of the adhesive foam pad.
- Excessive salivation may adversely affect the attachment of this device or cause the device to fail during use.
- Discontinue use immediately if skin irritation occurs.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection,

including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

PLEASE NOTE: If a serious incident is suspected from using the NEO-fit, report the details of the incident to CooperSurgical via phone number +1 203-601-5200 Ext 3100 or by email at ProductSurveillance@coopersurgical.com and to the local Health Authority in your country. A serious incident may have caused or contributed to a death, a delay in a procedure which resulted in death or serious injury, or a malfunction that could have caused an adverse event.

PATIENT POPULATION

The NEO-fit is intended to be used on patients in both emergency situations or under controlled conditions where quick and simple stabilization and support are necessary, including but not limited to the hospital, ER, ICU, or point of care. The specific patient population is relative to the product in use; the NEO-fit device is designed for neonate, infant, or pediatric patient use.

PRECAUTIONS

- Not for use with 2.0 mm ID or smaller ET tube.
- Do not use as sole method of support for ventilator tubes; use a tubing support device.
- Patient must be under professional supervision while intubated.
- Do not wrap the tube holder too tightly or the ET tube may be deformed and prevent suction catheter passage.
- Do not wrap the tube holder too loosely or the ET tube may move in the tube holder.

ADVERSE EVENTS

- Skin irritation
- ET tube extubation

Instructions for Use • English (Continued)

INSTRUCTIONS FOR USE


1. Trial fit the NEO-fit™ on patient prior to removing the adhesive liner, (adhesive pad should not overlap the mouth or eyes). Usage of a pectin-based skin preparation may minimize the risk of contact dermatitis.

2. Remove the adhesive liners from the bottom of the adhesive pads.

See Figure 2.


Caution: Do not peel the adhesive pad away from the plastic base during removal of the adhesive liner or the adhesive may be adversely affected. Use a new NEO-fit if this occurs.

Figure 2

3. Affix the adhesive pads to the patient adjacent to the ET tube. See Figure 3.


Figure 3

4. Wrap the tube holder strap snugly around the ET tube by pulling from the end of the strap. Pushing the strap around the tube may not secure the tube adequately. Do not wrap too tightly as this may interfere with the passage of a suction catheter through the ET tube.


See Figure 4.

Figure 4

5. Secure the strap using the hook and loop material of the strap. If needed, secure using "lollipop" foam adhesive tape, included, or other tape as appropriate. See Figure 5.

Figure 5

Instructions for Use • English (Continued)

MR Conditional

The Neo-fit™ Neonatal Endotracheal Tube Grip/Holder was determined to be MR-conditional.

Non-clinical testing demonstrated that the Neo-fit Neonatal Endotracheal Tube Grip/Holder is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the Neo-fit Neonatal Endotracheal Tube Grip/Holder produced the following temperature rise during MRI performed for 15-min of scanning (i.e., per pulse sequence) in the 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system:

Highest temperature change +1.7°C (35.06 °F)

Therefore, the MRI-related heating experiments for the Neo-fit Neonatal Endotracheal Tube Grip/Holder at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 -W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.7°C (35.06 °F).

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Neo-fit Neonatal Endotracheal Tube Grip/Holder. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 20-mm relative to the size and shape of the Neo-fit Neonatal Endotracheal Tube Grip/Holder.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	1,980-mm ²	1,453-mm ²	4,381-mm ²	3,678-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

GLOSSARY OF SYMBOLS

Marks bearing ™ are trademarks of CooperSurgical, Inc.
CooperSurgical is a registered trademark of CooperSurgical, Inc.

Ref #	Symbol	Title	Description	Standard Development Organization
5.1.6		Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1
5.1.5		Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1
5.1.4		Date of manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1
5.1.3		Use-by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1
5.4.5, Annex B.2		Not made with natural rubber latex	Indicates that natural rubber latex was not used in the manufacturing of the product, its container, or its packaging.	ISO 15223-1
5.4.2		Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1
5.4.4		Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	ISO 15223-1
5.4.3		Consult instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1
Table 2, Figure 6		MR Conditional	An item with demonstrated safety in the MR environment within defined conditions.	ISO 15223-1
5.2.7		Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.	ISO 15223-1
5.1.2		Authorized representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1
5.1.1		Manufacturer	Indicates the medical device manufacturer, as defined in EU Directive 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1
n.a.		Prescription device	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.	21 CFR 801.109
n.a.		Product conforms to the Medical Device Directive 93/42/EEC	Signifies European technical conformity.	n.a.

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I marchi contrassegnati con ™ sono marchi di CooperSurgical, Inc.
CooperSurgical è un marchio registrato di CooperSurgical, Inc.
CooperSurgical Distribution B.V.