SM400 • Select IUI®

Intrauterine Insemination Catheter

For intrauterine artificial insemination

WARNING: NEVER ATTEMPT TO RESTERILIZE.

MOUSE EMBRYO ASSAY (MEA) 2-CELL SURVIVAL >80%

This device is not intended for use for intrafallopian tube procedures,
In vitro fertilization (IVF) and Gamete intrafallopian transfer (GIFT).

DEVICE DESCRIPTION

The Select IUI® is a single-use, sterile, disposable, flexible catheter designed for use in intrauterine artificial insemination. Two lateral openings opposite one another, of a 1.2 mm diameter each, are located at 1.5 mm distance from the extreme distal, rounded, smooth, closed tip of the catheter.

The device has an overall length of 20 cm including a 3 cm long polyethylene Luer-lock hub which is molded onto its proximal end. The effective portion of the catheter (intended for insertion) comprises a 17 cm long clear polyethylene tube with a 2.1 mm O.D. (outer diameter) and a 1.6 mm I.D. (inner diameter).

The proximal 12.5 cm length of the catheter is tightly surrounded by an outer, clear polypropylene sheath with a 2.6 mm O.D. and a 2.15 mm I.D. for all but the distal 1 mm length, which joins the inner tube to the outer sheath. At this 1 mm joining length, the O.D. tapers inward from 2.6 mm to 2.1 mm. This leaves 4.5 cm of the inner catheter tube exposed at the distal end.

The proximal inner tube and the outer sheath circumferences are molded directly onto the distal end of the Luer-lock hub so that only the lumen of the inner tube can provide a channel for flow.

This design results in a single-channel catheter with two segments each having different degrees of flexibility: (i) a highly flexible 2.1 mm O.D., 4.5 cm long distal portion and (ii) a more rigid 2.6 mm O.D., 12.5 cm long proximal portion.

Graduation markings are at 5, 6, 7, 8, 9, and 10 cm distance from the extreme distal tip.

The two wings of the Luer-lock hub are positioned in line with the two catheter distal openings to enable orientation of these openings with tubal ostia while the catheter is in place within the uterine cavity.

INDICATIONS FOR USE

Intrauterine artificial insemination

CONTRAINDICATIONS

The Select IUI should not be used:

- in pregnant patients or in patients suspected of being pregnant
- in the presence of acute or chronic cervical infection
- in the presence of, or after recent Pelvic Inflammatory Disease (PID)

WARNINGS

Always wash the spermatozoa before performing intrauterine artificial insemination (including intracervical artificial insemination). The introduction of unwashed spermatozoa into

Instructions for Use (English)

Length: 17 cm (effective) 20 cm (overall)

O.D.: 2.6 mm stepped down to

2.1 mm for distal 4.5 cm of length

I.D.: 1.6 mm for the 17 cm effective length

Rx Only

CE2797

the uterine cavity with the Select IUI or by any other means will almost always result in severe adverse reactions which may include anaphylactic shock. Please refer to published medical literature for methods of preparing spermatozoa for intrauterine artificial insemination before performing this procedure.

Contents supplied sterile. Do not use if sterile barrier is damaged.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

PRECAUTIONS

Care should be taken prior to catheter insertion to ascertain the depth of the uterus and any uterine retroflexion or anteflexion present in individual patients. (This should be performed during a previous examination and not on the day of the artificial insemination when disturbance of the uterus should be minimized.) In no case should the device be forced against resistance. If passage of the Select IUI through the cervical canal cannot be achieved, the device should not be forced and the possibility of pathologic cervical stenosis should be considered.

ADVERSE REACTIONS

As with any instrument which is passed through the internal cervical os, mild cramping may be experienced by patients. In every case the INSTRUCTIONS FOR USE shall be followed; taking note of any CONTRAINDICATIONS, PRECAUTIONS and WARNINGS. ALWAYS USE WASHED SPERMATOZOA. REFER TO WARNINGS.

INSTRUCTIONS FOR USE

1. a. With the patient in the dorsal lithotomy position expose the uterine cervix to view with a vaginal speculum.
- b. Firmly attach a syringe (preferably 3 ml) with the plunger fully depressed, to the Luer-lock hub at the proximal end of the Select IUI.
- c. Place the distal tip of the Select IUI into the specimen of WASHED SPERMATOZOA so that the two lateral openings of the catheter are immersed in the specimen. Aspirate the desired quantity of the specimen (usually 0.5 ml) into the lumen of the catheter by retraction of the syringe plunger.

2. ALTERNATIVELY:

- a. A syringe filled with washed spermatozoa may be attached to the Luer-lock hub.
- b. Fill the catheter to the end with washed spermatozoa prior to inserting.
- c. Inseminate 0.5 ml of WASHED SPERMATOZOA.

3. Having previously determined the depth and direction of the patient's uterus, gently introduce the distal tip of the Select IUI into the external cervical os and advance it through the cervical canal into the uterine cavity until resistance is felt. The distal tip of the catheter should then be at the uterine fundus. Now retract the Select IUI by approximately 1 cm. The markings on the catheter provide guidance in estimating the length of catheter insertion and retraction.

NOTE: If any difficulty of insertion is encountered and bending of the catheter occurs, this can frequently be overcome by rotating the catheter gently while advancing it, or by grasping it close to the cervix with a ring forceps which is then used for advancement. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction of the cervix with a tenaculum may be required to perform full insertion. IN NO INSTANCE should force be used against resistance to achieve insertion.

4. With the distal tip of the Select IUI now positioned 1 cm from the uterine fundus, rotate the catheter to bring the Luer-lock hubwings to a horizontal position so as to orient the two lateral openings of the catheter within the uterine cavity. Then, SLOWLY INJECT THE WASHED SPERMATOZOA into the uterine cavity by depressing the syringe plunger fully.

NOTE: In no case should air be injected into the uterine cavity. Should there be any doubt about the complete expulsion of the specimen, detach the syringe from the Luer-lock hub; load the syringe with about 0.1 ml of sterile sodium chloride 0.9% (normal saline) solution; reattach the syringe to the Luer-lock hub; SLOWLY inject the sterile sodium chloride 0.9% solution to expel any remaining specimen into the uterine cavity.

5. Following injection of the specimen, slowly withdraw the Select IUI from the uterus and discard in accordance with Federal, State and local Medical / Hazardous waste practices. Remove the vaginal speculum and allow the patient to rest in the supine position for at least 10 minutes following the insemination.

REFERENCES

- Allen N.C., Herbert C.M., Maxson W.S., et al. Intrauterine insemination: a critical review. Fertility and Sterility 1985;44:5:569-580
- Martinez A.R., Bernardu R.E., Vermeiden J.P.W., Schoemaker J. Basic Questions on Intrauterine Insemination: An Update. Obstetrical & Gynecological Survey 1993;48:12:811-828
- Sertar R.T., Rufo S., Seibel M.M. Minimal Endometriosis and Intrauterine Insemination: Does Controlled Ovarian Hyperstimulation Improve Pregnancy Rates? Obstetrics & Gynecology 1992;80:1:37-40
- Stone S.C. Complications and Pitfalls of Artificial Insemination. Clin Obstet Gynecol 1980;23:4:667-682
- The American Fertility Society. Guideline for Practice: Intrauterine Insemination. 1991

EXPLANATION OF SYMBOLS

[REF] Reorder number **[STERILE EO]** Sterilized using ethylene oxide

[LOT] Batch code

Use-by date

Caution

Rx Only CAUTION: U.S.

Federal law restricts this device to sale by or on the order of a physician.

[EC REP] Authorized Representative in the European Community.

Not made with natural rubber latex

Manufacturer

Do not use if the package is damaged or opened

Consult instructions for use

Do not re-use

Do not re-sterilize

Select IUI® is a registered trademark of CooperSurgical, Inc.

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REFERENZEN

- Allen N.C., Herbert C.M., Maxson W.S., et al. Intrauterine insemination: a critical review. *Fertility and Sterility* 1985;44:5:569-580
- Martinez A.R., Bernardu R.E., Vermeiden J.P.W., Schoemaker J. Basic Questions on Intrauterine Insemination: An Update. *Obstetrical & Gynecological Survey* 1993;48:12:811-828
- Serta R.T., Rufo S., Seibel M.M. Minimal Endometriosis and Intrauterine Insemination: Does Controlled Ovarian Hyperstimulation Improve Pregnancy Rates? *Obstetrics & Gynecology* 1992;80:1:37-40
- Stone S.C. Complications and Pitfalls of Artificial Insemination. *Clin Obstet Gynecol* 1980;23:4:667-682
- The American Fertility Society. Guideline for Practice: Intrauterine Insemination. 1991

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© 2021 CooperSurgical, Inc.
Select IUI® è un marchio depositato di CooperSurgical, Inc.