TRAKE-fit™ Tracheostomy Tube Holder • Instructions for Use (English)
43-1850 Double Collar-Wide Tabs • 43-1851 Single Collar-Wide Tabs • 43-1862 Double Collar-Narrow Tabs

WARNINGS

- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.
- Unrelieved pressure may cause injury to the underlying tissues.
- Patient should be constantly monitored while this product is in use.

PLEASE NOTE: If a serious incident is suspected from using the TRAKE-fit, report the details of the incident to CooperSurgical via phone number +1 203-601-5200 Ext 3100 or by email at ProductSurveillance@coopersurgical.com and to the local Health Authority in your country. A serious incident may have caused or contributed to a death, a delay in a procedure which resulted in death or serious injury, or a malfunction that could have caused an adverse event.

INTENDED USERS

The TRAKE-fit is intended to be used by medical professionals only (including physicians, respiratory therapists, nurses, and emergency medical staff); they may be used in both emergency situations or under controlled conditions where quick and simple stabilization and support are necessary.


PATIENT POPULATION

The TRAKE-fit is intended to be used on patients in both emergency situations or under controlled conditions where quick and simple stabilization and support are necessary, including but not limited to the hospital, ER, ICU, or point of care. The specific patient population is relative to the product in use; the TRAKE-fit with narrow tabs will accommodate pediatric patients and the wide-tabbed TRAKE-fit should be used with larger patients and adults.

INSTRUCTIONS FOR USE

TRAKE-fit is supplied non-sterile

1. Place the neckband around the patient's neck with the hook tabs facing out. Remove the tabs from the neckband.
2. Insert the tabs, from the outside, through the slots in the flange of the tracheostomy tube so that the hook (rough) sides face the neckband.
3. Attach the hook tabs at a location on the neckband that is comfortable for the patient. Excess length can be cut off. Gauze pads can be used in the usual way.
4. For product 43-1850 and 43-1862, when the neckband becomes soiled, replace it with the extra one from the package.

GLOSSARY OF SYMBOLS
TRAKE-fit™ is a trademark of CooperSurgical, Inc.
CooperSurgical is a registered trademark of CooperSurgical, Inc.

Ref #	Symbol	Title	Description	Standard Development Organization
5.1.6		Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1
5.1.5		Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1
5.1.4		Date of manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1
5.1.3		Use-by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1
5.4.5, Annex B.2		Not made with natural rubber latex	Indicates that natural rubber latex was not used in the manufacturing of the product, its container, or its packaging.	ISO 15223-1
5.4.2		Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1
5.4.4		Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	ISO 15223-1
5.4.3		Consult instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1
5.2.7		Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.	ISO 15223-1
5.1.2		Authorized representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1
5.1.1		Manufacturer	Indicates the medical device manufacturer, as defined in EU Directive 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1
n.a.	RxOnly	Prescription device	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.	21 CFR 801.109
n.a.		Product conforms to the Medical Device Directive 93/42/EEC	Signifies European technical conformity.	n.a.

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© 2021 CooperSurgical, Inc.
GLOSSARIO DEI SIMBOLI TRAKE-fit™ è un marchio di CooperSurgical, Inc.
CooperSurgical è un marchio registrato di CooperSurgical, Inc.
CooperSurgical Distribution B.V.