Tampa Catheter • 5 French • 33 cm long • for Hysterosonography

61-2005

Instructions for Use (English)

CE2797

**CAUTION:** The Tampa Catheter Set is not intended to be used for aspiration of oocytes, in-vitro fertilization, gamete or zygote intrafallopian tube transfer, or any intratubal procedures.

**CAUTION:** U.S. Federal law restricts this device to sale by or on the order of a physician.

**WARNINGS:**

- Contents supplied sterile. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

Not made with natural rubber latex.

**CONTRAINDICATIONS:** Suspected infection, suspected pregnancy.

**INSTRUCTIONS FOR USE**

1. Prepare and position the patient in a manner consistent with the procedure for hysterosonography.
2. Fill the catheter with saline to expel air and adjust so that the distal end of the catheter protrudes one-half to one centimeter beyond the distal end of the sheath.
3. With the aid of a speculum, place the catheter and sheath so that the distal rounded end of the sheath gently contacts the external os of the cervix.
4. Hold the sheath stationary and advance the catheter through the cervical canal and well into the uterine cavity. **IMPORTANT: The tip of the catheter should be well into the uterus with the band at or near the external os of the cervix.**
5. Carefully remove the speculum and insert the transvaginal transducer into the proper position adjacent to the catheter.
6. After activation of ultrasound equipment, inject saline as required for the procedure. Typically, 3 to 20 cc (ml) of saline is sufficient for proper diagnosis.
7. Complete the study in a routine manner.

**EXPLANATION OF SYMBOLS**

|                                                                                                                 |                                                                                                       |                                                                                                                        |                                                                                                 |
|-----------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
|  Reorder Number                |  Do not re-use       |  Consult instructions for use         | <b>RxOnly</b> U.S. Federal law restricts this device to sale by or on the order of a physician. |
|  Batch Code                    |  Do not resterilize |  Do not use if package is damaged    | <b>EC REP</b> Authorized Representative in the European Community                               |
|  Use-by date                  |  Caution           |  Not made with natural rubber latex |                                                                                                 |
|  Sterilized using Irradiation |                                                                                                       |                                                                                                                        |                                                                                                 |
|  Manufacturer                |                                                                                                       |                                                                                                                        |                                                                                                 |

CooperSurgical is a registered trademark of CooperSurgical, Inc.

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© 2020 CooperSurgical, Inc.
CooperSurgical è un marchio depositato di CooperSurgical, Inc.