
Model 7033
Uni-Sem™ Intrauterine Catheter
for Intrauterine Artificial Insemination
Instructions for Use (English)

CE₂₇₉₇

DEVICE DESCRIPTION

The Uni-Sem™ intrauterine catheter is designed for use in intrauterine artificial insemination. It is a single-use, sterile, disposable, flexible step-down catheter with two, opposing, lateral ports, each located 2.0 mm from its rounded, smooth and closed distal tip.

The device has an overall length of 20.6 cm including a 3.5 cm long, polypropylene Luer lock adapter that is permanently affixed at its proximal end. The insertion portion of the device is comprised of a 17.1 cm long, clear polyethylene tube with a 2.3 mm OD (Outer Diameter) and a 1.7 mm ID (Inner Diameter) that is uniform throughout its length and is capable of containing 0.4 cc (ml) of fluid.

The proximal 12.6 cm of this tube is tightly surrounded by an outer, clear polypropylene sheath with a 2.6 mm OD for all but its distal 0.1 cm of length within which it gradually tapers to accomplish a smooth transition with the 2.3 mm OD of the inner tube; thereby leaving the distal 4.5 cm of the inner tube exposed. The proximal circumferences of both the inner tube and the outer sheath are molded directly into the distal end of the Luer lock adapter so that only the lumen of the inner tube can provide a channel for flow.

This design results in a single-channel catheter with two segments each having different degrees of flexibility:

- 1) a highly flexible, 2.3 mm OD, 4.5 cm long distal segment; and
- 2) a more rigid, 2.6 mm OD, 12.6 cm long proximal segment.

Colored graduation markings are located at 5.5 cm and 6.5 cm distance from the extreme distal tip of the catheter. Two wings on the luer lock adapter are positioned on the same plane as the catheter's two distal lateral ports to enable orientation of those ports within the uterus when they cannot be seen.

WARNINGS

- Always use washed spermatozoa when performing intrauterine (including intracervical) artificial insemination.
- The introduction of unwashed spermatozoa into the uterus with the Uni-Sem Intrauterine Catheter or by whatever other means will almost always result in severe adverse reaction which may include anaphylactic shock.
- Contents supplied sterile. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

CAUTION

U.S. Federal law restricts this device to sale by or on the order of a physician.

INDICATIONS FOR USE

Intrauterine Artificial Insemination

CONTRAINDICATIONS

- **NOT FOR USE IN INTRA-FALLOPIAN PROCEDURES.**
- The catheter should not be used in the presence of chronic cervical infection.
- The catheter should not be used in the presence of or after recent Pelvic Inflammatory Disease.

ADVERSE REACTIONS

- No adverse reactions to the use of the catheter have been reported. HOWEVER, as with any instrument which is passed through the internal cervical os, mild cramping may be expected.
- In every case, the INSTRUCTIONS FOR USE should be followed; taking note of any CONTRAINDICATION, PRECAUTIONS and WARNINGS.
ALWAYS USE WASHED SPERMATOZOA - see **WARNINGS**.

PRECAUTIONS

- By virtue of its flexibility and small diameter, the catheter can be introduced into the uterus with minimal danger of uterine wall perforation. However, care should be taken prior to its insertion to ascertain the depth of the uterus and any uterine retroflexion or antelexion present in individual patients. (This is preferably done during a previous examination and not on the day of artificial insemination when disturbance of the uterus should be minimized.) In NO CASE should the instrument be forced against resistance.
- If passage of the catheter through the cervical canal is not found to be easily achievable, the instrument should not be forced and the possibility of the presence of pathologic cervical stenosis considered.

INSTRUCTIONS FOR USE

- 1a. With the patient in the dorsal lithotomy position, expose the uterine cervix to view with a vaginal speculum.
b. Firmly attach a syringe (preferable 5 cc [ml] or 3 cc [ml]), with the plunger fully depressed, to the Luer lock adapter at the proximal end of the Uni-Sem Intrauterine Catheter.
c. Draw 1.0 cc (ml) of air into the syringe. Then place the distal tip of the catheter into the specimen of WASHED SPERMATOZOA so that both the lateral ports are immersed in the specimen. Now aspirate the desired quantity of the specimen into the lumen of the catheter by retraction of the syringe plunger. (The clear tubular portion of the catheter has the capacity of approximately 0.4 cc [ml] of fluid.)
2. **Alternatively**, a prewashed syringe of spermatozoa can be attached to the Luer lock of the Uni-Sem Intrauterine Catheter. The washed sperm should then be used to fill the entire catheter prior to insemination.

3. Having previously determined the depth and direction of the patient's uterus, gently introduce the distal tip of the catheter into the external cervical os and advance it through the cervical canal into the uterine cavity until resistance is felt. The distal tip of the catheter should then be at the uterine fundus. Now retract the instrument approximately 1.0 cm. (At the point that the narrower distal segment of the catheter steps up in size from 2.3 mm to 2.6 mm a total of 4.5 cm of the catheter's length will have been inserted. The markings at 5.5 cm and 6.5 cm on the catheter tube provide a further guide in judging the depth of insertion and retraction.) **NOTE:** The small OD of the distal segment of the catheter and its more rigid proximal segment usually permits its easy insertion. If any difficulty of insertion is encountered and bending of the catheter occurs, this frequently can be overcome by rotating the catheter gently while advancing it or by grasping it close to the cervix with a ring forceps which is then used for advancement. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to accomplish full insertion. In no event should force be used against resistance to achieve insertion.
4. With the distal tip of the catheter now positioned 1.0 cm from the uterine fundus, rotate the catheter to bring the wings on the Luer lock adapter to a horizontal position parallel to the examining table. This orients the two lateral ports at the distal end of the catheter within the uterine cavity toward the uterine ostium.
5. Now, SLOWLY INJECT NO MORE THAN 0.5 cc (ml) OF THE WASHED SPERMATOZOA into the uterine cavity.
6. Following placement of the specimen, slowly withdraw the Intrauterine Catheter from the uterus. Remove the speculum and allow the patient to rest.
7. Dispose of the catheter in accordance with all Federal, State and local Medical/Hazardous waste practices.

EXPLANATION OF SYMBOLS

REF Reorder number

LOT Batch code

 Use-by date

 Caution

 Do not re-use

 Do not resterilize

 Consult instructions for use

 Do not use if package is damaged

 Not made with natural rubber latex

STERILE EO Sterilized using ethylene oxide

R_x Only CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician

EC REP Authorized Representative in the European Community

 Manufacturer

Uni-Sem™ is a trademark of CooperSurgical, Inc.
CooperSurgical is a registered trademark of CooperSurgical, Inc.

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Uni-Sem™ è un marchio di fabbrica di CooperSurgical, Inc.
CooperSurgical è un marchio depositato di CooperSurgical, Inc.