# MX120 & MX124 • Uterine Explora™ Curette • Model I

## Curette with Sharp Randall-Type Cutting Edge and Vacu-Lok Syringe

### Instructions for Use (English)

CE 2797

#### DEVICE DESCRIPTION

The MX120 (3 mm) and the MX124 (4 mm) curettes are used to obtain clearly differentiated endometrial tissue, usually without anesthesia, in an office setting.

The flexibility of the nylon curette, the small diameter, and the sharp Randall-Type cutting edge allow better sampling from all areas of the endometrial cavity, including the cornual areas. The rounded end of this curette minimizes the possibility of perforation.

#### WARNINGS

- Contents supplied sterile. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.

#### CAUTION

U.S. Federal law restricts this device to sale by or on the order of a physician.

#### PRECAUTIONS

- Sounding of the uterus and the ensuing curettage should be performed with due care to avoid perforation of the uterine wall.
- NEVER use force to introduce the curette into the cervical canal.
- In cases where there is a marked cervical stenosis, Dilateria™ (*Laminaria japonica*) should be considered for use (to soften and dilate the cervical os) prior to performing the procedure.
- Any cervical manipulation may cause a vaso-vagal reaction. Patient should be watched for evidence of unusual pallor, nausea, vertigo or weakness. These symptoms will generally respond to about 3 to 10 minutes of rest and/or mild analgesic.
- In rare cases there may be spotting and mild cramping after the procedure has been performed. The patient should be instructed to contact her physician if bleeding persists or if a low grade fever is present.
- DO NOT DEPRESS SYRINGE PISTON WHILE CANNULA-CURETTE IS IN THE UTERUS.

#### INDICATIONS FOR USE

- Detection of endometrial hyperplasia
- Detection of endometrial cancer
- Initial endometrial biopsy prior to prescribing estrogen replacement therapy to ensure endometrial lining is normal
- Periodic monitoring of endometrial tissue in women receiving estrogen replacement therapy
- Endometrial dating
- Examination of glandular epithelium to determine ovulation and endometrial response to hormonal excretion and therapy
- Detection of pathology resulting in infertility
- Secondary Amenorrhea
- Monitor patient receiving TAMOXIFEN therapy for breast cancer

#### CONTRAINDICATIONS

- Patients where pregnancy is suspected
- Patients with, or recently recovered from, pelvic inflammatory disease
- Patients with any cervical or pelvic infections (*infections should be treated and cured before performing any endometrial sampling procedure*)
- Patients suffering from any diseases or conditions which could under all circumstances contraindicate outpatient surgical treatment – e.g. severe anemia, heart disease, or clotting mechanism deficiencies


#### INSTRUCTIONS FOR USE

1. Prepare vaginal area and cervix as you would for any sterile intrauterine procedure.
2. With a vaginal speculum in place and the cervical os visible, carefully sound the uterus with the sound curved in the direction of the canal, to determine both the position and the depth of the uterine cavity. In only an extremely small percentage of cases is anesthesia and/or premedication necessary. Advise patient that there may be a minimum amount of discomfort when the curette is first introduced.
3. It may be necessary to use a tenaculum to grasp the cervix and apply gentle traction to straighten the endocervical canal. To stabilize the cervix in those women with an anteverted uterus, grasp the anterior lip of the cervix with a tenaculum; if the uterus is retroverted, grasp the posterior lip of the cervix and apply gentle traction to straighten the cervical curvature.
4. After determining the depth of the uterine cavity, adjust the yellow sliding collar on the curette.

**NOTE:** Cutting edge of Explora Curette is aligned with the imprinted graduated markings and the index point (flattened surface at the base of curette).


# Uterine Explora™ Curette • Model I - Instructions for Use • English (Continued)

- Gently insert curette up to the yellow collar. (If cervical canal is very dry, apply a small amount of sterile water-soluble gel to the distal end of the curette.) Remove the stylet. Release traction on the tenaculum.
- Remove cap on syringe. Depress plunger to the bottom of the syringe barrel (this is the starting position of the syringe when attached to the curette). Attach syringe to curette by twisting.
- Withdraw plunger to desired level of negative pressure. Be careful not to twist plunger until ready to lock into position. Hold plunger as illustrated in **Figure A** to maintain suction. Use the following values as a guideline:


- 6 cc (ml) gives you 15 in. (38.1 cm) H<sub>2</sub>O negative pressure (approximate)
- 9 cc (ml) gives you 18 in. (45.7 cm) H<sub>2</sub>O negative pressure (approximate)
- 12 cc (ml) gives you 20 in. (50.8 cm) H<sub>2</sub>O negative pressure (approximate)

- To lock plunger in place:** grasp plunger close to barrel, as shown in **Figure B**. Twist plunger clockwise (to the right) 1/4 of a turn until you feel it catch. This locks plunger in place to maintain negative pressure throughout the procedure (test to be sure it is secure). You now have one-hand control of curette/syringe.


- Use long slow strokes to sample the four quadrants of the endometrium. For endometrial dating, one or two swipes are usually sufficient.
- Should the syringe fill with tissue and/or blood and you need additional negative pressure, disconnect syringe from curette, leaving the curette in the uterine cavity. Unlock syringe by twisting plunger to the left (counter-clockwise) 1/4 turn. Empty contents of syringe into collection jar containing fixative. Now reattach syringe to curette as directed in Steps 6 and 7.

**NOTE:** This avoids the need to reintroduce this or a second Explora™ into the same patient.

**WARNING:** Never introduce an instrument that has been exposed to fixative into the uterus.

- When you have completed the curettage, withdraw the curette from uterus. Hold over container of fixative. To expel tissue, twist plunger to the left (counter-clockwise) to unlock plunger. Expel tissue from curette into container of formalin or other acceptable fixative. Measure distance from yellow collar to distal (entering) end of curette. Enter this measurement into patient's chart.
- After the procedure is completed, cut distal end of curette (approximately 1 in. [2.5 cm]; see **Figure C**) and drop into container of fixative to be included with tissue expressed. See **Figure D**. (There may be residual tissue that pathologists would like to have available in this portion of the curette.)
- Withdraw fixative into curette and syringe and flush into specimen container. This will wash through any tissue that may have lodged in syringe or curette.
- Allow patient to remain in prone position on examining table for about 10-15 minutes.

**NOTE: Fill specimen container with fixative. Endometrial tissue is very hygroscopic and should be placed in fixative at least 10 times its volume to maintain tissue in satisfactory condition for pathologist.**

Inadequate tissue collection may occur in the presence of large uterine polyps or fibroids, an enlarged uterine cavity and/or an atrophic endometrium.

## EXPLANATION OF SYMBOLS

**REF** Reorder Number

**LOT** Batch Code

Use-by date

Caution

Do not use if package is damaged

Consult Instructions for Use

Do not re-use

**R<sub>x</sub> Only** **CAUTION:** U.S. Federal law restricts this device to sale by or on the order of a physician.

**EC/REP** Authorized Representative in the European Community.

Not made with natural rubber latex.

**STERILE EO** Sterilized using Ethylene Oxide

Manufacturer

Marks bearing ™ are trademarks of CooperSurgical, Inc.

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© 2022 CooperSurgical, Inc.
Cooper Surgical, Inc.