wallace®

Amniocentesis Needles

Rx Only

€ 2797

en

INSTRUCTIONS FOR USE

fr

MODE D'EMPLOI

de

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it

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es

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pt

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nl

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sv

BRUKSANVISNING

da

BRUGSANVISNING

no

BRUKSANVISNING

fi

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el

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2514 AP The Hague

The Netherlands

CooperSurgical, Inc.
95 Corporate Drive
Trumbull, CT 06611 USA

Wallace® Amniocentesis Needles

INSTRUCTIONS FOR USE

These instructions apply to the use of the following Wallace® products:

9022 Wallace® Amniocentesis Needles 22G x 90mm

12021 Wallace® Amniocentesis Needles 21G x 120mm

15021 Wallace® Amniocentesis Needles 21G x 150mm

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

DESCRIPTION

Wallace® Amniocentesis Needles (product codes 12021, 15021 and 9022) are sterile, single-use, devices provided sterile for the removal of amniotic fluid, trans-abdominally under ultrasound guidance. The devices consist of a stainless steel, 21 Gauge needle of differing lengths (120mm and 150mm) and a 22 Gauge 90mm long needle, all with self-locating stylets. The distal tip is engraved for 1cm of its length to aid location under trans-abdominal ultrasound. The proximal end of the needle is fitted with a Luer connection to ensure easy and safe connections.

INDICATIONS FOR USE

Amniocentesis is indicated for

- karyotyping to exclude Down's syndrome and other chromosomal abnormalities
- determining fetal sex in sex-linked disorders
- detection of numerous inherited disorders of metabolism
- assessment of severity of rhesus iso-immunisation by determining the optical density of the liquor, as an alternative to cordocentesis
- determining fetal lung maturity prior to delivery by measuring the lecithin/sphingomyelin ratio
- detection of raised alpha-fetoprotein in suspected cases of neural tube defect

The procedure is most commonly performed in women aged 35 years and over, or women who are felt to be high-risk for fetal abnormality. Amniocentesis should only be performed around 16 weeks of gestation^[1] and a detailed ultrasound study of the gestational sac and uterus is essential to confirm the gestational age and check the position of the uterus and placenta.

CONTRAINDICATIONS

Amniocentesis should not be performed in patients with:

- active vaginal bleeding
- coagulation disorders

This product should not be used for early amniocentesis (9-15 weeks) which may carry additional risk^[2].

INSTRUCTIONS FOR USE

- 1 Genetic counselling must be provided to all patients prior to the procedure.
- 2 Ensure the patient has an empty bladder before embarking on the procedure.
- 3 Perform a detailed ultrasound examination to visualise the fetus and the placenta, and to identify the most suitable liquor pool for sampling.
- 4 Prepare the entry site using antiseptic techniques.
- 5 If local anaesthetic is to be used, raise a bleb in the skin with a fine needle. Inject further anaesthetic through the abdominal and uterine wall (1% lignocaine may be used if no other contraindications exist).
- 6 Remove the needle from its packaging and discard the protective outer sheath.
- 7 With the stylet locked in position insert the needle, through needle guide if applicable, into the amniotic cavity under ultrasound guidance. The echomarking on the tip of the needle aids accurate placement in the liquor pool. Slowly remove the stylet whilst holding the needle firmly in place, and attach a syringe to the Luer fitting of the needle.
- 8 Withdraw the amniotic fluid, it may be advisable to discard the first 1ml to avoid contamination with maternal cells if karyotyping is to be undertaken. The volume of fluid withdrawn will depend upon the length of gestation, with 15-20ml recommended from the 16th week of pregnancy.
- 9 Withdraw the needle and dispose in accordance with local medical hazardous waste practices.
- 10 Expel the sample into a sterile, universal container labeled with the correct patient's details. Send sample to that laboratory as soon as possible.

PRECAUTIONS/CAUTION

- To minimize the risk of infection, aseptic procedures should be adhered to.
- Amniocentesis in multiple pregnancies requires special consideration and carries increased risk of abnormal results and spontaneous abortion^[3].
- The patient should be given thorough genetic counselling and be encouraged to ask questions prior to consenting to the procedure.
- To confirm gestational age, fetal viability, the site of the placenta and to exclude multiple pregnancy, an ultrasound examination should always precede any attempts at sampling. Fetal viability and age should be confirmed prior to detailed counselling and scheduling of amniocentesis.
- In the rhesus-negative, unsensitized patient anti-D immunoglobulin should be administered after the procedure.
- Patients should be advised to rest for 48 hours following the procedure.
- To be used by, or under the direction of a qualified person familiar with the technique of Amniocentesis.
- Discard if product or packaging damaged.

WARNINGS:

- Placental trauma penetration has been reported to increase amniocentesis-related loss rates^[4]. Penetration of the placenta and damage to the umbilical cord should be avoided.
- A new needle should be used for each sampling attempt. If a suitable sample has not been obtained after 2 consecutive attempts the procedure should be abandoned.

PATIENT INFORMATION

Patients should be informed of the risks associated with the procedure which should include:

- failure to withdraw amniotic fluid
- failure to culture the fluid after sampling
- an artifact in results of the test
- a spontaneous abortion rate of 0.3-1.5% within 4 weeks of the procedure^[5] and that this risk may increase with maternal age
- an increased risk of premature labour
- an increased risk *abruptio placentae*
- suggestions of an increased incidence of respiratory distress syndrome and pneumonia in the newborn^[5].

Patients should be made aware of possible post-amniocentesis symptoms which include a 3% chance of amniotic fluid leak^[6], vaginal spotting and cramping.

Informed consent should be obtained prior to the procedure and documented in the patient's notes.

STERILIZATION

This product is sterilized by ethylene oxide.

REFERENCES

- 1 Hanson FW et al. Ultrasonography-guided early amniocentesis in singleton pregnancies. *Am J Obstet Gynecol* 1990, **162**: 1376-83.
- 2 Smidt-Jensen S and Sundberg K. Early Amniocentesis. *Curr Opinion in Obs Gynaecol* 1995, **7**:117-121
- 3 Palle C et al. Increased risk of abortion after genetic amniocentesis in twin pregnancies. *Prenat Diagn* 1983; 83-89.
- 4 Hess LW et al. Significance of opaque discoloured amniotic fluid at second trimester amniocentesis. *Obstet Gynecol* 1986, **67**: 44.
- 5 Tabor A et al. Randomised controlled trial of genetic amniocentesis in 4606 low-risk women. *Lancet* 1986, **I**: 1287-1293.
- 6 Smidt-Jensen S et al. Randomised comparison of amniocenteses and transabdominal and transcervical chorionic villi sampling. *Lancet* 1992, **340**:1237-1244.

Wallace® is a registered trademark of CooperSurgical, Inc. and its subsidiaries.

CooperSurgical is a registered trademark of CooperSurgical, Inc.

© 2021 CooperSurgical, Inc. All rights reserved.

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RÉFÉRENCES

- 1** Hanson FW *et al.* Ultrasonography-guided early amniocentesis in singleton pregnancies. *Am J Obstet Gynecol* 1990, **162**: 1376-83.
- 2** Smidt-Jensen S and Sundberg K. Early Amniocentesis. *Curr Opinion in Obs Gynae* 1995, **7**:117-121
- 3** Palle C *et al.* Increased risk of abortion after genetic amniocentesis in twin pregnancies. *Prenat Diagn* 1983: 83-89.
- 4** Hess LW *et al.* Significance of opaque discoloured amniotic fluid at second trimester amniocentesis. *Obstet Gynecol* 1986, **67**: 44.
- 5** Tabor A *et al.* Randomised controlled trial of genetic amniocentesis in 4606 low-risk women. *Lancet* 1986, **I**: 1287-1293.
- 6** Smidt-Jensen S *et al.* Randomised comparison of amniocenteses and transabdominal and transcervical chorionic villi sampling. *Lancet* 1992, **340**:1237-1244.

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LITERATURHINWEISE

- 1 Hanson FW et al. Ultrasonography-guided early amniocentesis in singleton pregnancies. *Am J Obstet Gynecol* 1990, **162**: 1376-83.
- 2 Smidt-Jensen S and Sundberg K. Early Amniocentesis. *Curr Opinion in Obs Gynaecol* 1995, **7**:117-121
- 3 Palle C et al. Increased risk of abortion after genetic amniocentesis in twin pregnancies. *Prenat Diagn* 1983: 83-89.
- 4 Hess LW et al. Significance of opaque discoloured amniotic fluid at second trimester amniocentesis. *Obstet Gynecol* 1986, **67**: 44.
- 5 Tabor A et al. Randomised controlled trial of genetic amniocentesis in 4606 low-risk women. *Lancet* 1986, **I**: 1287-1293.
- 6 Smidt-Jensen S et al. Randomised comparison of amniocenteses and transabdominal and transcervical chorionic villi sampling. *Lancet* 1992, **340**:1237-1244.

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BIBLIOGRAFIA

- 1 Hanson FW *et al.* Ultrasonography-guided early amniocentesis in singleton pregnancies. *Am J Obstet Gynecol* 1990, **162**: 1376-83.
- 2 Smidt-Jensen S and Sundberg K. Early Amniocentesis. *Curr Opinion in Obs Gynae* 1995, **7**:117-121
- 3 Palle C *et al.* Increased risk of abortion after genetic amniocentesis in twin pregnancies. *Prenat Diagn* 1983: 83-89.
- 4 Hess LW *et al.* Significance of opaque discoloured amniotic fluid at second trimester amniocentesis. *Obstet Gynecol* 1986, **67**: 44.
- 5 Tabor A *et al.* Randomised controlled trial of genetic amniocentesis in 4606 low-risk women. *Lancet* 1986, **I**: 1287-1293.
- 6 Smidt-Jensen S *et al.* Randomised comparison of amniocenteses and transabdominal and transcervical chorionic villi sampling. *Lancet* 1992, **340**:1237-1244.

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REFERENCIAS

- 1 Hanson FW *et al.* Ultrasonography-guided early amniocentesis in singleton pregnancies. *Am J Obstet Gynecol* 1990, **162**: 1376-83.
- 2 Smidt-Jensen S and Sundberg K. Early Amniocentesis. *Curr Opinion in Obs Gynaecol* 1995, **7**:117-121
- 3 Palle C *et al.* Increased risk of abortion after genetic amniocentesis in twin pregnancies. *Prenat Diagn* 1983; 83-89.
- 4 Hess LW *et al.* Significance of opaque discoloured amniotic fluid at second trimester amniocentesis. *Obstet Gynecol* 1986, **67**: 44.
- 5 Tabor A *et al.* Randomised controlled trial of genetic amniocentesis in 4606 low-risk women. *Lancet* 1986, **I**: 1287-1293.
- 6 Smidt-Jensen S *et al.* Randomised comparison of amniocenteses and transabdominal and transcervical chorionic villi sampling. *Lancet* 1992, **340**:1237-1244.

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LITERATUUR

- 1 Hanson FW et al. Ultrasonography-guided early amniocentesis in singleton pregnancies. *Am J Obstet Gynecol* 1990, **162**: 1376-83.
- 2 Smidt-Jensen S and Sundberg K. Early Amniocentesis. *Curr Opin in Obs Gynae* 1995, **7**:117-121
- 3 Palle C et al. Increased risk of abortion after genetic amniocentesis in twin pregnancies. *Prenat Diagn* 1983: 83-89.
- 4 Hess LW et al. Significance of opaque discoloured amniotic fluid at second trimester amniocentesis. *Obstet Gynecol* 1986, **67**: 44.
- 5 Tabor A et al. Randomised controlled trial of genetic amniocentesis in 4606 low-risk women. *Lancet* 1986, **I**: 1287-1293.
- 6 Smidt-Jensen S et al. Randomised comparison of amniocenteses and transabdominal and transcervical chorionic villi sampling. *Lancet* 1992, **340**:1237-1244.

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- 1 Hanson FW et al. Ultrasonography-guided early amniocentesis in singleton pregnancies. *Am J Obstet Gynecol* 1990, **162**: 1376-83.
- 2 Smidt-Jensen S and Sundberg K. Early Amniocentesis. *Curr Opinion in Obs Gynae* 1995, **7**:117-121
- 3 Palle C et al. Increased risk of abortion after genetic amniocentesis in twin pregnancies. *Prenat Diagn* 1983: 83-89.
- 4 Hess LW et al. Significance of opaque discoloured amniotic fluid at second trimester amniocentesis. *Obstet Gynecol* 1986, **67**: 44.
- 5 Tabor A et al. Randomised controlled trial of genetic amniocentesis in 4606 low-risk women. *Lancet* 1986, **I**: 1287-1293.
- 6 Smidt-Jensen S et al. Randomised comparison of amniocenteses and transabdominal and transcervical chorionic villi sampling. *Lancet* 1992, **340**:1237-1244.

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REFERENSSIT

- 1 Hanson FW *et al.* Ultrasonography-guided early amniocentesis in singleton pregnancies. *Am J Obstet Gynecol* 1990, **162**: 1376-83.
- 2 Smidt-Jensen S and Sundberg K. Early Amniocentesis. *Curr Opinion in Obs Gynaecol* 1995, **7**:117-121
- 3 Palle C *et al.* Increased risk of abortion after genetic amniocentesis in twin pregnancies. *Prenat Diagn* 1983: 83-89.
- 4 Hess LW *et al.* Significance of opaque discoloured amniotic fluid at second trimester amniocentesis. *Obstet Gynecol* 1986, **67**: 44.
- 5 Tabor A *et al.* Randomised controlled trial of genetic amniocentesis in 4606 low-risk women. *Lancet* 1986, **I**: 1287-1293.
- 6 Smidt-Jensen S *et al.* Randomised comparison of amniocenteses and transabdominal and transcervical chorionic villi sampling. *Lancet* 1992, **340**:1237-1244.

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en**sv**

Single use. Do Not Reuse: Medical devices require specific material characteristics to perform as intended. These characteristics have been verified for single use only. Any attempt to re-process the device for subsequent re-use may adversely affect the integrity of the device or lead to deterioration in performance.

Sterile unless unit container is opened or damaged. Destroy after single use. Do not resterilize.

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