# SM507 • TheCurve®

## Intrauterine Insemination Catheter

For intrauterine artificial insemination

## Instructions for Use (English)

Length: 17.2 cm (effective) 19.9 cm (overall)

O.D.: 2.6 mm stepped down to

1.6 mm for distal 4.5 cm of length

I.D.: 1.1 mm for the 17.2 cm effective length

### WARNING: NEVER ATTEMPT TO RESTERILIZE.

MOUSE EMBRYO ASSAY (MEA) 2-CELL SURVIVAL >80%

**This device is not intended for use for intrafallopian tube procedures, in vitro fertilization (IVF) and gamete intrafallopian transfer (GIFT).**

**R<sub>x</sub> Only**  
**CE 2797**

### DEVICE DESCRIPTION

TheCurve® is a single-use, sterile, disposable, flexible catheter designed for use in intrauterine insemination. The opening is at the tip of the catheter which is rounded and beveled. The distal part of the catheter is curved to ease its insertion.

The device has an over all length of 19.9 cm including a 2.7 cm long polyethylene Luer-lock hub which is permanently affixed to its proximal end.

The effective portion of the catheter (intended for insertion) comprises a 17.2 cm long clear polyethylene tube with a 1.6 mm O.D. (outer diameter) and a 1.1 mm I.D. (inner diameter). When filled, TheCurve contains 0.2 ml of fluid. The proximal 12.7 cm length of the catheter is tightly surrounded by a clear polypropylene sheath with a 2.6 mm O.D. for all but the distal 1 mm length, which joins the inner tube to the outer sheath. At this 1 mm joining length, the O.D. tapers inward from 2.6 mm to 1.6 mm over the 1 mm spindle-shaped portion. This leaves 4.5 cm of the inner catheter tube exposed at the distal end.

The proximal inner tube and the outer sheath circumferences are molded directly onto the distal end of the Luer-lock hub so that only the lumen of the inner tube can provide a channel for flow. This design results in a single-channel catheter with two segments, each having different degrees of flexibility (i) a highly flexible 1.6 mm O.D., 4.5 cm long distal portion and (ii) a more rigid 2.6 mm O.D., 12.7 cm long proximal portion.

Graduation markings are placed at 5, 6, 7, 8, 9, and 10 cm distance from the extreme distal tip. The distal portion is curved to ease its passage into the cervical os and insertion through the cervical canal up to the uterine fundus.

### INDICATIONS FOR USE

Intrauterine artificial insemination

### CONTRAINDICATIONS

TheCurve should not be used in:

- pregnant patients or in patients suspected of being pregnant
- the presence of acute or chronic cervical infection
- the presence of, or after recent Pelvic Inflammatory Disease (PID)

### WARNINGS

The Instructions for Use are for guidance only. The clinical user has full responsibility of the technique to be used.

Always wash the spermatozoa before performing Intrauterine Artificial Insemination (including Intracervical Artificial Insemination). The direct introduction of unwashed spermatozoa into the uterine cavity will almost always result in severe adverse reactions which may include anaphylactic shock. Refer to published medical literature for methods of washing spermatozoa for Intrauterine Artificial Insemination before performing this procedure.

Contents supplied sterile. Do not use if sterile barrier is damaged.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

### PRECAUTIONS

Care should be taken prior to catheter insertion to ascertain the depth of the uterus and any uterine retroflexion or antelexion present in individual patients. (This should be performed during a previous examination and not on the day of the artificial insemination when disturbance of the uterus should be minimized.) In no case should the device be forced against resistance. If passage of TheCurve through the cervical canal cannot be achieved, the device should not be forced and the possibility of pathologic cervical stenosis should be considered.

### ADVERSE REACTIONS

As with any instrument which is passed through the internal cervical os, mild cramping may be experienced by patients. In every case the INSTRUCTIONS FOR USE should be followed; taking note of any CONTRAINDICATIONS, PRECAUTIONS and WARNINGS.

ALWAYS USE WASHED SPERMATOZOA.  
REFER TO WARNINGS.

## INSTRUCTIONS FOR USE

1. With the patient in the dorsal lithotomy position, expose the uterine cervix to view with a vaginal speculum.
2. According to the method used to wash and prepare the spermatozoa, EITHER firmly attach a Luer-lock fitted syringe (preferably 3 ml) containing the WASHED SPERMATOZOA to the hub at the proximal end of TheCurve and load the catheter with the washed spermatozoa OR firmly attach an empty 3 ml Luer-lock fitted syringe with the plunger fully depressed to the hub at the proximal end of TheCurve. Immerse the distal tip of TheCurve into the specimen of WASHED SPERMATOZOA and aspirate the desired quantity of the specimen (usually 0.5 ml) into the lumen by retraction of the syringe plunger.
3. Having previously determined the depth and direction of the patient's uterus, gently introduce the distal tip of the catheter, the curved tip at 12 o'clock, into the external cervical os and advance it through the cervical canal into the uterine cavity until resistance is felt. The curved distal tip should then be at the uterine fundus. Now retract the catheter by approximately 1 cm. The markings on the catheter tube provide guidance in estimating the length of catheter insertion and retraction.  
**NOTE:** If any difficulty of insertion occurs, this can frequently be overcome by grasping the catheter close to the cervix with a ring forceps which is then used for advancement. If the uterus is significantly flexed anteriorly or posteriorly, gentle traction on the cervix with a tenaculum may be required to perform full insertion. IN NO INSTANCE should force be used against resistance to achieve insertion.
4. With the distal curved tip of TheCurve now positioned 1 cm from the uterine fundus, SLOWLY INJECT THE WASHED SPERMATOZOA into the uterine cavity by depressing the syringe plunger fully. Should there be any doubt about the complete expulsion of the specimen in the uterus, detach the syringe from the Luer-lock catheter hub, load the syringe with about 0.1 ml of sterile sodium chloride 0.9% (normal saline) solution, reattach the syringe to the Luer-lock catheter hub and, SLOWLY inject the sterile sodium chloride 0.9% solution in the insemination catheter.
5. Following injection of the specimen, slowly withdraw TheCurve from the uterus and discard in accordance with Federal, State and local Medical / Hazardous waste practices. Remove the vaginal speculum and allow the patient to rest in the supine position for at least 10 minutes following the insemination.

## EXPLANATION OF SYMBOLS

 Reorder number	 Sterilized using ethylene oxide
 Batch code	<b>R<sub>x</sub>Only</b> <b>CAUTION:</b> U.S. Federal law restricts this device to sale by or on the order of a physician.
 Use-by date	
 Caution	 Authorized Representative in the European Community.
 Do not use if the package is damaged or opened	 Not made with natural rubber latex
 Consult instructions for use	 Manufacturer
 Do not re-use	
 Do not re-sterilize	

TheCurve® is a registered trademark of CooperSurgical, Inc.  
© 2021 CooperSurgical, Inc.

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TheCurve® è un marchio depositato di CooperSurgical, Inc.