
The Koh Colpotomizer™ System

Koh Cup™ Vaginal Fornices Delineator and Colpo-Pneumo Occluder™ (CPO-6)

Instructions for Use (English)

DESCRIPTION

Koh Cup™ Vaginal Fornices Delineator (KCS-Kit and KCP-Kit) is designed for use with The RUMI® System and Colpo-Pneumo Occluder™. The Koh Colpotomizer™ System is intended to be used in laparoscopic procedures where it is desirable to locate the position of the vaginal fornices. The Koh Cups are available in sizes 3.0 cm, 3.5 cm and 4.0 cm. The Koh Cup is available in Ultem® Polymeric Resin for use during electrosurgery. The stainless steel Koh Cup is available for use with laser or harmonic scalpel. These devices are non-sterile and must be sterilized according to standard procedure prior to use.

The Colpo-Pneumo Occluder (CPO-6) is a sterile single-use medical grade silicone device designed for use with The RUMI System and the Koh Cup Vaginal Fornices Delineator. The Koh Colpotomizer System is intended for use in laparoscopic procedures where it is desirable to minimize the loss of pneumoperitoneum after a colpotomy incision has been made.

WARNINGS

The Colpo-Pneumo Occluder is designed to accommodate no more than 180 cc of inflation. It is supplied as a sterile single-use device and should never be reused. The Koh Cup Vaginal Fornices Delineator is non-sterile and should be sterilized prior to use.

The Ultem Polymeric Resin Koh Cup is not intended to be used with harmonic or other types of ultrasonic energy sources. The stainless steel Koh Cup is not intended to be used with electrocautery. Use of electrocautery may result in inadvertent injury to adjacent structures.

Insure that placement of the Koh Cup Vaginal Fornices Delineator and Colpo-Pneumo Occluder on the RUMI Uterine Manipulator does not interfere or restrict the operation of the uterine manipulator in any way. The RUMI Uterine Manipulator with Koh Colpotomizer System devices attached should always be inserted into the vagina slowly and under direct vision. The cervix should be adequately dilated prior to insertion of the RUMI Uterine Manipulator Tip. Follow the RUMI Uterine Manipulator Instructions for Use for a detailed description of cervix dilation.

CAUTION

U.S. Federal law restricts this device to sale by or on the order of a physician.

INDICATIONS FOR USE

The Koh Colpotomizer System is indicated for use in laparoscopic procedures where use of a uterine manipulator is appropriate and the surgeon intends to remove or access intraperitoneal tissue through the vagina by use of a colpotomy or culdotomy incision.

CONTRAINDICATIONS

The Koh Colpotomizer System should not be used in patients who are pregnant or who are suspected of being pregnant, in patients that have an IUD in place, in patients planning gamete intrafallopian transfer, and in cases where the surgeon deems it inadvisable or finds it difficult to insert the silicone tip into the cervix or uterus.

ADVERSE REACTIONS

During surgery, the Koh Colpotomizer System with Koh Cup Vaginal Fornices Delineator and the Colpo-Pneumo Occluder become part of The RUMI System Uterine Manipulator. Adverse reactions which have been suspected or reported to be associated with the uterine manipulators are cramping, infection, uterine perforation and uterine spasm with accompanying temporary physiologic blockage of patent fallopian tubes. (The order of listing does not indicate frequency or severity.)

INSTRUCTIONS FOR USE

1. Prepare the patient for surgery according to standard surgical practices. Prepare the RUMI Uterine Manipulator for use in accordance with the RUMI Tip Instructions for Use.

2. Prior to patient insertion, slide the Colpo-Pneumo Occluder over the RUMI disposable tip to extend over the distal shaft of the uterine manipulator. (See Figures 1 and 2.)

3. Select the appropriate size Koh Cup by visually inspecting the diameter of the patient's cervix. Choose a Koh Cup with sufficient diameter to readily accept the cervix. If in question, it is recommended that the selected Koh Cup be tried on the cervix prior to attachment to the uterine manipulator. To do so, simply grasp the cup with an appropriate grasper and place it over the cervix to ensure that it fits correctly. (It is further recommended that a Koh Cup that is slightly too large be selected over one that is slightly too small.)

4. Slide the selected Koh Cup, small end first, over the RUMI disposable tip and onto the tip base until it is snugly attached to the tip base. The distal end of the Colpo-Pneumo Occluder will be between the Koh Cup and the tip base when properly installed. (See Figures 3 and 4.)

5. Test inflate the Colpo-Pneumo Occluder with 20-60 cc of sterile saline prior to patient insertion in order to verify proper operation of the occluding balloon and then deflate prior to patient insertion. Test inflate the RUMI Tip balloon according to the RUMI Instructions for Use.

6. With assembly complete, dilate the cervix to a sufficient diameter to easily accommodate the insertion of the uterine manipulator tip. Use a speculum or a combination of vaginal retractors to gain sufficient exposure of the cervix for device insertion. Slowly insert the manipulator device into the vagina and manipulator tip into the uterus under direct vision and in accordance with the RUMI Instructions for Use. Once insertion is complete and the Koh Cup is fully engaged on the cervix, the RUMI Uterine Manipulator may be operated as normal.

7. Prior to making colpotomy incisions, inflate the Colpo-Pneumo Occluder balloon with sufficient sterile saline or water (do not use air) to cause the inflated balloon to make continuous radial contact with the vaginal wall. 90 cc to 120 cc of inflation will generally be sufficient. With the Koh Colpotomizer System attached to the RUMI Uterine Manipulator and properly inserted in the patient, inflation of the Colpo-Pneumo Occluder will not impede movement of the RUMI's articulating tip and will prevent leakage of the pneumoperitoneum from the vagina during and after any colpotomy incisions for as long as the Colpo-Pneumo Occluder is inflated.

8. Laparoscopic identification and delineation of the vaginal fornices is aided by using the uterine manipulator to articulate the uterus into a position opposite the point of desired delineation while applying cephalad pressure to the uterine manipulator handle. This will increase tissue traction in the desired area of delineation so that the location of the rim of the Koh Cup which marks the position of the vaginal fornices) may easily be determined.

9. After dissection is complete, deflate the Colpo-Pneumo Occluder only. Do not deflate the uterine manipulator tip balloon.

10. Extend a tenaculum into the vagina through one of the windows on the Koh Cup and grasp the cervix.

11. Using the manipulator handle and tenaculum, pull the uterus into the vagina.

12. Removal of the uterus and the Koh Colpotomizer System:
 - a. For vaginal closure, continue to use the tenaculum and manipulator to pull the uterus through the vaginal canal until it is removed from the patient. The Koh Colpotomizer System will be removed during this process along with the uterus and uterine manipulator.
 - b. For laparoscopic closure, with the Colpo-Pneumo Occluder deflated, leave the uterus lodged in the vaginal canal so that it may serve to occlude the vagina. Once the pneumoperitoneum is reestablished, closure may proceed in the usual manner. The Koh Colpotomizer System will be removed along with the uterus and uterine manipulator during this process.
 - c. Check carefully to be sure the system is intact and no parts remain in the patient.

13. After surgery, discard the disposable components of the Koh Colpotomizer System according to standard hospital practices.

REPROCESSING INSTRUCTIONS KOH CUPS

1. INTRODUCTION

This section is intended to provide detailed instructions for effectively processing reusable Koh Cups. All reusable instruments should be thoroughly cleaned and sterilized to prepare them for use.

Instruments which can be disassembled should be disassembled for cleaning to ensure maximum contact between surfaces and cleaning tools and solutions (see specific instructions for each instrument). CooperSurgical, Inc. has validated the processes provided in these instructions to be capable of effective instrument processing. Equipment, operators, cleaning agents and procedures all contribute to the efficacy of processing. Healthcare facilities should ensure that selected processing steps are safe and effective within their systems.

Alternative methods of processing these instruments outside the instructions described in this document may be suitable for reprocessing; however, these have not been evaluated by CooperSurgical, Inc. Operators and healthcare facilities which choose to perform processes outside the instructions described in this document must validate these processes before use. In the event that national or regional government requirements conflict with the recommendations provided here, these shall override the recommendations of CooperSurgical, Inc., but must be validated before use. Reusable instruments must be rinsed thoroughly after use and after cleaning to ensure removal of gross soil and residual cleaners and solvents. Gross soil is damaging to instrument surfaces and inhibits thorough cleaning and subsequent sterilization. Residual cleaners and solvents can impede further processing and may affect the surface finish of the instrument over time.

CLEANING

There are two methods for cleaning that have been validated by CooperSurgical, Inc. Any healthcare facility should be able to perform the manual cleaning process. For those facilities that use automated washer-disinfectors, an automated method has also been validated.

Personnel are cautioned to use personal protective apparel due to the unknown and potentially hazardous nature of biological fluids and soils present. Specific detergents or cleaning solutions vary in requirements for concentrations and temperatures required for optimal cleaning performance. Use temperatures in this document as guidance, to be superseded by the detergent/cleaner manufacturer's instructions for use. Any processes outside the recommendations in this document must be validated before use.

The quality of water used to prepare solutions and rinse reusable instruments can affect the efficacy of the process and the instrument reuse life. CooperSurgical, Inc. recommends the use of freshly prepared purified water or sterile water for dilutions of solution and rinsing of instruments. Non-purified water can add mineral deposits and recontaminate instruments with microorganisms. Mineral deposits can impede sterilization and affect the condition of the instrument, resulting in staining, corrosion, and/or other damage.

Equipment and Materials Required for Manual Cleaning

- Personal protective apparel/equipment as recommended by cleaning solution supplier.
- Cleaning bath, sink, or other vessel large enough to accommodate full immersion of instruments.
- Freshly prepared cleaning solution intended for manual cleaning (use enzymatic, neutral, or alkaline detergents).
- Soft-bristled brushes and sponges.
- Clean, low-lint or lint-free cloth.

Equipment and Materials Required for Automated Cleaning

- Personal protective apparel/equipment as recommended by cleaning solution supplier.
- Legally marketed (FDA-cleared and/or CE-marked) medical device washer or washer-disinfector.
- Cleaning solution intended for automated cleaning (use enzymatic, or neutral, alkaline detergents).
- Racks and accessories to hold and support instruments during processing.

POINT OF USE

During and immediately after the clinical procedure, instruments should be treated to reduce and remove gross soil materials. Instruments should be wiped with low-lint or lint-free absorbent materials, such as low-lint or lint-free wipes to remove solid tissues and the majority of bodily fluids. Immediately after use, all disposable single-use accessories should be removed and discarded into appropriate biohazardous waste containers. Instruments should be placed into a transport container that will protect them from damage during transport and should be kept moist.

TRANSPORT TO PROCESSING AREA

Transport the instruments to the processing area as soon as practical. Minimize holding time before removing organic debris. Ensure that lighter instruments are not mixed with heavier instruments to avoid damage to the instruments due to mechanical abrasion. Seal or enclose the transport container to prevent fluid loss and the potential for cross-contamination of other areas.

PREPARATION

- Remove Koh Cup from RUMI or RUMI II Handle.
- Remove all single-use accessories and discard in an appropriate waste container designed for biohazardous materials if this has not been performed at point of use.
- The only devices used with the Koh Cup that are reusable are the RUMI or RUMI II Handles and the Advincula Arch™ Handles. Specific reprocessing instructions for the RUMI, RUMI II Handles and the Advincula Arch Handles are included with packaging for these devices.


Figure 1

Figure 2

Figure 3


Figure 4

CLEANING

Cleaning should take place as soon after use as practical. Soils should not be allowed to dry onto the instrument surfaces. The manual process and the automated process were validated independently and are considered equivalently effective when performed according to the instructions in this document. Do not use metal-bristled brushes on the instruments.

MANUAL CLEANING PROCESSES

1. Prepare a fresh solution of neutral or alkaline detergent specifically designed for manual cleaning of medical instruments. Prepare the solution at the detergent manufacturer's recommended concentration and temperature using purified or sterile water.
2. Submerge the Koh Cup in the freshly prepared detergent solution, ensuring that it is completely immersed.
3. Scrub all surfaces of the Koh Cup with a soft-bristled brush.
4. Rinse with warm purified or sterile water (100-120°F/38-49°C) for one (1) minute.
5. Prepare a neutral or alkaline enzymatic detergent specifically designed for manual cleaning of medical instruments. Prepare a fresh solution of detergent at the detergent manufacturer's recommended concentration using water at the detergent manufacturer's recommended temperature.
6. Soak the instrument for a minimum of ten (10) minutes at room temperature (68-77°F/20-25°C) or longer if indicated on the detergent manufacturer's instructions for use. Use the detergent manufacturer's recommended temperature if it is different from room temperature.
7. Rinse with warm purified or sterile water (100-120°F/ 38-49°C) for one (1) minute.
8. Examine the instrument for signs of residual soil. If there are any signs of residual soil on the instrument, repeat the cleaning process.
9. Dry the instrument with a low-lint or lint-free towel or wipe.

Automated Cleaning Process

1. Ensure that all gross soil has been wiped or rinsed from the surfaces of the instruments. All instruments within the scope of this document may be processed together in an automated washer or washer-disinfector.
2. Use only a legally marketed (FDA-cleared and/or CE-marked) medical device washer or medical device washer-disinfector. Industrial use and household use dishwashers are not acceptable for the cleaning of medical devices.
3. Place the instruments loosely within a rack designed for the medical device washer or medical device washer-disinfector, using a hold-down screen if indicated by the washer/washer-disinfector manufacturer.
4. Load the rack with the instruments into the washer/washer-disinfector and close the door.
5. Select a cycle intended for instruments that has the following parameters at a minimum:

PARAMETER	TOLERANCE
Detergent	Neutral or Alkaline medical device detergent
Pre-wash rinse	Minimum of one
	Minimum of 00:15 (mm:ss)
	Cold tap water or better quality
Wash	Minimum of one repetition
	Minimum of 02:00 (mm:ss)
	Warm tap or purified water (minimum of 140°F/60°C)
	Minimum concentration per detergent manufacturer's IFU
Post-wash rinse	Minimum of one
	Minimum of 00:15 (mm:ss)
	Warm tap or purified water (minimum of 110°F/43.3°C)
Thermal rinse*	Minimum of one
	Minimum of 01:00 (mm:ss)
	Purified water (minimum of 180°F/82.2°C)
Drying**	Minimum of default cycle setting

** Only heat-based medical device washer-disinfectors will include a thermal rinse. Medical device washers that do not include a thermal disinfection phase should be set to have a minimum of a 01:00 (mm:ss) rinse following the post-wash rinse using purified water.*
*** Drying time will be variable dependent on the size and composition of the load, altitude, environmental conditions, and air temperature and source characteristics.*

NOTE: automated wash processes that include additional phases or longer phase durations are anticipated to deliver equivalent or better cleaning efficacy as the minimal validated parameters. Additional phases and longer phase durations may be added if they are within the healthcare facility's normal procedures.

6. Start the cycle and allow it to run through conclusion.
7. Medical washer-disinfector default cycles often use high heat during the final stages of the processing cycle. Use caution and heat-protective gloves when opening and removing the wash rack at the end of the cycle.
8. Examine the instruments for signs of residual soil. If there are any signs of residual soil on the instrument, repeat the cleaning process.
9. If the instruments are not thoroughly dry following the automated process, dry the instruments with low-lint or lint-free towels or wipes.

2. INSPECTION

All instruments must be inspected for signs of damage, wear or residual soil each time they are reprocessed. If residual soil is observed during inspection, the cleaning process must be repeated for the instrument. Remaining soils can impede effective sterilization by shielding remaining microorganisms and can lead to tissue reactions in subsequent patients.

Surface Finish

Inspect the surface of the instrument. The surface should be smooth and evenly colored. Some staining and discoloration may occur depending on the exact formulation of the detergents used and the water quality. Surface staining and darkening may occur with the Ultem plastic models in cases of excessive heat, particularly if the device is in close contact with metal instrument surfaces during sterilization by moist heat. Staining will not impact performance of the Koh Cups as long as there are no texture changes in the surface. Corrosion on the surface of the stainless steel models of the Koh Cups may be a sign of a weakened material that could fracture under normal loads.

Cracks/Crazing

The Ultem plastic material is smooth and resistant to scratches. If cracks, scratches, gouges, crazing, or cuts are observed in the surface, the material may be damaged in a manner that could weaken its ability to support normal loads. The device should be replaced if these visual cues are present.

Debris

Look carefully at any areas that appear to be a different color than the overall instrument as these may be signs of residual soil adhered to the surface.

Attachment Ring

Pay attention to the shape and appearance of the inner ring area where the Koh Cup attaches to the snap drum of the handles. Gouges, distortions, cracks, corrosion, and other signs of surface damage in this location may impede proper attachment of the Koh Cup to the handle. Damage in this region indicates that the instrument should be replaced.

3. WRAPPING

Prior to terminal sterilization, the instruments must be packaged to preserve sterility after processing during storage. Legally marketed sterilization packaging (e.g. CSR wrap) that is compatible with and labeled for use with the chosen sterilization method should be used. CooperSurgical, Inc. has not evaluated the sterilization efficacy of the processes recommended in this document using containment systems (e.g. procedure/case trays, cassettes, rigid reusable sterilization containers), therefore, users who choose to use a containment system in the sterilization of these instruments should verify that the process is effective under these conditions. Containment systems can impede the ingress and egress of sterilant (moist heat). CooperSurgical, Inc. recommends individually wrapping the instruments prior to sterilization. The instruments should be wrapped using one of the common closure techniques, using properly sized material and secured with tape or adhesive designed and labeled for use with the sterilization process chosen. Instruments should be labeled with identifying information and sterilization date to facilitate selection for use after processing.

4. STERILIZATION - MOIST HEAT (STEAM)

Individually packaged instruments using legally marketed sterile barriers with labeling consistent with the following parameters may be sterilized prior to use and between subsequent uses. Legally marketed (FDA-cleared and/or CE-marked) moist heat (steam) sterilizers should be used. The following parameters have been validated for use with the instr ments within the scope of this document:

STERILIZATION PROCESS	EXPOSURE TEMPERATURE	EXPOSURE TIME	DRY TIME
Gravity Displacement	250°F/121°C	30 minutes	30 minutes
Pre-vacuum	270°F/132°C	4 minutes	30 minutes
Pre-vacuum	273°F/134°C	3 minutes	30 minutes

5. STORAGE

Individually packaged instruments should be cooled to a safe handling temperature following moist heat sterilization prior to final storage. Individually packaged instruments should be transported to a clean, temperature and humidity controlled, sterile supply storage area. If instrument packaging is compromised during storage (becomes wet, torn, or seals are open), the instrument must be reprocessed through cleaning, inspection, and sterilization prior to use. Do not stack instruments or allow other product or instruments to be stacked on top of them. Protection from dust and moisture is recommended. Ensure that traceability and identity is maintained throughout storage and up to subsequent use.

6. VALIDATION OF PROCESSES

CooperSurgical, Inc. performed validations for the cleaning of the instruments. The following detergents were used in the validations demonstrating efficacy:

Klenzyme™, Manu-Klenz®, ENZOL®, Neutrad®, Prolystica® Ultra Concentrate Neutral Detergent and Prolystica® Ultra-Concentrate Alkaline Detergent

CooperSurgical, Inc. does not endorse the use of the identified products in lieu of other similar products designed for use with medical devices, however, this information may provide a basis for comparing the formulations to select an appropriate alternate in facilities or locations where these detergents are not readily available. The parameters listed in the automated cleaning section as the minimum parameters for automated cleaning are identical to the parameters used in the validation, with the exception of the dry time. Post-cleaning dry time was eliminated for the cleaning process validation to facilitate detection of soil residuals, if present, in order to reduce the potential for a false negative result.

EXPLANATION OF SYMBOLS

Reorder number

Batch code

Use-by date
(CPO-6 Only)

Consult
instructions for use

Non-sterile

Do not re-use
(CPO-6 Only)

Do not use if
package is damaged

Caution

Not made with
natural rubber latex

Manufacturer

Sterilized using ethylene oxide (CPO-6 Only)

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician

Authorized Representative in the European Community

CooperSurgical is a registered trademark of CooperSurgical, Inc.

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RUMI® is a registered trademark of CooperSurgical, Inc.
Koh Colpotomizer™, Koh Cup™ and Colpo-Pneumo Occluder™ are trademarks of CooperSurgical, Inc.
Utem® is a registered trademark of SABIC Innovative Plastics.
Klenzyme™, Manu-Klenz® and Polystica® are trademarks of STERIS Corporation.
ENZOL® is a registered trademark of Johnson and Johnson.
Neutrad® is a registered trademark of Decon Labs, Inc.

95 Corporate Drive
Trumbull, CT 06611 USA
Phone: (800) 243-2974
Fax: (800) 262-0105

34145-IFU-C • Rev. A • 06/20

International
Phone: +1 (203) 601-9818
Fax: +1 (203) 601-4747
www.coopersurgical.com

© 2020 CooperSurgical, Inc

[EC]REP EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Made in the USA

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CooperSurgical is een gedeponeed handelsmerk van CooperSurgical, Inc.