Disposable Os Finder

920019

Instructions for Use (English)

CE2797

Read all safety information and instructions before using this product.

WARNINGS:

- Do not force passage
- For single use only. Do not re-use, reprocess or resterilize.
Re-use, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Re-use, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

CAUTION:

- Insert gently
Caution: U.S. Federal law restricts this device to sale by or on the order of a physician

ADVERSE REACTIONS:

- None recorded

CONTRAINDICATIONS:

- Pregnancy or possible pregnancy
- Unexplained Uterine Bleeding

DESCRIPTION/INDICATION FOR USE by physicians in clinics and hospitals

Total Length of each dilator is 21 cm (8") with a maximum usable length of 2.3 cm (1") for Os Locator, 4.5 cm (1 3/4") for Canal Finder, 8.5 cm (3.34") Fundus Finder. Useable width increases gradually to 3.8 mm (0.15") for each. The Os Finder set (Os Locator, Canal Finder, and Fundus Finder) is intended to gradually dilate the Os and is used for cervical stenosis and the post-menopausal patient.

INSTRUCTIONS FOR USE:

Three steps are to be followed:

- Os Locator is identified by a single ring around the handle. It is useful for locating the external Os by gently pushing the tip against the cervix.
- Canal Finder has a double ring around the handle. The slightly longer tip will insinuate the canal while effecting dilation.
- Fundus Finder has a triple ring around the handle. The longest of the three tips, this instrument can be passed through the entire length of the canal and will open the internal Os.

EXPLANATION OF SYMBOLS

 Reorder Number

 Batch Code

 Use-by date

 Consult instructions for use

 Do not re-use

 Do not re-sterilize

 Do not use if package is damaged

 Caution

 Not made with natural rubber latex

 U.S. Federal law restricts this device to sale by or on the order of a physician

 Sterilized Using Ethylene Oxide

 Authorized Representative in the European Community

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