GELLHORN

Multiple Drain Pessaries

Support for third-degree prolapse/procedentia
(Stage III or Stage IV) **€2797**

AVAILABLE MODELS

GELLHORN, FLEXIBLE	GELLHORN, 95% RIGID
REF MXPGE (Pessary Only)	REF MXPGRS (Pessary Only)
MXXPGPE (Pessary, TRIMO-SAN™, Jel Jector)	MXXPGRS (Pessary, TRIMO-SAN™, Jel Jector)

GELLHORN, SHORT STEM

REF MXPGS (Pessary Only)	REF MXPGRSS (Pessary, TRIMO-SAN™, Jel Jector)
---------------------------------------	--

Single-Patient Use RxOnly

Wash pessary with mild soap and thoroughly rinse prior to initial use.

DESCRIPTION

The Milex® GELLHORN Pessary is a medical device made from silicone that is inserted into the vagina to function as a supportive structure of the uterus, bladder or rectum.

WARNINGS

• Do not leave pessaries in place for long periods of time, as serious complications can occur which may require surgical intervention. Do not use these pessaries on a patient with a known silicone allergy.

CoperSurgical 38651-B • Rev. A • 11/2020

- Chemicals in various vaginal preparations can interact with the pessary material, resulting in discoloration or deterioration of the pessary. TRIMO-SAN™ does not interact with the pessary material. Use only TRIMO-SAN. Other materials have not been tested for compatibility.

INDICATIONS FOR USE

For effective support of third-degree prolapse or procedentia. The cervix rests behind the flat base of the pessary and only the stem shows in the vaginal entrance when the patient does a Valsalva maneuver. The GELLHORN pessary requires a relatively capacious vagina and an intact perineum.

CONTRAINDICATIONS

- The presence of pelvic infections and/or lacerations
- A noncompliant patient
- Endometriosis
- Pregnant patient
- Sexually active patient

RECOMMENDED CLEANING INSTRUCTIONS*

- Prepare a cleaning solution by mixing a mild soap (such as Dawn® or equivalent) with tap water using the soap manufacturer's recommended concentration. Prepare this solution in a container large enough to fully submerge the device.
- Soak and Scrub
 - Soak the device in the container of prepared soap solution for a minimum of 5 minutes.
 - Following the 5-minute soak period, scrub the device for a minimum of 15 seconds with a soft-bristled brush, such as a tooth brush and/or pipe brush. Scrub device below water line to prevent aerosolization of contaminants.
 - Following scrub, inspect device for visible soil residue.
- Rinse
 - Remove the device from the soap solution and thoroughly rinse under flowing tap water for a minimum of 30 seconds.
 - Allow the device to dry.

*New pessaries are powdered with food-grade powder.

PATIENT INSTRUCTIONS FOR PESSARY USE

CAUTION

COOPERSURGICAL RECOMMENDS PESSARIES BE INSERTED AND REMOVED BY THE PHYSICIAN OR OTHER HEALTHCARE PROFESSIONAL UNLESS OTHERWISE DIRECTED.

- To ensure the desired correction of your condition, your healthcare professional needs your full cooperation.
- It is essential that your healthcare professional inspect your vagina at frequent intervals for evidence of pressure and/or allergic reaction.
- A gradual increase in the interval of inspection may be considered at the discretion of the treating practitioner.
- It is important to evaluate the effects of the vaginal gel TRIMO-SAN.

TRIMO-SAN is available at the prescription counter of your pharmacy. A prescription is required.

REPORT ANY OF THE FOLLOWING SYMPTOMS TO YOUR PHYSICIAN

- Any difficulty in urinating
- Any discomfort
- Any changes in the color or consistency of vaginal discharge
- Any increase in the amount of vaginal discharge or vaginal bleeding
- Any foul odor associated with vaginal discharge
- Vaginal itching
- If the pessary falls out

For medical emergencies and for all medically related advice, consult your healthcare professional.

FOR THE PHYSICIAN / HEALTHCARE PROFESSIONAL

Review these instructions with the patient to establish use regimen.

INSTRUCTIONS

- Wear dry gloves. When necessary, lubricate only the entering end of the pessary with TRIMO-SAN. Hold the pessary as shown in Figure 1.
- Use one finger to depress the perineum.
- Guide the pessary, inserting it edgewise almost parallel to the


Figure 1 Figure 2

introitus (see Figure 1), avoiding the urethral opening while the perineum is strongly pushed downward. Use a corkscrew motion while introducing the GELLHORN into the vagina.

- Once the large flat disc is past the introitus, push the pessary upward until only the end of the stem shows in the vaginal entrance. The cervix rests behind the flat disc (see Figure 2).
- Have the patient sit, stand and bear down. Examine the patient while she is in the standing position to ensure the pessary has not shifted position. The patient should not feel the pessary once it is in position. The pessary should not be too loose as it may turn or be expelled and it should not be too tight as it may cause discomfort.
- The healthcare professional should be able to sweep one finger between the pessary and vaginal walls. If there is not enough space to do this, the next smaller size should be tried. If excessive space exists, the pessary will not be effective and may rotate or even be expelled.
- It may be necessary to refit the patient with a different size or type of pessary after a period of time. Do not assume that a replacement will always be the same size as the previous one. Check the fitting to ensure continued patient comfort and relief of symptoms. The useful life of a pessary is limited. Examine frequently for signs of deterioration (such as cracks or breaks in silicone outer surface). A pessary should be replaced if damaged.

Note: If the patient is unable to urinate with the pessary in position, remove it and fit her with the next smaller size. Repeat as necessary.

- If the patient can void without difficulty, pessary remains in position upon re-examination, and the patient is comfortable with the pessary in place, this is a good indication that the correct size may have been selected. Patient experience may vary.
- Ulcerations and erosions frequently occur in cases of complete prolapse due to irritation of the exteriorized cervix. Whenever possible, reducing the mass and treating the irritation are primary

steps before using a pessary. Prolapse reduction may resolve cervical vaginal irritation. Verification of cervical cytology (pap) and/or biopsy as clinically indicated.

- During each visit, the vagina should be carefully inspected for evidence of pressure or allergic reaction. The patient should be questioned concerning douching, discharge, disturbance of bowel function or urination. It may be necessary to fit another size or an entirely different type of pessary.

- At the physician's discretion, the patient can be instructed in the proper removal, cleaning and reinsertion techniques for her own pessary. This process can be performed nightly or even weekly by the patient under ideal circumstances.

TO REMOVE

- Use one finger to depress the perineum.
- Use other hand to grasp the "knob", pulling the pessary away from the cervix, turning the pessary so that the disc is almost parallel to the introitus. Using corkscrew motion ease the pessary out. See Figure 3.
- Removal may be facilitated by passing a finger along the stem and behind the disk and folding it against the stem for removal.

Figure 3

RECOMMENDED FOLLOW-UP

- Have patient return within 24 hours for first examination.
- Have patient return for second examination within 3 days.
- Have patient return for examination every few months.

Note: The above schedule of follow-up examinations may be altered to fit the needs of the individual patient at the discretion of the healthcare provider.

EXPLANATION OF SYMBOLS

REF Reorder number

LOT Batch code

Use-by date

Consult instructions for use

Not made with natural rubber latex

RxOnly CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician

EC/REP Authorized Representative in the European Community

Manufacturer

Dawn® is a registered trademark of Procter & Gamble, Cincinnati, Ohio 45202

TRIMO-SAN™ is a trademark of CooperSurgical, Inc.

Milex® and CooperSurgical are registered trademarks of CooperSurgical, Inc.

© 2020 CooperSurgical, Inc.

Made in the USA

CoperSurgical
95 Corporate Drive
Trumbull, CT 06611 USA
Phone: (800) 243-2974
Fax: (800) 262-0105
International
Phone: +1 (203) 601-9818
Fax: +1 (203) 601-4747
www.coopersurgical.com

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TRIMO-SAN™ est une marque de CooperSurgical, Inc.
Milex® et CooperSurgical sont des marques déposées de CooperSurgical, Inc.
TRIMO-SAN™ är ett varumärke som tillhör CooperSurgical, Inc.
Milex® och CooperSurgical är registrerade varumärken som tillhör CooperSurgical, Inc.