H/S Catheter Set with Integrated Stylet for Hysterosonography and Hysterosalpingography

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H/S Catheter Set with Integrated Stylet for Hysterosonography and Hysterosalpingography • 61-3005 and 61-3007

Instructions for Use (English)

DEVICE DESCRIPTION

The H/S Catheter Set consists of a latex-free balloon catheter, insertion sheath and a 1.5 cc (ml) syringe (5F catheter) or 3.0 cc (ml) syringe (7F catheter). (See Figure 1.) The catheter can be used with aqueous based contrast media. Not made with natural rubber latex.


Figure 1

WARNINGS

- Contents supplied sterile. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

INTENDED USE/INDICATIONS

For administering contrast media during Hysterosalpingography or Hysterosonography procedures to detect uterine pathology such as polyps, fibroids, adhesions or endometrial thickening, and/or patency of fallopian tubes.

CONTRAINDICATIONS

Suspected infection, suspected pregnancy, profuse bleeding or sexually transmitted disease.

PRECAUTIONS

Do not exceed the recommended balloon inflation volume of 1.5 cc (ml) (5F catheter) or 3.0 cc (ml) (7F catheter) or the balloon may burst. The use of OIL-BASED contrast media such as ethyl esters may interact with the balloon of the catheter, causing rupture. The use of oil-based contrast media is not recommended.

ADVERSE EVENTS

Some patients may have a hypersensitivity to contrast media.

INSTRUCTIONS FOR USE

Catheter Preparation

1. Grasp the insertion sheath connected to the trapezoidal hub.
2. Remove and discard the crimped protective cover, exposing the balloon catheter tip.
3. Test the balloon integrity by inflating with air, saline or water using the syringe in the set. If performing Hysterosonography remove as much air as possible from the balloon. Deflate completely by pulling back on the syringe plunger and closing the stopcock.
4. Attach a contrast media-filled syringe (not supplied) to the blue luer hub and fill the catheter with contrast media to expel air. Advance the insertion sheath so that the distal end of the catheter protrudes slightly from the distal end of the sheath. (See Figure 2.)

Figure 2

Catheter Placement

5. Visualize the external cervical os with the aid of a speculum and advance the sheath and catheter so that the tip of the catheter enters the cervical canal.
6. The catheter may be shaped/curved to facilitate insertion.
7. Hold the sheath stationary and advance the catheter into the cervical canal and into the uterine cavity.
8. Open the stopcock and slowly inflate the balloon with up to 1.5 cc (ml) (5F catheter) or 3.0 cc (ml) (7F catheter) of air, saline or water. Turn the stopcock off allowing the balloon to remain inflated.
9. Withdraw the catheter until it occludes the internal os of the cervix.

Injection of Contrast Media

10. Inject the contrast media, close the clamp and complete the study in a routine manner.
11. Occasionally, when access to the uterine cavity is difficult, the balloon must be inflated within the endo-cervical canal.

Catheter Removal

12. Open the stopcock and deflate the balloon by pulling back on the syringe plunger. Withdraw the catheter.
13. Dispose of in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.

EXPLANATION OF SYMBOLS

REF Reorder Number

LOT Batch Code

 Use-by date

 Do not re-use

 Do not resterilize

 Consult instructions for use

 Do not use if package is damaged

 Caution

 Not made with natural rubber latex

 Manufacturer

STERILE EO Sterilized using ethylene oxide

EC REP Authorized Representative in the European Community

Rx Only U.S. Federal law restricts this device to sale by or on the order of a physician

CooperSurgical is a registered trademark of CooperSurgical, Inc.

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© 2020 CooperSurgical, Inc.
CooperSurgical è un marchio depositato di CooperSurgical, Inc.