H/S Elliptosphere™ Catheter Set
for Hysterosonography and Hysterosalpingography
Instructions for Use (English)

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H/S Elliptosphere™ Catheter Set for Hysterosonography and Hysterosalpingography • 61-4005

Instructions for Use (English)

DEVICE DESCRIPTION

The H/S Elliptosphere™ Catheter Set consists of a latex free balloon catheter, insertion sheath and a 1.5 cc (ml) syringe. (See Figure 1.) The catheter can be used with aqueous or oil-based contrast media. Not made with natural rubber latex.


Figure 1

WARNINGS

- Contents supplied sterile. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

INTENDED USE/INDICATIONS

For administering contrast media during Hysterosalpingography or Hysterosonography procedures to detect uterine pathology such as polyps, fibroids, adhesions or endometrial thickening, and/or patency of fallopian tubes.

CONTRAINDICATIONS

Suspected infection, suspected pregnancy, profuse bleeding or sexually transmitted disease.

PRECAUTIONS

Do not exceed the recommended balloon inflation volume of 1.5 cc (ml) or the balloon may burst.

ADVERSE EVENTS

Some patients may have a hypersensitivity to contrast media.

INSTRUCTIONS FOR USE

Catheter Preparation

1. Grasp the insertion sheath connected to the trapezoidal hub.
2. Remove and discard the crimped protective cover, exposing the balloon catheter tip.
3. Test the balloon integrity by inflating with air or contrast media using the syringe in the set. Deflate completely by pulling back on the syringe plunger and closing the stopcock.
4. Attach a contrast media-filled syringe (not supplied) to the blue luer hub and fill the catheter with contrast media to expel air. Advance the insertion sheath so that the distal end of the catheter protrudes slightly from the distal end of the sheath. (See Figure 2.)

Figure 2

Catheter Placement

5. Visualize the external cervical os with the aid of a speculum and advance the sheath and catheter so that the tip of the catheter enters the cervical canal.
6. Hold the sheath stationary and advance the catheter into the cervical canal 0.4-0.8 in. (1-2 cm).
7. Open the stopcock and slowly inflate the balloon with up to 1.5 cc (ml) of either air or saline to provide occlusion without causing discomfort to the patient. Turn the stopcock off allowing the balloon to remain inflated.

Injection of Contrast Media

8. Inject the contrast media, close the clamp and complete the study in a routine manner.
9. If uterine placement is desired, deflate the balloon by opening the stopcock and pulling back on the syringe plunger. Advance the catheter into the uterus, fully inflate the balloon using the enclosed syringe, close the stopcock and withdraw the catheter until the balloon occludes the internal os.

Catheter Removal

10. Open the stopcock and deflate the balloon by pulling back on the syringe plunger. Withdraw the catheter.
11. Dispose of in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.

EXPLANATION OF SYMBOLS

- REF** Reorder Number
- LOT** Batch Code
- Use-by date
- Do not re-use
- Do not resterilize
- Consult instructions for use
- Do not use if package is damaged
- Caution
- Not made with natural rubber latex
- Manufacturer

STERILE EO Sterilized using ethylene oxide

EC REP Authorized Representative in the European Community

R_x Only U.S. Federal law restricts this device to sale by or on the order of a physician

Marks bearing™ are trademarks of CooperSurgical, Inc
CooperSurgical is a registered trademark of CooperSurgical, Inc

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© 2020 CooperSurgical, Inc.