H/STYLET™ FOR USE WITH H/S CATHETERS

61-0020

CE 2797

Instructions for Use (English)

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

CAUTION: This device is intended for use in the infusion lumen. Use in the balloon inflation lumen may result in difficulty to use or damage to the catheter assembly. For use with H/S Catheters 61-5005, 61-5007 and 61-4005 only.

WARNING: DO NOT USE if the stylet is protruding from the distal tip of the catheter assembly, this may cause patient discomfort or injury.

- Contents supplied sterile. Do not use if sterile barrier is damaged.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

Not made with natural rubber latex.

INSTRUCTIONS FOR USE

Use when additional stiffness is considered necessary.

1. Remove the protective cap from distal end of stylet and discard.
2. Insert the stylet into the blue luer infusion hub of the H/S Catheter and advance it.
3. Secure the stylet to the blue luer by twisting the plastic rotating luer lock of the stylet onto blue luer.
4. The catheter may be shaped/curved to facilitate insertion.
5. Insert the H/S Catheter with stylet into the cervical canal; position the balloon as per the specific H/S Catheter instructions.
6. Slowly remove the stylet from the H/S Catheter.

EXPLANATION OF SYMBOLS

REF	Reorder Number		Do not re-use		Consult instructions for use	Rx Only	U.S. Federal law restricts this device to sale by or on the order of a physician
LOT	Batch Code		Do not re-sterilize		Do not use if package is damaged	EC REP	Authorized Representative in the European Community
	Use-by date		Sterile		Manufacturer		
STERILE R	Sterilized using Irradiation		Caution		Not made with natural rubber latex		

CooperSurgical is a registered trademark of CooperSurgical, Inc.

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© 2020 CooperSurgical, Inc.
CooperSurgical è un marchio depositato di CooperSurgical, Inc.