
LEEP Disposable Universal Electrodes

STERILE	Unless package has been opened or damaged. Ethylene Oxide Gas Sterilized.
CAUTION	U.S. Federal law restricts this device to sale by or on the order of a physician.
CAUTION	Electrodes may be used with Electrosurgical Generators with a maximum HF Voltage rating of 5000 volts.

CE 2797

Instructions for Use (English)

DESCRIPTION

CooperSurgical's LEEP Disposable Universal Electrodes are offered in two basic types: one is designed for electro-excisional procedures. The second is used for electro-fulguration procedures for diagnosing non-invasive HPV-related or unrelated lesions of the lower genital tract and for treating non-invasive HPV-related or unrelated lesions of the lower genital tract (see Indications for Use).

DIAGNOSTIC USAGES (Electroexcision)

- Intraepithelial neoplasia of lower genital tract
- Condylomata acuminata
- Lesions such as polyps, molluscum, contagiosum, nevi, seborrheic keratosis, achricordon, etc.

THERAPEUTIC USAGES (Electrofulguration)

- Intraepithelial neoplasia of lower genital tract
- Condylomata acuminata
- Symptomatic conservative treatment-resistant ectropion
- Lesions such as polyps, molluscum, contagiosum, nevi, seborrheic keratosis, achricordon, etc.

INDICATIONS FOR USE

For electroexcisional procedures, loop and needle electrodes are used with BLEND-1 cut current whereas electrofulguration is performed with ball-shaped or needle-type electrodes using coagulation power outputs. All electrodes are disposable and have an insulated shaft measuring 5.5 cm to 12 cm in length, and loop electrodes, a T-shaped crossbar. The latter is essential to prevent heat transfer to tissue and thermo-coagulating the lesional epithelium during excision.

Loop and needle electrodes are made of 0.2 mm hard tungsten wire with high melting temperature and excellent conductiveness of alternating current. To prevent excessive tissue removal, loop electrodes are shallow, the largest measures 1.0 cm deep x 2.0 cm large, and smallest 0.4 cm deep x 1.0 cm large. The former is designed for LEEP of the cervix, whereas the latter for LEEP of vaginal and external anogenital lesions. It is not recommended to use larger than a 1.0 cm deep loop electrode because they may inadvertently lead to excess tissue removal, and in the dense fibrous tissue, tend to bend. When wire bends, the crossbar precedes the wire and the electrode has to be removed prematurely once the crossbar has reached the vaginal wall. This results in an incomplete excision of the lesional tissue. The maximum depth of tissue that can be removed with a single excisional pass using a 1.0 cm deep electrode is 0.9 cm because the crossbar of the electrode prevents deeper penetration, with a 0.8 cm deep electrode, the depth is 0.7 cm, etc.

Loop-shaped electrodes are used for non-invasive cervical lesions,

whereas those with square or rectangular configurations are generally used for HPV-related condylomata acuminata and intraepithelial neoplasia of the vagina, vulva and perineal area. The flat shape of the square/rectangular electrodes allows for better and even depth control than the loop or radius type electrode.

The 8 mm long thin needle electrode with a 12 cm long shaft is designed for electro-conization of deeply located endocervical lesions step-by-step rather than in one stroke when using loop electrodes. The thin needle electrode with a 5.5 cm long shaft may be used for excising larger condylomata and lesions suspected of invasive disease in the external anogenital tissue.

Hemostasis of bleeding sites and/or ablation of subclinical lesional tissues are achieved by the 5 mm large ball or the 0.8 cm thick macro-needle electrode using coagulation power outputs. Those with 12 cm long shafts are designed for the vagina and cervix whereas those with a 5.5 cm shaft are used for external anogenital lesions. The macro-needles are particularly useful to coagulate bleeding vessels with minimal destruction of adjacent, normal tissues and virtually bloodless excision (using coagulation current) of very large condylomata and thrombosed hemorrhoidal tissues.

CONTRAINDICATIONS

Absolute contraindications include hypovolemic shock, cardiac decompensation (Class IV) and coagulopathy/blood dyscrasia; acute/chronic severe, active infection caused by Herpes, Chlamydia or Gonococcus and other aerobic and anaerobic organisms; intraepitheal disease extending to or close to lower uterine segment; severe allergic reaction to local anesthetic and iodine solutions; and a refusal to sign an informed consent form.

Relative contraindications include clinically apparent invasive cancer (punch biopsy is sufficient for histologic diagnosis); pregnancy (except in expert hands when invasion is suspected); less than 3 months post-partum and hypoplastic DES-cervix.

PRECAUTIONS

- LEEP is a powerful surgical procedure with occasional intra-operative and post-operative complications. You MUST obtain informed signed consent prior to the procedure.
- Severe bleeding may rarely develop after deep excisional procedures and may be difficult to impossible to control by direct pressure, electro-fulguration and application of AstrinGyn® gel or Monsel's paste. You must have the appropriate instruments for suturing blood vessels in the office.
- Adequate anesthesia and complete visualization of the lesional area should be obtained before initiating LEEP.

CooperSurgical

- Use only specially designed insulated instruments for performing LEEP procedures, to avoid burns to patient or physician.
- For maximum safety, the vaginal walls may be protected from activated electrodes and inadvertent cuts and burns by a latex condom placed on any of the LEEP specula. The use of condoms or LEEP vaginal retractor (*F400H*) is of great help to retract vaginal walls allowing for good visualization of the cervix and vaginal fornix. Check first for patient allergy to latex before using condoms.
- Prior to any loop electrosurgical procedures, make a “dry pass” with the electrode under colposcopic guidance. This will determine that the path of the activated electrode will not be obstructed or come in contact with non-insulated or anesthetized tissue.
- Electrosurgical excision with loop electrodes must be relatively slow but continuous. If the electrode must be stopped or stalled, remove electrode from the tissue and begin from the opposite side and attempt at meeting the original line of excision. In this way, excessive thermal damage to tissue will be avoided.
- When excising at the periphery of large exocervical lesions, do not extend electrode deeper than 3 mm to avoid cutting into the cervical branches of the uterine artery and producing important bleeding.
- To increase good visualization during LEEP, it is necessary to have an adequate fume evacuation vacuum-suction system. Make sure it functions properly before initiating LEEP.
- Make sure that the electrosurgical unit (ESU) is maintained in the “Stand-by” mode until electrosurgery begins.
- Physicians with limited experience with LEEP will obtain a stable hand by placing the shaft of the loop electrode on the posterior blade of the LEEP speculum. Similarly, when excising external ano-genital lesions, good depth control is achieved by resting the 5th finger of the hand that holds the square-shaped electrode on the index finger of the other hand which is used to spread the skin at the excision site.
- These electrodes are for single use **ONLY**. After use, discard per local hazardous waste procedures.

DIRECTIONS FOR USE • GENERAL OPERATION

Loop Electrosurgical Excision Procedure (LEEP) of the cervix should not be used unless presence of an unequivocal lesion is ascertained by cytology or colposcopy. If in doubt, proceed with traditional diagnostic triage using endocervical curettage and cervical biopsies.

LEEP can be performed in the office for localized lesions; however, general anesthesia or epidural/spinal anesthesia is needed for extensive external anogenital lesions and deeply located endocervical intraepithelial neoplasia.

USE IN THE CERVIX

- Determine size and distribution of lesion.
- Provide local anesthetic.
- Select electrode that will allow a one pass excision.
- If the lesion is large, excise central area first (that is most likely to be the seat of occult invasion), then proceed with additional passes as required.
- Select the appropriate power settings in BLEND-1 cutting mode. These depend on the size, composition of the electrode and the electrosurgical generator being used.
- Turn on the vacuum suction, activate the ESU and proceed with continuous, slow motion excisions.
- After excision, inspect crater base and endocervical canal using an endocervical specula to insure all lesional tissue has been removed.
- Bleeding vessels are electrofulgurated using the 5 mm ball or macro-needle electrode in coagulation power output.

AstrinGyn® gel is applied to crater beds.

- Provide written post-operative instructions to the patient that should include a return to the office if persistent bleeding and/or pelvic infection develops; if not, in 3 to 6 months for a follow-up visit.

USE IN THE ENDOCERVICAL SQUAMOUS INTRAEPITHELIAL

There are 2 approaches to this type of lesion —

1. Excise endocervical portion first with either a large 2.0 cm x 1.0 cm loop or a small square-shaped 1.0 cm x 1.0 cm loop electrode and repeat procedure once. The resulting depth of the crater generally reaches 1.8 cm. Additional depth may be obtained if required.
2. As an alternative to loop electrodes, use the 8 mm thin needle electrode, excise lesional tissue step-by-step (as with a scalpel blade) and detach base of cone with the needle electrode or, to avoid thermocoagulation of the endocervical margins, with a tonsillar snare with a #8 wire.

USE IN THE VAGINA

Identify lesions; elevate lesions to 1 cm above the adjacent, normal epithelium by infiltrating the sub-epithelial connective tissue stromat with anesthetic solution. Place the shaft of the vaginal electrode on the posterior blade of the speculum (this provides for a steady hand). Excise large lesions with a shallow 1.0 cm x 0.4 cm square-shaped wire bent to 15° (*S1004*). Electrofulgurate bleeding sites and flat, subclinical lesions with the 5 mm ball electrode in coagulation power.

USE IN THE EXTERNAL ANOGENITAL

(vulva, perineum, anus, penis)

- Elevate lesions using local anesthesia as described for vaginal lesions. Excise raised lesions with a shallow short-shafted (1.0 cm x 0.4 cm) square-shaped electrode (*S1004*).
- Depth control is achieved by properly placing hand close to lesional tissue.
- Step-by-step electroexcisional procedures can be achieved with the 0.8 mm thin needle electrode (*N0550*).

Flattened subclinical lesions are best ablated with ball or macro needle electrodes using coagulation power outputs. These electrodes are also used for providing hemostasis to bleeding sites; to ascertain whether all lesional tissue has been destroyed, and to obtain good healing, the hard eschar which is produced by fulguration is removed by a 3 mm or 5 mm large Fox curette or saline-moistened cotton-tipped applicators. The treatment fields are covered with 1% silver sulphadiazine cream. Patients are instructed to return at six weeks after LEEP for external anogenital disease or earlier if complications such as bleeding or infection develop.

EXPLANATION OF SYMBOLS

 REF	Catalog/Reorder Number		Shock Hazard Warning
 LOT	Lot number		Consult instructions for use.
	Use-by date		Do not use if package is damaged
 STERILE EO	Sterilized using ethylene oxide		Manufacturer
R_x Only	Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.	 EC REP	Authorized Representative in the European Community.
	Do not re-use		Not made with natural rubber latex
	Do not re-sterilize		BF applied part
	Caution		

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© 2020 CooperSurgical, Inc.
AstrinGyn® is a registered trademark of CooperSurgical, Inc.