Milex® Pessary Fitting Kit

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Milex[®] Pessary Fitting Kit

Fitting Instructions for the Healthcare Professional

GENERAL INFORMATION

PURPOSE

The Milex[®] Pessary Fitting Kit was developed to help you as a healthcare professional determine accurate pessary type and size prior to prescribing a pessary for your patient.

INDICATIONS FOR USE

See individual Pessary sections in this document for indications.

CONTENTS

The pessaries included in this kit are: INCONTINENCE DISH with Support, RING with Support, DONUT, GELLHORN and CUBE. Also included in this kit are a foam model, carrying case and pessary size chart.

TRIMO-SAN[™]* and a Universal Applicator with instructions are only included in kits sold in the U.S.

WARNING: Fitting pessaries are not to be worn more than 15 minutes.

FITTING INSTRUCTIONS

IMPORTANT

- Pessaries are fitted by trial and error. There are no mechanical devices available that can accurately determine the size or type of pessary your patient requires to obtain the desired results.
- Even before fitting a pessary, the patient should be informed that it is not uncommon to have to change the size or type of pessary more than once after being originally fitted. This is why it is so important that your patient be instructed to return within 24 hours of the initial fitting and again in 72 hours. Thereafter, re-examination every few months is recommended to ensure that a proper fit is maintained as long as the patient is wearing the pessary.
- At each visit the pessary should be removed and the vaginal vault inspected for signs of allergic reaction or undue pressure.
- At the physician's discretion, the patient can be instructed in the proper removal, cleaning and reinsertion techniques for her own pessary. This process can be performed nightly by the patient under ideal circumstances.

Note: When you purchase the pessaries for your patients, detailed fitting instructions are included in the box for your use along with pessary instructions for the patients.

RECOMMENDED CLEANING, HIGH LEVEL DISINFECTION AND STERILIZATION INSTRUCTIONS**

CLEANING

- Prepare a cleaning solution by mixing a mild soap (such as Dawn[®] or equivalent) with tap water using the soap manufacturer's recommended concentration. Prepare this solution in a container large enough to fully submerge the device.

2. Soak and Scrub

- Soak the device in the container of prepared soap solution for a minimum of 5 minutes.
- Following the 5-minute soak period, scrub the device for a minimum of 15 seconds with a soft-bristled brush, such as a tooth brush and/or pipe brush. Scrub device below water line to prevent aerosolization of contaminants.
- Following scrub, inspect device for visible soil residue.

3. Rinse

- Remove the device from the soap solution and thoroughly rinse under flowing tap water for a minimum of 30 seconds.
- Allow the device to dry.

HIGH LEVEL DISINFECTION

- Prepare a container of Revital-Ox[™] RESERT[®] STERIS, Cidex[®] OPA or equivalent disinfectant in a container large enough to fully submerge the device.
- Transfer the device to the container to soak. Ensure that the device is completely immersed in the disinfecting solution. Soak the device for at least 9 minutes in RESERT or 12 minutes in Cidex OPA.
- Rinse the device thoroughly in sterile water.
- Place the cleaned/disinfected device in an appropriate dry storage area.

STERILIZATION

Sterilization Process	Exposure Temperature	Exposure Time	Dry Time
Pre-Vacuum	273 °F / 134 °C	3 minutes	30 minutes
Pre-Vacuum	270 °F / 132 °C	4 minutes	30 minutes
Gravity	250 °F / 121 °C	30 minutes	30 minutes

INCONTINENCE DISH with Support Fitting Pessary

For stress urinary incontinence with mild prolapse

Fitting pessary sizes available in kit and corresponding pessary sizes

INCONTINENCE DISH with Support				
Fitting Kit Replacement Pessary Part No.	Size	Product Part No.	inches	mm
MXFIT0801	2	MXPCONDS02	2-1/2"	65 mm
MXFIT0802	3	MXPCONDS03	2-3/4"	70 mm
MXFIT0803	4	MXPCONDS04	3"	75 mm
MXFIT0804	5	MXPCONDS05	3-1/8"	80 mm

* TRIMO-SAN is not sold outside the USA.

**New pessaries are powdered with food-grade powder.

DESCRIPTION

The Millex® INCONTINENCE DISH Pessary is a medical device made from silicone that is inserted into the vagina to function as a supportive structure of the uterus, bladder and/or rectum.

WARNINGS

- Do not use these pessaries on a patient with a known silicone allergy.
- Chemicals in various vaginal preparations can interact with the pessary material, resulting in discoloration or deterioration of the pessary. TRIMO-SAN™* does not interact with the pessary material. Use only TRIMO-SAN*. Other materials have not been tested for compatibility.

INDICATIONS FOR USE

Stress urinary incontinence

- The INCONTINENCE DISH with Support is indicated if stress urinary incontinence with a mild uterine prolapse is complicated by a mild cystocele.

CONTRAINDICATIONS

- The presence of pelvic infections and/or lacerations
- A noncompliant patient
- Endometriosis
- Pregnant patient

INSTRUCTIONS

Review these instructions with the patient to establish use regimen.

1. Wear dry gloves. When necessary, lubricate only the entering end of the pessary with TRIMO-SAN* or equivalent vaginal jelly. Compress the pessary (bringing sides together) as shown in Figure 1.
2. Use one finger of the opposite hand to depress the perineum. Hold the pessary almost parallel with the introitus (see Figure 2). Direct the entering end of the pessary past the cervix into the posterior fornix. Allow the pessary to open into shape after passing the introitus.
3. Use the index finger to guide the pessary along the lower vaginal wall, behind the cervix, and into the posterior fornix.
4. Use the index finger to bring the knob up behind the symphysis pubis (see Figure 3).

Note: If the patient is unable to urinate with the pessary in position, remove it and fit her with the next smaller size. Repeat as necessary.

5. If the patient can void without difficulty, and the pessary remains in position upon re-examination, and the patient is comfortable with the pessary in place, this is a good indication that the correct size may have been selected. Patient experience may vary.
6. Examine the patient while she is in the standing position to ensure the pessary has not shifted position. The patient should not feel the pessary once it is in position. The pessary

Figure 1


Figure 2

Figure 3

should not be too loose as it may turn or be expelled and it should not be too tight as it may cause discomfort.

7. The healthcare professional should be able to sweep one finger between the pessary and the vaginal walls. If there is not enough space to do this, the next smaller size should be tried. If excessive space exists, the pessary will not be effective and may rotate or even be expelled.
8. It is sometimes necessary to refit the patient with a different size or type of pessary after a period of time. Do not assume that a replacement will always be the same size as the previous one. Check the fitting to ensure continued patient comfort and relief of symptoms. The useful life of a pessary is limited. Examine frequently for signs of deterioration (such as cracks or breaks in silicone outer surface). A pessary should be replaced if damaged.

TO REMOVE

Use the index finger to depress the perineum. Hook other index finger under the knob of the INCONTINENCE DISH pessary and pull down. Fold the pessary by bringing the sides together, angling it so that it is almost parallel to the introitus, and gently ease the pessary out.

RING with Support Fitting Pessary

For uterine prolapse (Stage I and Stage II)

Fitting pessary sizes available in kit and corresponding pessary sizes

RING with Support				
Fitting Kit Replacement Pessary Part No.	Size	Product Part No.	inches	mm
MXFIT0020	2	MXPRS02	2-1/4"	57 mm
MXFIT0021	3	MXPRS03	2-1/2"	64 mm
MXFIT0022	4	MXPRS04	2-3/4"	70 mm
MXFIT0023	5	MXPRS05	3"	76 mm
MXFIT0024	6	MXPRS06	3-1/4"	83 mm
MXFIT0025	7	MXPRS07	3-1/2"	89 mm

DESCRIPTION

The Millex® RING Pessary is a medical device made from silicone that is inserted into the vagina to function as a supportive structure of the uterus, bladder and/or rectum.

WARNINGS

- Do not leave pessaries in place for long periods of time, as serious complications can occur which may require surgical intervention.
- Do not use these pessaries on a patient with a known silicone allergy.
- Chemicals in various vaginal preparations can interact with the pessary material, resulting in discoloration or deterioration of the pessary. TRIMO-SAN™* does not interact with the pessary material. Use only TRIMO-SAN*. Other materials have not been tested for compatibility.

INDICATIONS FOR USE

- The RING with Support pessary is indicated for support in Stage I and Stage II prolapse complicated by a mild cystocele.

* TRIMO-SAN is not sold outside the USA.

CONTRAINDICATIONS

- The presence of pelvic infections or lacerations
- A noncompliant patient
- Endometriosis
- Pregnant patient

INSTRUCTIONS

Review these instructions with the patient to establish use regimen.

1. Wear dry gloves. When necessary, lubricate only the entering end of the pessary with TRIMO-SAN* or equivalent vaginal jelly. Hold as illustrated in Figure 1.
2. The pessary is folded along the axis of the bigger outer holes by bringing the small round holes together (see Figure 2). The arch formed points downward as shown with the RING pessary.
3. Direct the pessary past the cervix into the posterior fornix. Allow the pessary to open again into the ring shape after passing the introitus.
4. The index finger is inserted deep into the vagina to turn the pessary approximately 90° (see Figure 3). The RING pessary in this position cannot be folded and pushed out.
5. Ask the patient to sit, stand and bear down slightly. If there is no leakage, and the patient is comfortable with the pessary in position, have her empty her bladder. A properly fitted pessary takes up slack in redundant tissue, holding the uterus higher in the vagina.
6. If the patient can void without difficulty, and the pessary remains in position upon re-examination, and the patient is comfortable with the pessary in place, this is a good indication that the correct size may have been selected. Patient experience may vary.

Note: If the patient is unable to urinate with the pessary in position, remove it and fit her with the next smaller size. Repeat as necessary.

7. Examine the patient while she is in the standing position to ensure the pessary has not shifted position. The patient should not feel the pessary once it is in position. The pessary should not be too loose as it may turn or be expelled and it should not be too tight as it may cause discomfort.
8. The healthcare professional should be able to sweep one finger between the pessary and the vaginal walls. If there is not enough space to do this, the next smaller size should be tried. If excessive space exists, the pessary will not be effective and may rotate or even be expelled.
9. It is sometimes necessary to refit the patient with a different size or type of pessary after a period of time. Do not assume that a replacement will always be the same size as the previous one. Check the fitting to ensure continued patient comfort and relief of symptoms. The useful life of a pessary is limited. Examine frequently for signs of deterioration (such as cracks or breaks in silicone outer surface). A pessary should be replaced if damaged.

TO REMOVE

- Use one finger to depress the perineum. Turn the pessary until the notches face the introitus. Fold the pessary and gently ease it out.

Figure 1

Figure 2

Figure 3

DONUT Fitting Pessary

Support for third-degree prolapse/procidentia (Stage III + Stage IV)

Fitting pessary sizes available in kit and corresponding pessary sizes

DONUT				
Fitting Kit Replacement Pessary Part No.	Size	Product Part No.	inches	mm
MXFIT0805	2	MXPDO02	2-1/2"	64 mm
MXFIT0806	3	MXPDO03	2-3/4"	70 mm
MXFIT0807	4	MXPDO04	3"	76 mm
MXFIT0808	5	MXPDO05	3-1/4"	83 mm

DESCRIPTION

The Milex® Donut Pessary is a medical device made from silicone that is inserted into the vagina to function as a supportive structure of the uterus, bladder and/or rectum.

WARNINGS

- Do not leave pessaries in place for long periods of time, as serious complications can occur which may require surgical intervention.
- Do not use these pessaries on a patient with a known silicone allergy.
- Chemicals in various vaginal preparations can interact with the pessary material, resulting in discoloration or deterioration of the pessary. TRIMO-SAN™* does not interact with the pessary material.
- Use only TRIMO-SAN*. Other materials have not been tested for compatibility.

INDICATIONS FOR USE

For effective support of third-degree prolapse or procidentia.

CONTRAINDICATIONS

- The presence of pelvic infections and/or lacerations
- Sexually active patient
- A noncompliant patient
- Endometriosis
- Pregnant patient

INSTRUCTIONS

Review these instructions with the patient to establish use regimen.

1. Wear dry gloves. When necessary, the entering end of the pessary can be coated with TRIMO-SAN* or equivalent vaginal jelly.
2. Use one finger to depress the perineum. Hold the DONUT almost parallel to the introitus. Guide the pessary into the vaginal vault using corkscrew motion (see Figure 1).
3. The cervix should rest behind the DONUT (see Figure 2).


Figure 1

Figure 2

* TRIMO-SAN is not sold outside the USA.

4. Ask the patient to sit, stand and bear down slightly. If there is no leakage, and the patient is comfortable with the pessary in position, have her empty her bladder. A properly fitted pessary takes up slack in redundant tissue, holding the uterus higher in the vagina.
5. If the patient can void without difficulty, pessary remains in position upon re-examination, and the patient is comfortable with the pessary in place, this is a good indication that the correct size may have been selected. Patient experience may vary.

Note: If the patient is unable to urinate with the pessary in position, remove it and fit her with the next smaller size. Repeat as necessary.

6. Examine the patient while she is in the standing position to ensure the pessary has not shifted position. The patient should not feel the pessary once it is in position. The pessary should not be too loose as it may turn or be expelled and it should not be too tight as it may cause discomfort.
7. The healthcare professional should be able to sweep one finger between the pessary and the vaginal walls. If there is not enough space to do this, the next smaller size should be tried. If excessive space exists, the pessary will not be effective and may rotate or even be expelled.
8. It is sometimes necessary to refit the patient with a different size or type of pessary after a period of time. Do not assume that a replacement will always be the same size as the previous one. Check the fitting to ensure continued patient comfort and relief of symptoms. The useful life of a pessary is limited. Examine frequently for signs of deterioration (such as cracks or breaks in silicone outer surface). A pessary should be replaced if damaged.

TO REMOVE

Hook one finger inside the center of the pessary. Use the thumb and middle finger to compress the side of the DONUT (see Figure 3). Use the finger of other hand to press down on the perineum. Angle the pessary and pull gently through the introitus.

Figure 3

GELLHORN Multiple Drain Fitting Pessary

Support for third-degree prolapse/procidentia (Stage III or Stage IV)

Fitting pessary sizes available in kit and corresponding pessary sizes

GELLHORN (Flexible)				
Fitting Kit Replacement Pessary Part No.	Size	Product Part No.	inches	mm
MXFIT0012	2-1/4"	MXPGE2-1/4	2-1/4"	57 mm
MXFIT0013	2-1/2"	MXPGE2-1/2	2-1/2"	64 mm
MXFIT0014	2-3/4"	MXPGE2-3/4	2-3/4"	70 mm
MXFIT0015	3"	MXPGE3-	3"	76 mm

DESCRIPTION

The Milex® GELLHORN Pessary is a medical device made from silicone that is inserted into the vagina to function as a supportive structure of the uterus, bladder and/or rectum.

WARNINGS

- Do not leave pessaries in place for long periods of time, as serious complications can occur which may require surgical intervention.
- Do not use these pessaries on a patient with a known silicone allergy.
- Chemicals in various vaginal preparations can interact with the pessary material, resulting in discoloration or deterioration of the pessary. TRIMO-SAN™* does not interact with the pessary material. Use only TRIMO-SAN*. Other materials have not been tested for compatibility.

INDICATIONS FOR USE

For effective support of third-degree prolapse or procidentia. The cervix rests behind the flat base of the pessary and only the stem shows in the vaginal entrance when the patient does a Valsalva maneuver. The GELLHORN pessary requires a relatively capacious vagina and an intact perineum.

CONTRAINDICATIONS

- The presence of pelvic infections and/or lacerations
- A noncompliant patient
- Endometriosis
- Pregnant patient
- Sexually active patient

INSTRUCTIONS

Review these instructions with the patient to establish use regimen.

1. Wear dry gloves. When necessary, lubricate only the entering end of the pessary with TRIMO-SAN* or equivalent vaginal jelly. Hold the pessary as shown in Figure 1.
2. Use one finger to depress the perineum.
3. Guide the pessary, inserting it edgewise almost parallel to the introitus (see Figure 1), avoiding the urethral opening while the perineum is strongly pushed downward. Use a corkscrew motion while introducing the GELLHORN into the vagina.
4. Once the large flat disc is past the introitus, push the pessary upward until only the end of the stem shows in the vaginal entrance. The cervix rests behind the flat disc (see Figure 2).
5. Have the patient sit, stand and bear down. Examine the patient while she is in the standing position to ensure the pessary has not shifted position. The patient should not feel the pessary once it is in position. The pessary should not be too loose as it may turn or be expelled and it should not be too tight as it may cause discomfort.
6. The healthcare professional should be able to sweep one finger between the pessary and vaginal walls. If there is not enough space to do this, the next smaller size should be tried. If excessive space exists, the pessary will not be effective and may rotate or even be expelled.
7. It may be necessary to refit the patient with a different size or type of pessary after a period of time. Do not assume that a replacement will always be the same size as the previous one. Check the fitting to ensure continued patient comfort and relief of symptoms. The useful life of a pessary is limited.


Figure 1


Figure 2

* TRIMO-SAN is not sold outside the USA.

Examine frequently for signs of deterioration (such as cracks or breaks in silicone outer surface). A pessary should be replaced if damaged.

Note: If the patient is unable to urinate with the pessary in position, remove it and fit her with the next smaller size. Repeat as necessary.

8. If the patient can void without difficulty, pessary remains in position upon re-examination, and the patient is comfortable with the pessary in place, this is a good indication that the correct size may have been selected. Patient experience may vary.

TO REMOVE

1. Use one finger to depress the perineum.
2. Use other hand to grasp the "knob," pulling the pessary away from the cervix, turning the pessary so that the disc is almost parallel to the introitus. Using corkscrew motion ease the pessary out. See Figure 3.
3. Removal may be facilitated by passing a finger along the stem and behind the disk and folding in against the stem for removal.

Figure 3

CUBE Flexible Fitting Pessary

Support for Stage III prolapse/procidentia

Fitting pessary sizes available in kit and corresponding pessary sizes

CUBE				
Fitting Kit Replacement Pessary Part No.	Size	Product Part No.	inches	mm
MXFIT0809	2	MXPEC02	1-3/8"	35 mm
MXFIT0810	3	MXPEC03	1-1/2"	38 mm
MXFIT0811	4	MXPEC04	1-5/8"	41 mm
MXFIT0812	5	MXPEC05	1-3/4"	44 mm

DESCRIPTION

The Millex® CUBE Pessary is a medical device made from silicone that is inserted into the vagina to function as a supportive structure of the uterus, bladder and/or rectum.

WARNINGS

- Do not leave pessaries in place for long periods of time, as serious complications can occur which may require surgical intervention.
- Do not use these pessaries on a patient with a known silicone allergy.
- Chemicals in various vaginal preparations can interact with the pessary material, resulting in discoloration or deterioration of the pessary. TRIMO-SAN™* does not interact with the pessary material. Use only TRIMO-SAN*. Other materials have not been tested for compatibility.

INDICATIONS FOR USE

- The CUBE pessary is effective in the support of Stage III prolapse. The CUBE pessary will give additional support in those cases associated with complete lack of vaginal tonicity. The CUBE pessary is also effective in patients having procidentia combined with rectocele and/or cystocele.
- The CUBE pessary is also indicated in vaginal wall prolapse. The unique support is due to the suction action of the six concavities, which, by a slight negative pressure, supports the prolapse and the vaginal walls.
- The CUBE pessary should be removed nightly and washed with mild soap and water, then thoroughly rinsed. If possible, leave out overnight.

CONTRAINDICATIONS

- The presence of pelvic infections and/or lacerations
- A noncompliant patient
- Endometriosis
- Pregnant patient

INSTRUCTIONS

Review these instructions with the patient to establish use regimen.

1. Wear dry gloves. The pessary is compressed between the thumb and forefingers (see Figure 1). If necessary, the entering edge can be coated with TRIMO-SAN* or equivalent vaginal jelly.
2. Use the fingers of the other hand to spread the labia and guide the pessary into the vagina. Check the fit (see Figure 2).
3. Have the patient sit, stand and bear down. Examine the patient while she is in the standing position to ensure pessary has not shifted position. The patient should not feel the pessary once it is in position.
4. Have the patient (or caregiver) insert and remove the CUBE pessary several times before the patient leaves your office.

Figure 1

5. It is sometimes necessary to refit the patient with a different size or type of pessary after a period of time. Do not assume that a replacement will always be the same size as the previous one. Check the fitting to ensure continued patient comfort and relief of symptoms. The useful life of a pessary is limited. Examine frequently for signs of deterioration (such as cracks or breaks in silicone outer surface). A pessary should be replaced if damaged. The pessary should not be too loose as it may turn or be expelled and it should not be too tight as it may cause discomfort.

Figure 2

Note: If the patient is unable to urinate with the pessary in position, remove it and fit her with the next smaller size. Repeat as necessary.

6. If the patient can void without difficulty, and the pessary remains in position upon re-examination, and the patient is comfortable with the pessary in place, this is a good indication that the correct size may have been selected. Patient experience may vary.
7. Examine the patient while she is in the standing position to ensure the pessary has not shifted position. The patient should not feel the pessary once it is in position. The pessary should not be too loose as it may turn or be expelled and it should not be too tight as it may cause discomfort.

* TRIMO-SAN is not sold outside the USA.

8. The healthcare professional should be able to sweep one finger between the pessary and the vaginal walls. If there is not enough space to do this, the next smaller size should be tried. If excessive space exists, the pessary will not be effective and may rotate or even be expelled.

TO REMOVE

DO NOT PULL ON CORD TO REMOVE PESSARY!

The cord is to help locate the CUBE pessary. Pulling on the cord to remove the CUBE can damage vaginal mucosa. To remove the pessary, the suction with the vaginal walls must be broken. Work the finger tips between the vaginal mucosa and the pessary, compress the CUBE between the thumb and forefinger and remove (see Figure 1).

Keep the plunger pressed until the Universal Applicator is removed from the vagina.

NOTE: For those wearing an INFLATOBALL or CUBE pessary, remove the pessary per the physician's instructions before using TRIMO-SAN.

6. Clean the Universal Applicator immediately after using by flushing the Universal Applicator (pull and press plunger back and forth) several times, holding the Universal Applicator under water (see Figure 3).
Make sure that the Universal Applicator is clean and dry before reusing.
7. After each use, wipe the threads of the tube and inside the cap. Recap tightly.

TRIMO-SAN™ INFORMATION*

(TRIMO-SAN is not sold outside the USA)

Note: Advise your patient that TRIMO-SAN* is available at the prescription counter upon request (no Rx required).

WARNING

Chemicals in various vaginal preparations can interact with the pessary material, resulting in discoloration or deterioration of the pessary. TRIMO-SAN* does not interact with the pessary material. Use only TRIMO-SAN*; other materials have not been tested for compatibility.

PURPOSE OF TRIMO-SAN*

1. To help restore and maintain normal vaginal acidity
2. To coat the walls of the vagina with a lubricating film that helps reduce odor-causing bacteria

APPLICATION

- Pessary wearers should use 1/2 applicator (2g) of TRIMO-SAN* three (3) times the first week after original insertion of the pessary, unless otherwise directed by the healthcare professional
- Use 1/2 applicator (2g) of TRIMO-SAN* twice a week thereafter, unless otherwise directed by the healthcare professional
- Reaction or irritation caused by TRIMO-SAN* is very rare, but should it occur, discontinue use
- If symptoms persist or worsen, contact your healthcare professional

TRIMO-SAN* HAS NOT BEEN TESTED DURING PREGNANCY

- If the patient wishes to become pregnant, please advise her against using TRIMO-SAN* for 6 hours prior to or following intercourse

PROPER USE OF THE UNIVERSAL APPLICATOR

1. Unscrew the cap from the tube.
2. Pierce the end of the tube with top of cap (a small amount of TRIMO-SAN may escape when the top of the tube is pierced).
3. Attach the threaded end of the Universal Applicator to the TRIMO-SAN tube (see Figure 1).
4. While holding the Universal Applicator, squeeze the tube from the bottom until TRIMO-SAN reaches the specified dosage level (see Figure 2).
5. Detach the Universal Applicator from the tube. Gently insert the Universal Applicator approximately 1/2 inch into the vagina. Press in the plunger of the Universal Applicator to dispense the gel.

Figure 1


Figure 2

Figure 3

* TRIMO-SAN is not sold outside the USA.

EXPLANATION OF SYMBOLS

	Reorder number
	Batch code
	Use-by date
	Consult instructions for use
	Caution
	Non-Sterile
	Do not use if package is damaged or opened
	Not made with natural rubber latex
Rx Only	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician
	Authorized Representative in the European Community.
	Manufacturer

Dawn® is a registered trademark of Procter & Gamble, Cincinnati, Ohio 45202
CIDEX® OPA is a registered trademark of Johnson & Johnson.
Revital-Ox™ RESERT® are trademarks of STERIS Corporation.
TRIMO-SAN™ is a trademark of CooperSurgical, Inc.
Milex® and CooperSurgical are registered trademarks of CooperSurgical, Inc.

Assembled in the USA

* TRIMO-SAN is not sold outside the USA.

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