RUMI® II Koh-Efficient™

Vaginal Fornices Delineator & Colpo-Pneumo Occluder™ for use with the RUMI® II Handle

Koh-Efficient™ for Electrocautery

Koh-Efficient™ for Harmonic Scalpel and Laser

2.5 cm • 3.0 cm • 3.5 cm • 4.0 cm

Instructions For Use (English)

CE 2797

STERILE

SINGLE-USE

DEVICE DESCRIPTION

The RUMI® II Koh-Efficient™ is a sterile, single-use vaginal fornices delineator with an integrated Colpo-Pneumo Occluder™. It is made of medical grade plastics with an Ultem® Polymeric rigid KOH Cup™. The Harmonic Scalpel/laser compatible version (KCS-RUMI) is equipped with a stainless steel distal surface. Not made with natural rubber latex.

The RUMI II Koh-Efficient (Koh-Efficient) is compatible with the RUMI® II Handle (Handle) and the RUMI® Uterine Manipulator Tip (Tip). These three products, used in conjunction with the Cervical Sizer (Sizer), enable unobstructed placement of a standard Tip into the uterus. Once in place, the Koh-Efficient is advanced along the Arm of the Handle and locks into place around the cervix.

RUMI II Handle
(UMH650)

Cervical Sizer
(KC-SIZER)

**RUMI
Uterine
Manipulator Tip**

**RUMI II
Koh-Efficient**
(KC-RUMI or
KCS-RUMI)

WARNINGS

- Contents supplied sterile. Do not use if sterile barrier is damaged.
- The integrated Colpo-Pneumo Occluder is designed to accommodate no more than 180 cc of inflation.
- The Koh-Efficient is single-use only, and should NEVER be reused.
- Using an incorrect size KOH Cup could result in patient injury.
- The Ultem KOH Cup is not intended to be used with Harmonic or other types of ultrasonic or laser energy sources.
- The stainless steel KOH Cup is not intended to be used with electrocautery. Use of electrocautery may result in inadvertent injury to adjacent structures.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not

limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Dispose of in accordance with all applicable Federal, State, and local Medical/Hazardous waste practices.

CAUTION

- U.S. Federal law restricts this device to sale by or on the order of a physician.
- The RUMI II Koh-Efficient is to be used ONLY with the RUMI II Handle (UMH650).

INTENDED USE/INDICATIONS FOR USE

The Koh-Efficient is a combination KOH Cup and Colpo-Pneumo Occluder that is indicated for use in laparoscopic procedures where it is desirable to locate the position of the vaginal fornices, typically during a total laparoscopic hysterectomy or supra-cervical laparoscopic hysterectomy.

CONTRAINDICATIONS


The Koh-Efficient should not be used in patients who are pregnant or who are suspected of being pregnant, in patients that have an IUD in place, in patients with suspected pelvic infections or in cases where the surgeon deems it inadvisable and finds it difficult to insert the silicone tip into the cervix or uterus.

PRECAUTIONS

Always confirm that the Koh-Efficient is locked in place prior to starting the colpotomy.

INSTRUCTIONS FOR USE

RUMI II System (Assembled)

For complete information on how to prepare The RUMI® II System, please refer to the Instructions for Use for the RUMI II Handle (UMH650), the RUMI Uterine Manipulator Tips and the Cervical Sizer (KC-SIZER).

1. Prepare the Handle for use by following the Handle's Instructions for Use.

2. Attach the proper size Tip to the Handle. **Fully press the Tip Catheters into the three channel positions on the Arm** (see Figure 1).

3. Using the Sizer, determine the proper size Koh-Efficient.


Figure 1

4. With the Handle in the neutral position (see Figure 2), and the Koh-Efficient in the open position, slide the Koh-Efficient over the Tip and Arm (see Figure 3) to where the arrow on the Koh-Efficient meets the arrow on the Arm (see Figure 4).

Figure 2


Figure 3


Figure 4

5. Close the Koh-Efficient (see Figure 3) by firmly pressing it onto the Arm (see Figure 5). Slide the Koh-Efficient back to the start position (see Figure 6).


Figure 5


Figure 6

6. Test inflate the Colpo-Pneumo Occluder with 20 to 60 cc of sterile saline, then deflate prior to patient insertion.
7. Insert the Tip into the uterus and inflate as per Tip Instructions for Use.
8. Press on the Push Pad and slide the Koh-Efficient forward until locked in place (see Figure 7) and a “click” is heard. Confirm locked position by pulling back on the Push Pad. Visually confirm that the indicator marks on the Koh-Efficient are lined up with the indicator line on the Arm (see Figures 7 & 8) and the Push Pad is in line with the end of the arrow (see Figure 8). *Note: The last 1/4” of travel will have increased resistance.*


Figure 7


Figure 8

9. Confirm that the distal end of the cup resides in the fornices.
- Note: If it is necessary to move the Koh-Efficient, push the release tab to disengage fingers from Arm (see Figure 9) and slide Koh-Efficient back.
10. The Tip can now be articulated by pressing the trigger on the Grip and rotating the Grip.
- When the Trigger is released, the Tip will be locked in position.
11. Just prior to making colpotomy incisions, inflate the Colpo-Pneumo Occluder balloon with 60 to 120 cc of sterile saline or water (do not use air).
12. When it is time to remove the uterine manipulator from the patient, you must first deflate the Colpo-Pneumo Occluder balloon.
13. Remove the Koh-Efficient by spreading outward and down the Removal Flats (see Figure 9).


Figure 9

14. Dispose of in accordance with all applicable Federal, State and local Medical/Hazardous waste practices.

EXPLANATION OF SYMBOLS

Reorder Number

Batch Code

Use-by date

Ethylene Oxide Sterilization

Do not resterilize

Do not use if package is damaged

Do not re-use

Caution

Consult instructions for use

Not made with natural rubber latex

Manufacturer

Rx Only

U.S. Federal law restricts this device to sale by or on the order of a physician.

Authorized Representative in the European Community

The marks bearing the symbols “TM” and “®” are trademarks and registered trademarks of CooperSurgical, Inc. Ultem® is a registered trademark of SABIC Innovative Plastics CooperSurgical is a registered trademark of CooperSurgical, Inc.

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© 2020 CooperSurgical, Inc.
As marcas que apresentam os símbolos “TM” e “®” são marcas comerciais e marcas registadas da CooperSurgical, Inc.
CooperSurgical é uma marca registada da Coopersurgical, Inc.