



CerviPATH™

Patient Management Systems

Visualization and Documentation at its Best

Features and Benefits

- Image capture, annotation and documentation for tracking lesion regression/progression
- Create, store and export images and reports
- Live video feed for educational purposes
- Added security with password protection and data hide
- Available exclusively on two LEISEGANG OptiK™ Model 2 colposcope models, the colposcopy "Gold Standard"



CooperSurgical

Keeping You at the Forefront
of Women's Health Care™

Perform the most comprehensive exam possible with CerviPATH™ Digital Colposcopy Systems

document, compare and track with CerviPATH™

Specifications

- LEISEGANG OptiK™ Model 2 Colposcopes
 - Parfocal three-step drum magnification (3.75X, 7.5X and 15X)
 - Field of view 19 mm to 76 mm
 - Depth of field 2 mm to 8 mm
 - Bright white LED illumination system
 - Individually adjusting diopters for correction of ametropia
 - Variable gross and fine adjustments for optimum height and focus
 - Green filter
- Laptop, tray, and infrared printer
- CerviPATH™ Patient Management software
- 1.3 Megapixel Infinity 1.1 camera or integrated camera
- Available in either Swing or Balance-O-Matic base
- LCD monitor accessory available
- Complimentary system setup and in-service
- One-year standard warranty. Additional service options available under the Customer Assurance Program (CAP)

Documented regression rates of up to 60% for CIN I and up to 40% for CIN II are a fact.* Clinical best practices suggest patient monitoring with documentation for optimal care. CerviPATH™ allows you to:

- Document findings and capture images
- Annotate images with text, measurements and drawings
- Compare images within current visit or against previous visit
- Visually track the cervical health of your patients

CerviPATH™ is only available with LEISEGANG OptiK Model 2 photo colposcopes. Select from either a fully integrated system or the external camera variation. Both offer high quality photography for your documentation and the performance you've come to expect from LEISEGANG.

Features

- Specially selected optics in a convergent path for pristine 3D visualization
- High resolution image capture and live video feed
- Annotate images with text, measurement, voice and drawing tools
- Compare images from current and previous visits
- Track lesion progression and regression
- Client/server capable
- Export image (jpegs) and report (PDFs) to disk, e-mail or network
- Detailed patient demographic fields
- Comprehensive exam documentation, including findings and follow-up
- Single-image and multiple-image print options
- Four comprehensive report templates
- Patient image slide show
- Record 15-60 second video clips
- Multiple file cabinets for patient record storage

To find out more about the CerviPATH™ Digital Colposcopy Systems, or to schedule a demonstration, call 800.243.2974, or visit www.coopersurgical.com

*Adapted from Arends, J, et al, *J Clin Pathol*. 1998;51:96, Holowaty, P, et al, *J Nat Cancer Inst*. 1999;91:252, McIndoe, WA, et al, *Obstet Gynecol*. 1984;64:451, and McInikow, J, et al, *Obstet Gynecol*. 1998;92:727.

*Adapted from Arends, J, et al, *J Clin Pathol*. 1998;51:96, Holowaty, P, et al, *J Nat Cancer Inst*. 1999;91:252, McIndoe, WA, et al, *Obstet Gynecol*. 1984;64:451, and McInikow, J, et al, *Obstet Gynecol*. 1998;92:727.

Ordering Information

DESCRIPTION	ORDER NUMBER	BASE CONFIGURATION
OptiK Model 2 0.5" CCD Optimized with 1.1 MP Camera and CerviPATH™ Patient Management System	OPTIK2-06-004	Swing Base
	OPTIK2-07-004	Balance-O-Matic Base
OptiK Model 2 Integrated 0.5 MP Camera and CerviPATH™ Patient Management System	OPTIK2-13-003	Swing Base
	OPTIK2-14-003†	Tilt Base
	OPTIK2-15-003†	Upright Base
CerviPATH Upgrade: Software, laptop, Tray, Infinity Camera Video Adapter	53321 08706	

† Does not include laptop tray