



Endometrial Samplers

The best options for the best results

Products intended for endometrial tissue sampling

- SelectCells Mini®
- Pipet Curet
- Uterine Explora Model I

Features and benefits

- Safe
- Sterile/Disposable
- Provides reliable results
- Latex free



@operSurgical

**Keeping You at the Forefront
of Women's Health Care™**

CooperSurgical provides you with the best product options designed to maximize histological sample retrieval

SelectCells Mini® - SM203

Experience the full-size advantage of a small-diameter endometrial sampler

- Smallest diameter of any endometrial sampling device
- Narrow profile for greater patient comfort
- Eases insertion through stenotic os
- Sterile/Disposable (One-time Use)



Milex™ Pipet Curet - MX140

Doctor preferred for smoother sampling

- Thicker stylet with double "O-Ring" maximizes suction for a smoother pull
- Centimeter markings clearly marked for better visibility
- Larger opening facilitates maximum tissue retrieval
- Sterile/Disposable (One-time Use)



Uterine Explora Model I - MX120

Utilizes Syringe Technique for exceptional control and recovery

- Nylon curette for better patient comfort
- Vacu-Lok Syringe for one-handed operation
- Randall-type cutting edge allows for superior sampling
- Sterile/Disposable (One-time Use)



Product Specifications	Milex™ Pipet Curet	SelectCells Mini®	Explora Model I
Part Number	MX140	SM203	MX120
Outer Diameter	3 mm	1.95 mm	3 mm
Length	23.97 cm	25.6 cm	19.69 cm
Suction Technique	Piston	Piston	Syringe
Key Features	<ul style="list-style-type: none"> • Double O-Ring • Distinct Centimeter Markings 	<ul style="list-style-type: none"> • Centimeter Markings • Non-traumatic Design 	<ul style="list-style-type: none"> • Vacu-Lok Syringe • Nylon Curette



www.coopersurgical.com

95 Corporate Drive, Trumbull, CT 06611
800.243.2974 • 203.601.5200 • Fax 800.262.0105

Form # 81257 03/08
© Copyright 2008 CooperSurgical, Inc.