

## Uterine Manipujector



## The Kronner Manipujector®

Intra-operatively Adjustable Uterine Manipulator-Injector

### Features and Benefits

- Affords rapid intra-operative repositioning
- Provides excellent manipulative control
- Securely seals internal os for no-leak injection
- Safe, simple and economical

Designed For Maximum  
Procedure Versatility




**CoperSurgical**

*Keeping You at the Forefront  
of Women's Health Care™*

The latest advance in balloon-tipped uterine manipulator-injectors. Designed for patient safety, ease-of-use, convenience and maximum procedural versatility.

### **Unique Pinch-Grip Handle - For Rapid, Effortless Intra-Operative Adjustment**

The Kronner Manipujector's unique glide action, pinch-grip handle is designed to eliminate the time consuming tasks required to remove, reinsert and reinflate when intra-operative repositioning is required. Now, just pinch the grip and glide the handle to the new position. Then release the grip, which sets the handle firmly, and you're ready to continue your procedure.



### **Anatomically Curved Shape - Provides Excellent Manipulative Control**

Radial design optimizes cervical visualization, allowing easy insertion and maneuverability of either anteverted or retroverted uteri. Curvature affords superb elevation for ease of surgical dissection and visualization.

### **Positive Balloon Positioning - Reduces Potential For Expulsion**

Elimination of spring-controlled automatic retraction obviates the need to maintain forward pressure during balloon inflation, and helps avoid inadvertent balloon inflation within the cervical canal – a common cause of spontaneous balloon expulsion requiring device reinsertion.

***The Kronner Manipujector is safe, simple and economical to use. It is designed to afford maximum operator convenience and minimal uterine trauma.***

**CooperSurgical**

[www.coopersurgical.com](http://www.coopersurgical.com)

To find out more about the Kronner Manipujector®, or to schedule a demonstration call 800.243.2974.

US Patent No. 4,775,362  
Kronner Manipujector® is a trademark of CooperSurgical Inc.

**Kronner Manipujector®**

95 Corporate Drive, Trumbull, CT 06611  
800.243.2974 • 203.601.5200  
International:  
203.601.9818 • fax 203.601.4747

Form # 81357 02/08  
© Copyright 2008 CooperSurgical, Inc.