


# INTRODUCING Mystic<sup>®</sup> II



Enhancing a winner —  
The Mityvac line  
of vacuum assisted  
delivery systems

*Remember, one size cup does  
not fit all deliveries!*

## Safety

Lower scalp lacerations by 68%

## Low Profile

Mushroom<sup>®</sup> Cup flexes up to 90° in any  
direction for easier insertion

## Superior Control

Handle can rotate 180° without transferring  
rotation to the cup

## Versatile

An array of cups to choose from

CooperSurgical

95 Corporate Drive, Trumbull, CT 06611 • 1(800)MITYVAC • [www.coopersurgical.com](http://www.coopersurgical.com)  
©2007 CooperSurgical, Inc.

Mystic<sup>®</sup> is a product of CooperSurgical, Inc.  
Form # 81237OBGYN 06/07