

Pessary Material Information Sheet

Product Code	Product Description	Size Range	Silicone	Latex	Silicone and Metal	Other
PCOND PCONDS	Incontinence Dish with Support without Support	55 mm — 85 mm				
PDO	Donut	2" — 2 ¾"				
PEC	Cube	0 — 7				
PECH	Cube with Drainage	0 — 7				
PECT	Tandem Cube	2/0 — 7/5				
PGE	Silicone Flexible Gellhorn	1 ½" — 3 ½"				
PGRS	95% Rigid Silicone Gellhorn	1 ½" — 3 ½"				
PINF	Inflatoball	S — XL		Latex		
PSH	Shatz	1 ½" — 3 ½"				
PRS	Ring with Support	0 — 10	With nylon pegs			
PRS	Ring with Support	11, 12, 13			Remove before MRI*	
PER	Ring	0 — 10	With nylon pegs			
PER	Ring	11, 12, 13			Remove before MRI*	
PRSK	Ring with Support and Knob	0 — 8	With nylon pegs			
	Ring with Support and Knob	9, 10, 11, 12, 13			Remove before MRI*	
PERK	Ring with Knob	0 — 8	With nylon pegs			
PERK	Ring with Knob	9, 10, 11, 12, 13			Remove before MRI*	
PCON	Incontinence Ring	0 — 13			Remove before MRI*	
PGS	Gehrung	1 — 9			Remove before MRI*	
PGSK	Gehrung with Knob	1 — 9			Remove before MRI*	
PREG	Regula	2 — 8			Remove before MRI*	
PEH PHS PHSK PEHK	Hodge Hodge with Support Hodge with Support and Knob Hodge with Knob	0 — 9			Remove before MRI*	
PRI	Risser	0 — 9			Remove before MRI*	
PES	Smith	0 — 8			Remove before MRI*	

* Removal before x-ray may be warranted if pessary will obscure the target area.



CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products which aid clinicians in the management and treatment of commonly seen conditions.

203.601.5200 | 800.243.2974 | coopersurgical.com | 95 Corporate Dr. | Trumbull, CT 06611

©2011 CooperSurgical, Inc.

Form #81474 • Rev. 05/11