




# Lone Star **APS** Retractor™

See the difference unobstructed transvaginal surgery can make

## For Obstetric and Gynecologic Procedures

### Benefits

- Unobstructed visualization and access
- Designed specifically for modern pelvic and vaginal surgery
- Effectively manages the surgical site, simplifying the procedure
- VELCRO® backing for increased stability and ease of placement



**Better exposure.**  
**Better access.**  
**Better all around.**

# The unique Lone Star APS Retractor

Unprecedented access - for surgical excellence


The versatile APS system provides unparalleled exposure. For simpler, better managed surgeries, the innovative design of the self-retaining APS Retractor™ offers a clear pathway, unobstructed anatomical visibility, and dramatic improvement in site management with minimum effort. Adjustable clips, stable VELCRO® backing and ultra-precise hook placements quickly adapt the APS Retractor to present the surgical site exactly as needed. The system frees up OR personnel, saving time and operating expense. There's nothing like it – and nothing that performs as well.

## The APS Retractor Designed specifically for:

- Reconstructive pelvic surgery
  - Sling procedures
  - Anterior repair procedures
  - Posterior repair procedures
- Mesh excision
- Episiotomy repair
- Vaginal hysterectomy
- Vaginal biopsies

## Achieve unparalleled visualization and access...

*The APS Retractor combined with the recommended Stays provide excellent visual exposure of the vagina.*



- Convenient clip manages the Foley catheter
- VELCRO backing facilitates site set up and management
- Closely placed slots allow greater precision in hook/stay placement
- Custom holes for non-perforating towel clamps keep retractor stable and stationary
- Ideal width for a variety of instrument introduction
- Unique, open-end design facilitates use of weighted speculum, or access to the posterior compartment

## Specifications

### APS Retractor

(Packaged singly with 2 catheter clips)

Order # 3715	Recommended Stays	
Height 16.6 cm	<b>3311</b>	5 mm Sharp Hook
Width 16.6 cm	<b>3327</b>	7 mm Two-finger Sharp Hook
Weight 77 g	<b>3350</b>	2 mm Blunt Hook
	<b>3384</b>	12 mm Two-finger Blunt Hook

Products are sterile, single-use only



coopersurgical.com • also on YouTube at youtube.com/coopersurgical

95 Corporate Drive, Trumbull, CT 06611  
800.243.2974 • 203.601.5200

To find out more about the Lone Star APS Retractor or to place an order call 800-243-2974.

The marks bearing the symbols "TM" and "®" are trademarks and registered trademarks of CooperSurgical, Inc. VELCRO® is a registered trademark of Velcro Industries, B.V.

Form # 81710 • Rev. 03/10  
© Copyright 2010 CooperSurgical, Inc