



Quality and Innovation to Support Today's Surgeons

# Surgical Products Catalog







Quality and innovation to support

# today's surgeon



CooperSurgical®



Our portfolio of high-quality surgical and laparoscopic instruments was founded on the specialized needs of women's surgery, yet many of our surgical products are found throughout the hospital. Our products include devices for hysterectomy, pelvic surgery and C-section procedures, as well as surgical retractors and instruments for wound and port site closure.



# Surgical Brands

At CooperSurgical, we develop and market innovative solutions for surgical procedures. Within this catalog, you will find an extensive range of instruments and devices that facilitate surgery and improve patient care.



ALLY Uterine Positioning System®

Advincula Delineator™ with Koh-Efficient® Technology

Koh-Efficient®

RUMI® II System

Advincula Arch®

Kronner Manipujector®

INSORB® Subcuticular Skin Stapler

Carter-Thomason CloseSure System®

SeeClear® Laparoscopic Smoke Evacuation System

Lone Star® Retractor System

OR/Mobius®

# Table of Contents

**Uterine Manipulation Devices**  
Gynecologic Laparoscopic Product Selection ..... 6–7  
ALLY Uterine Positioning System® ..... 8  
Advincula Delineator™ ..... 9  
RUMI® II System ..... 10  
RUMI® II Long ..... 10  
RUMI® II / Koh-Efficient® System ..... 11  
Cervical Sizer ..... 11  
Advincula Arch® ..... 12–13  
Sacrocolpopexy and Sacrocervicopexy Tips ..... 14  
Hoyte Sacro Tips ..... 14  
Disposable Diagnostic Uterine Manipulators ..... 15

**Subcuticular Skin Stapler** ..... 16

**Laparoscopic Port Site Closure**  
Carter-Thomason CloseSure System® ..... 17  
Carter-Thomason® II Port Closure System ..... 18  
SureGrasp™ Suture Passer ..... 19

**Laparoscopic Smoke Evacuation**  
SeeClear® Smoke Evacuation System ..... 20

**Trocars** ..... 21–22

**Retractors**  
Lone Star® Retractor System ..... 23–24  
OR/Mobius® ..... 25

**Hysterectomy Instruments**  
Zeppelin® Hysterectomy Clamps ..... 26  
Zeppelin® Vaginal Hysterectomy Clamps ..... 27  
Z-Clamp™ ..... 27  
SuperCut Scissors ..... 28  
Z-Scissors™ ..... 28

**Laparoscopic Instruments** ..... 29

**Laparoscopic Disposables** ..... 29

**Index** ..... 30

**Warranty Information** ..... 31

**Customer Service** ..... 31

CE Available for sale where CE mark is required

CA Available for sale in Canada



# Gynecologic Laparoscopic Product Selection



Scan the QR Code to Access the In-Service App

## INDICATIONS

### Uterine Tips




Gynecologic Laparoscopic Procedures and Diagnostics

### Koh-Efficient®




Gynecologic Laparoscopic Procedures with RUMI II

### Koh-Efficient



Gynecologic Laparoscopic Procedures with Advincula Arch

### Sacrocolpopexy Tips



Vaginal Prolapse Repair Without Cervix

### Sacro cervicopexy Tips



Vaginal Prolapse Repair With Cervix

### Hoyte Sacrocolpopexy Tips



Vaginal Prolapse Repair Without Cervix

### Hoyte Sacro cervicopexy Tips



Vaginal Prolapse Repair With Cervix

### Disposable Manipulator With Integrated Koh-Efficient



TLH/LSH/LAVH

### Manipujector®



Diagnostics



# ALLY Uterine Positioning System®

The ALLY UPS, hands-free control for a full range of gynecologic laparoscopic procedures

## ALLY Uterine Positioning System

NUMBER	DESCRIPTION	
AU-UPS	ALLY Uterine Positioning System	
AU-CART	ALLY Uterine Positioning System Cart	€€
AU-AD	Adapter Drape for RUMI II and Advincula Arch (20/Box)	€€
AU-AD-DLNTR	Adapter Drape for Advincula Delineator (20/Box)	€€
AU-UPS-EXT	Extended Power Cable, 20 ft	

## International ALLY Uterine Positioning System

NUMBER	DESCRIPTION	
AU-INT	ALLY UPS, International	€€
AU-KITF	AU Accessory Kit F	€€
AU-KITG	AU Accessory Kit G	€€
AU-KITH	AU Accessory Kit H	
AU-KITI	AU Accessory Kit I	
AU-KITJ	AU Accessory Kit J	
AU-KITM	AU Accessory Kit M	

\*Kit contains power cord, plug and rail clamp.  
Please consult with your territory representative for the appropriate kit.


## ALLY Uterine Positioning System Cart

- Designed to bring the ALLY UPS directly to the OR table
- Offers storage and maneuverability throughout the hospital



# Advincula Delineator™

with Koh-Efficient®  
Technology



Disposable Uterine  
Manipulator



## Advincula Delineator™ with Ultem® Koh-Efficient®

NUMBER	DESCRIPTION	
AD750-KE25	Advincula Delineator 2.5CM Ultem Plastic cup (3/Box)	CA* CE
AD750-KE30	Advincula Delineator 3.0CM Ultem Plastic cup (3/Box)	CA* CE
AD750-KE35	Advincula Delineator 3.5CM Ultem Plastic cup (3/Box)	CA* CE
AD750-KE40	Advincula Delineator 4.0CM Ultem Plastic cup (3/Box)	CA* CE

## Advincula Delineator™ with Stainless Steel Koh-Efficient®

NUMBER	DESCRIPTION	
AD750-KES25	Advincula Delineator 2.5CM Stainless Steel cup (3/Box)	CA* CE
AD750-KES30	Advincula Delineator 3.0CM Stainless Steel cup (3/Box)	CA* CE
AD750-KES35	Advincula Delineator 3.5CM Stainless Steel cup (3/Box)	CA* CE
AD750-KES40	Advincula Delineator 4.0CM Stainless Steel cup (3/Box)	CA* CE

## Advincula Delineator™ with Flexible Soft Cup Koh-Efficient®

NUMBER	DESCRIPTION	
AD750SC-KE25	Advincula Delineator 2.5CM Flexible Soft cup (3/Box)	
AD750SC-KE30	Advincula Delineator 3.0CM Flexible Soft cup (3/Box)	
AD750SC-KE35	Advincula Delineator 3.5CM Flexible Soft cup (3/Box)	
AD750SC-KE40	Advincula Delineator 4.0CM Flexible Soft cup (3/Box)	



# RUMI® II System

## A. The RUMI® II System CA★ CE

The articulating RUMI II System, with locking handle, soft distal tip, silicone balloon, and integrated catheter for chromopertubation provides uterine manipulation for gynecologic laparoscopic procedures.

NUMBER	DESCRIPTION
UMH650	RUMI II Uterine Manipulator Handle
UMH650L	RUMI II Long Handle – 2” Longer




Locking handle



Soft distal end of tip and silicone balloon

### B.



6.7 mm x 12 cm 6.7 mm x 10 cm 6.7 mm x 8 cm 6.7 mm x 6 cm 5.1 mm x 6 cm 5.1 mm x 3.75 cm

## B. The RUMI II System Tips CA★ CE

RUMI II single-use intrauterine tips are available in six sizes to accommodate individual patient anatomy. Made from flexible silicone, RUMI II tips have a soft, flexible distal end to be used with RUMI II system.

NUMBER	DESCRIPTION
UMY514	RUMI II Tip Yellow, 5.1 mm x 3.75 cm (5/Box)
UML516	RUMI II Tip Lavender, 5.1 mm x 6 cm (5/Box)
UMW676	RUMI II Tip White, 6.7 mm x 6 cm (5/Box)
UMB678	RUMI II Tip Blue, 6.7 mm x 8 cm (5/Box)
UMG670	RUMI II Tip Green, 6.7 mm x 10 cm (5/Box)
UM0672	RUMI II Tip Orange, 6.7 mm x 12 cm (5/Box)

# RUMI® II with Koh-Efficient® Technology

The Koh-Efficient system combines the Koh Cup™ and pneumo-occlusion balloon into one device. This sterile, single-use device is available in four sizes.

## C. RUMI II Koh-Efficient System CA★ CE

### For Electrosurgery

NUMBER	DESCRIPTION
KC-RUMI-25	RUMI II / Koh-Efficient, 2.5 cm (5/Box)
KC-RUMI-30	RUMI II / Koh-Efficient, 3.0 cm (5/Box)
KC-RUMI-35	RUMI II / Koh-Efficient, 3.5 cm (5/Box)
KC-RUMI-40	RUMI II / Koh-Efficient, 4.0 cm (5/Box)

### For Harmonic Scalpel/Laser

NUMBER	DESCRIPTION
KCS-RUMI-25	Stainless Steel Capped RUMI II / Koh-Efficient, 2.5 cm (5/Box)
KCS-RUMI-30	Stainless Steel Capped RUMI II / Koh-Efficient, 3.0 cm (5/Box)
KCS-RUMI-35	Stainless Steel Capped RUMI II / Koh-Efficient, 3.5 cm (5/Box)
KCS-RUMI-40	Stainless Steel Capped RUMI II / Koh-Efficient, 4.0 cm (5/Box)



Stainless Steel Capped Cup available in all four sizes

## Colpo-Pneumo Occluder™ CA★ CE

NUMBER	DESCRIPTION
CPO-6	Colpo-Pneumo Occluder, Sterile, Single-Use (6/Box)

## Koh Cups, Ultem™ Resin CA★ CE

### For Electrosurgery

NUMBER	DESCRIPTION
KCP-30-2	Ultem Koh Cup, Nonsterile, 3 cm (2/Box)
KCP-35-2	Ultem Koh Cup, Nonsterile, 3.5 cm (2/Box)
KCP-40-2	Ultem Koh Cup, Nonsterile, 4 cm (2/Box)

## Koh Cups, Stainless Steel CA★ CE

### For Harmonic Scalpel/Laser

NUMBER	DESCRIPTION
KCS-30	Stainless Steel Koh Cup, One Nonsterile, 3 cm
KCS-35	Stainless Steel Koh Cup, One Nonsterile, 3.5 cm
KCS-40	Stainless Steel Koh Cup, One Nonsterile, 4 cm

## Cervical Sizer CA★ CE

Use the Cervical Sizer to select the appropriate Koh-Efficient size based on patient anatomy.

NUMBER	DESCRIPTION
KC-Sizer	Cervical Sizer (1/Each)





# Advincula Arch®

Reusable Lightweight Uterine Manipulator

NUMBER	DESCRIPTION	CA* CE
UMH750	Advincula Arch Uterine Manipulator Handle	

## Arch Koh-Efficient® System

The Arch Koh-Efficient combines the Koh Cup™ and the Colpo-Pneumo Occluder™ balloon into one device for use with the Advincula Arch® handle. The sterile, single-use device is available in four sizes to suit your patient’s anatomy, and is offered in both stainless steel (for harmonic scalpel and laser procedures) and Ultem® (for use with electro-surgical cutting devices).

### Arch Koh-Efficient® – For Electrosurgery CA\* CE

NUMBER	DESCRIPTION
KC-ARCH-25	Ultem™ Plastic Arch Koh-Efficient, 2.5 cm (5/Box)
KC-ARCH-30	Ultem Plastic Arch Koh-Efficient, 3.0 cm (5/Box)
KC-ARCH-35	Ultem Plastic Arch Koh-Efficient, 3.5 cm (5/Box)
KC-ARCH-40	Ultem Plastic Arch Koh-Efficient, 4.0 cm (5/Box)

### Arch Koh-Efficient® – For Harmonic Scalpel/Laser CA\* CE

NUMBER	DESCRIPTION
KCS-ARCH-25	Stainless Steel Capped Arch Koh-Efficient, 2.5 cm (5/Box)
KCS-ARCH-30	Stainless Steel Capped Arch Koh-Efficient, 3.0 cm (5/Box)
KCS-ARCH-35	Stainless Steel Capped Arch Koh-Efficient, 3.5 cm (5/Box)
KCS-ARCH-40	Stainless Steel Capped Arch Koh-Efficient, 4.0 cm (5/Box)

The Advincula Arch shares the same tips as the articulating RUMI II.



Stainless Steel Capped Cup  
Available in all four sizes





## Sacrocolpopexy and Sacrocervicopexy Tips

For Vaginal Prolapse Repair



### Sacrocolpopexy / Sacrocervicopexy Tips CA★ CE

The Sacrocolpopexy Tip can be used in conjunction with any RUMI or Advincula Arch handles.

NUMBER	DESCRIPTION
<b>SACRO-1</b>	Sacrocolpopexy Tip – Large 1.5" x 4.5" (3.8 cm x 11.4 cm) (3/Box)
<b>SACRO-1S</b>	Sacrocolpopexy Tip – Small 1.3" x 3.5" (3.2 cm x 9 cm) (3/Box)
<b>SACRO-1C</b>	Sacrocolpopexy Tip – Concave Distal End 1.5" x 4.5" (3.8 cm x 11.4 cm), Cavity 1.3" x 1.5" deep (3.3 cm x 3.8 cm) (3/Box)

### Hoyte Sacro Tips CA★ CE

The Hoyte Sacro Tips can be used in conjunction with the Advincula Arch® Handle.

NUMBER	DESCRIPTION
<b>LH-SCOLPO-LRG</b>	Hoyte Sacrocolpopexy Tip – Large, Single-use (3/Box)
<b>LH-SCOLPO-SML</b>	Hoyte Sacrocolpopexy Tip – Small, Single-use (3/Box)
<b>LH-SCERV-LRG</b>	Hoyte Sacrocervicopexy Tip – Large, Single-use (3/Box)
<b>LH-SCERV-SML</b>	Hoyte Sacrocervicopexy Tip – Small, Single-use (3/Box)
<b>UMH750</b>	Advincula Arch Uterine Manipulator Handle



## Disposable Diagnostic Uterine Manipulators



### Kronner Manipujector® CA\* CE

The disposable Kronner Manipujector used for laparoscopic procedures.

NUMBER	DESCRIPTION
6003	Kronner Manipujector (12/Box)



### HUMI® Harris-Kronner Uterine Manipulator Injector

The HUMI is a curved, rigid, double-lumened device that combines the capability to manipulate the uterus and to introduce dyes, contrast media or gas without changing instruments.


NUMBER	DESCRIPTION
6001	HUMI Harris-Kronner Uterine Manipulator Injector (12/Box)



### HUI® Harris Uterine Injector CA\*

For HSG contrast media injection or dye introduction during tuboplasty surgery.

NUMBER	DESCRIPTION
6002	HUI Harris Uterine Injector (12/Box)



### ZUMI™ Zinnanti Uterine Manipulator Injector

Ergonomically curved shaft allows manipulation of both anteverted and retroverted uteri.

NUMBER	DESCRIPTION
ZSI1151	ZUMI™ 4.5 mm Uterine Manipulator (12/Box)
ZSI1152	ZUI™ 4.0 ZSI Uterine Injector (12/Box)
ZSI1153	ZUI™ 2.0 ZSI Uterine Injector (12/Box)




# INSORB®

## Subcuticular Skin Stapler

The INSORB Absorbable Skin Stapler is a subcuticular, interrupted closure modality designed to create uniform, symmetric dermis-to-dermis closures.

### INSORB® Skin Stapler

NUMBER	DESCRIPTION
2030	INSORB Absorbable Skin Stapler – Box of 6
9100	Adson Forceps



# Carter-Thomason CloseSure System®

Port Site Closure

## Carter-Thomason CloseSure System® CA★ C€


NUMBER	DESCRIPTION
CTI-512N	Carter-Thomason CloseSure System – Suture Passer, Pilot Guides: (1) 10/12 mm, (1) 5 mm, sterile/single-use (5 Guides/Box)
CTI-1015N	Carter-Thomason CloseSure Suture Passer, Pilot Guides: (1) 10/12 mm, (1) 15 mm, sterile/single-use (5 Guides/Box)
CTXL	Carter-Thomason CloseSure Suture Passer XL, Pilot Guides: (1) 10/12 mm XL, (1) 15 mm XL, sterile/single-use (3 Systems/Box)
PG-15S	Carter-Thomason CloseSure 15 mm Pilot Guides (Standard): sterile/single-use (5 Guides/Box)
PG-5/10	Pilot Guides: (1) 5 mm, (1) 10/12 mm, sterile/single-use (5 Guides/Box)
CTI-SP	(1) Suture Passer



Carter-Thomason CloseSure System:  
Ideal for Average to Obese Patients



Carter-Thomason CloseSure System XL:  
Ideal for Morbidly Obese and Bariatric Patients




The Carter-Thomason Suture Passer features a recessed, hinged jaw, which is engineered to enable the tip to capture the suture

# Carter-Thomason® II Port Closure System


Port Site Closure

## Carter-Thomason II Port Closure System C€

NUMBER	DESCRIPTION
CTI-1012P	(1) Suture Passer, (1) 10/12 mm Suture Guide (5 Systems/Box)
CTI-1015P	(1) Suture Passer, (1) 10/12 mm Suture Guide, (1) 15 mm Suture Guide (5 Systems/Box)
CTXL-P	(1) XL Suture Passer, (1) 10/12 mm XL Suture Guide, (1) 15 mm XL Suture Guide (3 Systems/Box)
CTI-15P	(1) Suture Passer, (1) 15 mm Suture Guide, (1) 15 mm XL Suture Guide (5 Systems/Box)



Extendable and blunt grasping fingers



Retractable grasping fingers

Suture Passer

15 mm Suture Guide

10/12 mm Suture Guide



# SureGrasp™ Suture Passer


SureGrasp Suture Passer with retractable grasping fingers for laparoscopic procedures.

## Carter-Thomason SG with SureGrasp™ Technology

NUMBER	DESCRIPTION	CE
CTI-1015SG	(1) Carter-Thomason CloseSure SureGrasp Suture Passer, Pilot Guides: (1) 10/12 mm, (1) 15 mm Pilot Guide sterile/single-use (5 Systems/Box)	CE
CTXL-SG	(1) Carter-Thomason CloseSure SureGrasp Suture Passer XL Pilot Guides: (1) 10/12 mm XL, (1) 15 mm XL sterile/single-use (3 Systems/Box)	CE
CTSG	(1) Carter-Thomason II Port Closure System Suture Passer (5 Passers/Box)	CE



# SeeClear® Smoke Evacuation System



SeeClear® is a passive, disposable, multistage filter system that traps smoke and aerosolized particles.

## SeeClear® Smoke Evacuation System €€

NUMBER	DESCRIPTION
SC082500	SeeClear Maximum Flow Rate: 8.0 L/min. at 15 mm HG pressure differential (10/Box)



# Trocars

## Apple-Hunt Secondary Cannula® / Trocar


The trocar features a smooth cannula/trocar transition and is offered in both a pyramidal and conical tip.

### A. 5 mm Apple-Hunt Secondary Cannula / Trocar CA☼ CE

NUMBER	DESCRIPTION
900-800	5 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar (5/Box)
900-801	5 mm Apple-Hunt Secondary Cannula / Conical Tip Trocar (5/Box)
900-802	5 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar with Side Port (5/Box)
900-803	5 mm Apple-Hunt Secondary Cannula / Conical Tip Trocar with Side Port (5/Box)
900-820	5 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar (2-Pack) (5/Box)
900-821	5 mm Apple-Hunt Secondary Cannula / Conical Tip Trocar (2-Pack) (5/Box)
900-822	5 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar with Side Port (2-Pack) (5/Box)
900-823	5 mm Apple-Hunt Secondary Cannula / Conical Tip Trocar with Side Port (2-Pack) (5/Box)
900-840	5 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar – extended length (5/Box)
900-841	5 mm Apple-Hunt Secondary Cannula / Conical Tip Trocar – extended length (5/Box)
900-842	5 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar – extended length (2-Pack) (5/Box)
900-844	5 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar with Side Port – extended length (5/Box)
900-845	5 mm Apple-Hunt Secondary Cannula / Conical Tip Trocar with Side Port – extended length (5/Box)
900-846	5 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar with Side Port – extended length (2-Pack) (5/Box)
900-847	5 mm Apple-Hunt Secondary Cannula / Conical Tip Trocar with Side Port – extended length (2-Pack) (5/Box)

### B. 12 mm Apple-Hunt Secondary Cannula / Trocar CA☼ CE

NUMBER	DESCRIPTION
900-850	12 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar with 5 mm and 10 mm Reducers (5/Box)
900-851	12 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar with 5 mm and 10 mm Reducers – extended length (5/Box)



### C. 10 mm Apple-Hunt Secondary Cannula / Trocar CA☼ CE

NUMBER	DESCRIPTION
900-860	10 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar with 5 mm Reducer (5/Box)
900-861	10 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar with 5 mm Reducer – extended length (5/Box)
900-863	10 mm Apple-Hunt Secondary Blunt Tip Trocar / Cannula with 5 mm Reducer (5/Box)

### D. 7 mm Apple-Hunt Secondary Cannula / Trocar CA☼ CE

NUMBER	DESCRIPTION
900-870	7 mm Apple-Hunt Secondary Cannula / Pyramidal Trocar with 5 mm Reducer (5/Box)
900-871	7/8 mm Apple-Hunt Secondary Cannula / Pyramidal Tip Trocar with 5 mm Reducer – extended length (5/Box)



### E. Apple-Hunt Instrument Organizer

The Apple-Hunt Organizer has a "saddle bag" design that can be draped across the patient's leg to keep instruments within easy reach. Includes eight individual compartments.

NUMBER	DESCRIPTION
900-305	Apple-Hunt Instrument Organizer (50/Box)

## Insufflation Needles

The Pneumo-Matic Insufflation Needle connects to the OR air source.

### F. Pneumo-Matic Insufflation Needle

NUMBER	DESCRIPTION
900-200	Pneumo-Matic Needle, 120 mm (10/Box)
900-210	Pneumo-Matic Needle, 150 mm (10/Box)

# Lone Star® Retractor System

## Lone Star® Disposable Retractor Rings

Rings are available in a variety of flexible, articulating shapes appropriate to the surgical site. Disposable plastic models are pre-sterilized and packaged singly in trays. Slotted design keeps elastic stay hooks precisely where you need them.


### Lone Star® Disposable Retractor Rings CA\* CE

NUMBER	DESCRIPTION
3301G	Retractor, sterile plastic, 25.3 cm X 11.3 cm
3302G	Retractor, sterile plastic, 22.0 cm X 11.3 cm
3304G*	Retractor, sterile plastic, 32.5 cm X 18.3 cm Includes 2 each Catheter Clips
3307G	Retractor, sterile plastic, 14.1 cm x 14.1 cm
3308G*	Retractor, sterile plastic, 28.3 cm X 18.3 cm Includes 2 each Catheter Clips
3309G	Retractor, sterile plastic, 18.6 cm X 8.9 cm
3310G	Retractor, sterile plastic, 25.0 cm x 25.0 cm
3715*	APS Retractor, sterile plastic, VELCRO® backing, 16.6 cm x 16.6 cm. Includes 2 each Catheter Clips

\*Packaged singly with 2 Catheter Clips

### Lone Star Retractor Kits CA\* CE

NUMBER	DESCRIPTION
3704	Gynecology Kit: each kit contains one disposable retractor, 32.5 cm x 18.3 cm, Model #3304G, and one pack of eight 5 mm sharp stay hooks, Model #3311-8G (5/Box)
3707	Colorectal Kit: each kit contains one square disposable retractor ring, 14.1 cm x 14.1 cm, Model #3307G, and one pack of eight 5 mm sharp stay hooks, Model #3311-8G (5/Box)



## Lone Star® Elastic Stay Hooks

Elastic stays adapt to the surgical site anatomy with or without the retractor ring.

### Sharp Hooks CA\* CE

NUMBER	DESCRIPTION
3311-1G	5 mm sharp hook (50/Case)
3311-8G	5 mm sharp hook (8/Pack)
3314-1G	3 mm sharp hook (50/Case)
3314-8G	3 mm sharp hook (8/Pack)
3325-4G	5 mm two-finger (4/Pack)
3327-4G	7 mm two-finger (4/Pack)
3329-4G	9 mm two-finger (4/Pack)

### Rake and Blade Hooks CA\* CE

NUMBER	DESCRIPTION
3330-4G	13 mm x 16 mm solid blade (4/Pack)
3331-4G	13 mm x 20 mm Spira blunt hook (4/Pack)
3332-4G	14 mm x 16 mm two-finger (4/Pack)
3333-4G	16 mm x 20 mm Spira two-finger (4/Pack)
3334-4G	20 mm x 16 mm four-finger (4/Pack)
3335-4G	20 mm x 20 mm Spira three-finger (4/Pack)
3338-4G	6.5 mm x 16 mm solid blade (4/Pack)

### Blunt Hooks CA\* CE

NUMBER	DESCRIPTION
3350-1G	12 mm blunt hook (50/Case)
3350-8G	12 mm blunt hook (8/Pack)
3350L-4G	12 mm blunt hook, elongated lead (4/Pack)
3384-4	12 mm two-finger blunt hook (4/Pack)
3550-1G	12 mm blunt hook, dual lead (50/Case)
3316-1G	5 mm blunt hook (50/Case)
3316-8G	5 mm blunt hook (8/Pack)
3205-8G	5 mm semi-blunt PRO-Stay (8/Pack)



# OR/Mobius®

Self-Retaining Elastic  
Abdominal Retractor

## OR/Mobius® CA★ CE

NUMBER	DESCRIPTION
900-508	Mini Retractor, 8 cm I.D. ring, for incisions to 6 cm (5/Box)
900-500	Small Retractor, 11 cm I.D. ring, for incisions to 9 cm (5/Box)
900-515	Medium Retractor, 15 cm I.D. ring, for incisions to 14 cm (5/Box)

Also available in larger sizes for maximizing exposure during  
C-sections or procedures with larger incision lengths.

## OB/Mobius® CA★ CE

NUMBER	DESCRIPTION
900-517	Medium-Large Retractor, 17 cm I.D. ring, for incisions to 15 cm (5/Box)
900-520	Large Retractor, 22 cm I.D. ring, for incisions to 17 cm (5/Box)



# Zeppelin® Clamps

Zeppelin Hysterectomy Clamps are offered in straight, curved and angled shapes for abdominal procedures.

## A. Zeppelin® Clamp, Gold Plated, Straight

NUMBER	LENGTH
16070GY	21 cm (8.5")
16075GY	24 cm (9.5")
16375GY	30 cm (12")

## B. Zeppelin Clamp, Gold Plated, Slightly Curved


NUMBER	LENGTH
16071GY	21 cm (8.5")
16076GY	24 cm (9.5")
16376GY	30 cm (12")

## C. Zeppelin Clamp, Gold Plated, Strong Curve

NUMBER	LENGTH
16072GY	21 cm (8.5")
16077GY	24 cm (9.5")
16377GY	30 cm (12")

## D. Zeppelin Clamp, Gold Plated, Angled

NUMBER	LENGTH
16073GY	21 cm (8.5")
16078GY	24 cm (9.5")
16378GY	30 cm (12")



## Vaginal Hysterectomy Clamps

### A. Zeppelin® "S" Vaginal Hysterectomy Clamp

NUMBER	LENGTH
16577GY	24 cm (9.5")

### B. Sheth Adnexa Zeppelin® Hysterectomy Clamp

NUMBER	LENGTH
16372GY	26 cm (10")

### C. Z-Clamp™

The clamps are offered in both straight and offset handles, including extended-length instruments for deeper pelvic procedures and obese patients.

#### Straight Clamp

Straight Handle		Offset Handle	
NUMBER	LENGTH	NUMBER	LENGTH
Z16070GY	21 cm (8.25")	Z16175GY	24 cm (9.5")
Z16075GY	24 cm (9.5")	Z16175GYL	30.5 cm (12")
Z16075GYLX	35.5 cm (14")	Z16175GYLX	35.5 cm (14")

#### Strongly Curved Clamp

Straight Handle		Offset Handle	
NUMBER	LENGTH	NUMBER	LENGTH
Z16072GY	21 cm (8.25")	Z16167GY	21 cm (8.25")
Z16077GY	24 cm (9.5")	Z16177GY	24 cm (9.5")
Z16077GYL	30.5 cm (12")	Z16177GYL	30.5 cm (12")
Z16077GYLX	35.5 cm (14")	Z16177GYLX	35.5 cm (14")

#### Slightly Curved Clamp

Straight Handle		Offset Handle	
NUMBER	LENGTH	NUMBER	LENGTH
Z16071GY	21 cm (8.25")	Z16166GY	21 cm (8.25")
Z16076GY	24 cm (9.5")	Z16176GY	24 cm (9.5")
Z16076GYL	30.5 cm (12")	Z16176GYL	30.5 cm (12")
Z16076GYLX	35.5 cm (14")	Z16176GYLX	35.5 cm (14")

#### Angled Clamp

Straight Handle		Offset Handle	
NUMBER	LENGTH	NUMBER	LENGTH
Z16073GY	21 cm (8.25")	Z16168GY	21 cm (8.25")
Z16078GY	24 cm (9.5")	Z16178GY	24 cm (9.5")
Z16078GYL	30.5 cm (12")	Z16178GYL	30.5 cm (12")
Z16078GYLX	35.5 cm (14")	Z16178GYLX	35.5 cm (14")

### SuperCut Scissors

#### D. SuperCut Scissors, Slight Curve

NUMBER	DESCRIPTION
SC16060	9", Slight Curve
SC16260	10.5", Slight Curve

#### E. SuperCut Scissors, Strong Curve

NUMBER	DESCRIPTION
SC16061	9", Strong Curve
SC16261	10.5", Strong Curve

#### F. SuperCut Scissors, Angled

NUMBER	DESCRIPTION
SC16062	9", Angled
SC16262	10.5", Angled

### G. Z-Scissors™

The curvature of the Z-Scissors is designed to work with the Z-Clamps for improved accuracy during dissection. Available in three curve configurations.

#### Z-Scissors, Slight Curve

NUMBER	LENGTH
Z16060GY	23 cm (9")
Z16060GYL	26.5 cm (10.5")
Z16060GYLX	30.5 cm (12")

#### Z-Scissors, Curve

NUMBER	LENGTH
Z16061GY	23 cm (9")
Z16061GYL	26.5 cm (10.5")
Z16061GYLX	30.5 cm (12")

#### Z-Scissors, Fully Curved

NUMBER	LENGTH
Z16062GY	23 cm (9")
Z16062GYL	26.5 cm (10.5")
Z16062GYLX	30.5 cm (12")





## Laparoscopic Instruments

### A. Nu-Tip® Instruments

Nu-Tip Instruments feature a reusable handle with disposable stainless steel tips.

NUMBER	DESCRIPTION
M2100M	Reusable Handle/Shaft Assembly, Insulated, Monopolar, 32 cm
M2101D	Scissors, Curved Metzenbaum, Disposable Tips (10/Box)
M2102D	Scissors, Straight, Disposable Tips (10/Box)
M2103D	Scissors, Hook, Disposable Tips (10/Box)

M2101D  
Curved

A.



M2102D  
Straight



M2103D  
Hook

### B. Knot Pushers

NUMBER	DESCRIPTION
M4500P	Clarke-Reich Micro: 6-0 and smaller, 29 cm, 5 mm
M4501P	Clarke-Reich Macro: 5-0 and larger, 29 cm, 5 mm
M4502P	Reddick-Saye Macro: 5-0 and larger, angled, 29 cm, 5 mm
M4503P	Closed Macro: 5-0 and larger, 29 cm, 5 mm
M4601P	Clarke-Reich Macro: 5-0 and larger, 45 cm, 5 mm
M4603P	Closed Macro: 5-0 and larger, 45 cm, 5 mm

B.



<b>Advincula Arch System</b> .....	12-13	<b>Hysterectomy Scissors</b>		<b>RUMI II / Koh-Efficient System</b> .....	11
<b>Advincula Arch System</b>		<b>SuperCut Scissors</b> .....	28	<b>RUMI II Tips</b> .....	10
Advincula Arch Koh-Efficient for Electrosurgery.....	13	<b>Z-Scissors</b> .....	28	<b>Sacrocolpopexy Tips</b> .....	14
Advincula Arch Koh-Efficient for Harmonic Scalpel / Laser .....	13	<b>INSORB Subcuticular Skin Stapler</b> .....	16	<b>Sacrocolpopexy Tips</b> .....	14
Advincula Delineator.....	9	<b>Insufflation Needles</b>		<b>ALLY Uterine Positioning System</b> .....	8
Hoyte Sacro Tips.....	14	<b>Pneumo-Matic Insufflation Needle</b> .....	22	<b>SeeClear Smoke Evacuation System</b>	
Sacrocolpopexy Tips .....	14	<b>Kronner Manipulator</b> .....	15	Maximum Flow Rate 8.0 Liters/min. ....	20
Sacrocolpopexy Tips.....	14	<b>Laparoscopic Instruments</b>		<b>Smoke Evacuation (Laparoscopic)</b> .....	20
<b>ALLY Uterine Positioning System</b> .....	8	<b>Knot Pushers</b> .....	29	<b>SureGrasp Suture Passer</b>	
<b>Apple-Hunt Instrument Organizer</b> .....	22	<b>Nu-Tip Instruments</b> .....	29	10/12 mm, 15 mm System .....	19
<b>Apple-Hunt Trocars</b>		<b>Lone Star Retractor System</b>		10/12 mm, 15 mm XL System .....	19
5 mm Cannula / Trocar .....	21	<b>Retractor Kits</b> .....	23	<b>Trocars</b> .....	21-22
7 mm Cannula / Trocar .....	22	<b>Retractor Rings</b> .....	23	<b>Uterine Manipulation</b>	
10 mm Cannula / Trocar .....	21	<b>Stay Hooks</b> .....	24	Advincula Arch.....	12
12 mm Cannula / Trocar .....	21	<b>Mobius – OR and OB</b> .....	25	Advincula Delineator.....	9
<b>Carter-Thomason CloseSure System</b>		<b>OB/Mobius Retractor System</b>		Cervical Sizer.....	11
5 mm, 10/12 mm System .....	17	<b>15 cm Incisions</b> .....	25	Disposable Diagnostic Uterine Manipulators.....	15
10/12 mm, 15 mm System .....	17	<b>17 cm Incisions</b> .....	25	Hoyte Sacro Tips.....	14
10/12 mm, 15 mm XL System .....	17	<b>OR/Mobius Retractor System</b>		Product Selection Chart .....	6-7
Pilot Guides.....	17	<b>6 cm Incisions</b> .....	25	RUMI II / Koh-Efficient System.....	11
<b>Carter-Thomason II Port Closure System</b>		<b>9 cm Incisions</b> .....	25	RUMI II System.....	10-11
10/12 mm System.....	18	<b>14 cm Incisions</b> .....	25	Sacrocolpopexy and Sacrocolpopexy Tips .....	14
10/12 mm, 15 mm System .....	18	<b>Port Site Closure (Laparoscopic)</b> .....	17-19	ALLY Uterine Positioning System.....	8
10/12 mm, 15 mm XL System .....	18	<b>Retractors</b>		<b>Z-Clamps</b> .....	27
<b>Disposable Diagnostic Uterine Manipulation</b>		<b>Lone Star Retractor System</b> .....	23-24	<b>Zeppelin Clamps</b> .....	26
Advincula Delineator .....	9	<b>OB/Mobius Retractor System</b> .....	25	Sheth Adnexa Zeppelin Hysterectomy Clamp .....	27
HUI Harris Uterine Injector.....	15	<b>OR/Mobius Retractor System</b> .....	25	Zeppelin "S" Vaginal Hysterectomy Clamp.....	27
HUMI Harris-Kronner Uterine Manipulator Injector ..	15	<b>RUMI II System</b>		<b>ZUMI Zinnanti Uterine Manipulator Injector</b> .....	15
Kronner Manipulator .....	15	Cervical Sizer.....	11		
ZUMI Zinnanti Uterine Manipulator Injector .....	15	Colpo-Pneumo Occluder .....	11		
<b>HUI Harris Uterine Injector</b> .....	15	Koh Cups, Ultem Resin for Electrosurgery .....	11		
<b>HUMI Harris-Kronner Uterine Manipulator Injector</b> .....	15	Koh Cups, Stainless Steel			
<b>Hysterectomy Clamps</b>		for Harmonic Scalpel / Laser .....	11		
Sheth Adnexa Zeppelin Hysterectomy Clamp .....	27	Hoyte Sacro Tips.....	14		
Z-Clamps .....	27	RUMI II Uterine Manipulator Handle.....	10		
Zeppelin Clamps.....	26	RUMI II / Koh-Efficient for Electrosurgery.....	11		
Zeppelin "S" Vaginal Hysterectomy Clamp.....	27	RUMI II / Koh-Efficient			
		for Harmonic Scalpel / Laser .....	11		



## Committed to Our Customers

At CooperSurgical, we strive for the highest level of expertise and knowledge of our products and to develop an ongoing relationship with each of our customers that is built on trust, quality and consistency. We look to provide service that exceeds not only our customers' expectations, but industry standards as well. CooperSurgical's Customer Service Department provides full-service order entry with extended hours to facilitate order placement. As an added convenience we also now offer online order entry. Call now to have your Customer Service Representative set up your account today. Once you are registered for online ordering, you will be able to place orders, track shipments, review prior purchase history and schedule orders for future delivery.

☞ Available for sale where CE mark is required

CA★ Available for sale in Canada

CE and CA labeling is subject to change. Contact Customer Service or your local sales representative for the most up-to-date information.

## Warranty Information

CooperSurgical, Inc. ("CooperSurgical") warrants that its products offered in this catalog, **other than** single-use products and consumables, will be free from defects in materials and workmanship for one year from the date of purchase. CooperSurgical warrants that its single-use products and consumables offered in this catalog will, through the earlier of 90 days from the date of purchase of each such product and the expiration date of that product stated on that product's package, be free from defects in materials and workmanship that are discovered prior to or during the single use of that product.

If CooperSurgical determines that a product offered in this catalog fails to conform to the CooperSurgical warranty applicable to that product, CooperSurgical will, as the sole remedy for that failure to conform, repair or replace that product, at CooperSurgical's discretion, free of charge.

To return a product to CooperSurgical, a customer must comply with CooperSurgical's Returned Goods Policy described in this catalog. A customer will not have any remedy if a CooperSurgical product does not conform to the warranty for that product unless that product is returned to CooperSurgical in accordance with that Returned Goods Policy. CooperSurgical will ship returned products that it repaired or replaced under warranty to the customer who returned those products, at CooperSurgical's expense F.O.B. the customer's facility. Under all other circumstances, CooperSurgical will ship returned products to the customer who returned those products at the customer's expense F.O.B. CooperSurgical's facility.

CooperSurgical's warranties, including those in this catalog, do not cover damage caused by misuse, improper care, improper use of chemicals or cleaning methods, loss, theft, servicing or sharpening by non-authorized personnel or negligent or intentional conduct on the part of the owner or user of the product. CooperSurgical's warranties do not cover normal wear and tear or general maintenance, including cleaning and sharpening services required as a result of normal use. Any modifications or changes to a product will void that product's warranty.

THE ABOVE WARRANTIES ARE IN LIEU OF, AND COOPERSURGICAL HEREBY DISCLAIMS, ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, WRITTEN OR ORAL, WITH RESPECT TO COOPERSURGICAL'S PRODUCTS, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO TERMS, CONDITIONS, UNDERSTANDINGS OR AGREEMENTS THAT PURPORT TO MODIFY THE ABOVE WARRANTIES OR THAT MAKE ANY ADDITIONAL WARRANTIES FOR ANY COOPERSURGICAL PRODUCT SHALL HAVE ANY LEGAL EFFECT UNLESS MADE IN WRITING AND SIGNED BY AN AUTHORIZED COOPERSURGICAL CORPORATE OFFICER.

COOPERSURGICAL SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE FOR LOST PROFITS, DAMAGES FROM LOSS OF USE OR LOST DATA, OR INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES UNDER ITS WARRANTIES OR OTHERWISE FOR ANY CLAIM RELATED TO COOPERSURGICAL'S PRODUCTS, EVEN IF COOPERSURGICAL HAS BEEN ADVISED, KNEW OR SHOULD HAVE KNOWN OF THE POSSIBILITY OF SUCH DAMAGES. COOPERSURGICAL'S LIABILITY WITH RESPECT TO A PRODUCT COVERED BY A WARRANTY OR OTHERWISE SHALL BE LIMITED IN ALL CIRCUMSTANCES TO THE PURCHASE PRICE OF THAT PRODUCT.

## 5 Ways to Order



### Local Sales Representative

[www.coopersurgical.com/pages/salesrepsearch.aspx](http://www.coopersurgical.com/pages/salesrepsearch.aspx)



### Email

U.S.: [Orders@coopersurgical.com](mailto:Orders@coopersurgical.com)

Int'l: [InternationalCustomerService@coopersurgical.com](mailto:InternationalCustomerService@coopersurgical.com)



### E-Ordering (U.S. Only)

Call 800.243.2974 to set up your account today.



### Phone Orders

800.243.2974 (8:00 AM - 8:30 PM EST)

Int'l: 203.601.9818



### Fax Orders

800.262.0105 (24 hours a day)

Int'l: 203.601.4747

## Returned Goods Policy

### 1. Goods will be accepted for return for the following reasons:

- If shipment was made without Customer's authorization or order
- If incorrect goods were shipped
- If defective goods are covered by the standard warranty

### 2. To return goods, you must contact a Customer Service Representative by telephone (1.800.243.2974 or 1.203.601.5200) for a Returned Merchandise Authorization (RMA) number. Items will not be accepted without an RMA number.

#### Please have the following information:

- Reason you wish to return the goods
- Quantity, description, part number, serial number of the goods
- Date of Receipt of Order
- Customer's Purchase Order and the CooperSurgical Invoice Number

**All used instruments must be cleaned and sterilized prior to being shipped for return.**

### 3. Shipment must be sent prepaid by Customer. Freight collect shipments will not be accepted, and goods will be returned to sender.

### 4. If Customer intends to return equipment ordered in error, the following restocking charges and terms will apply:

- There is no restocking fee if goods are returned within 30 days of shipment
- There is a restocking fee of 25% if goods are returned within 31-60 days of shipment
- Goods must be returned unused, in the original carton, and in marketable condition
- Refurbishing and replacement charges will be added to the restocking charges for damaged or missing items
- No return after 60 days
- No refund on sterile, single-use disposable products

**Send All Returns to: CooperSurgical  
95 Corporate Drive  
Trumbull, CT 06611**

## Trademark Notices

The marks bearing the symbols "TM" and "®" are trademarks and registered trademarks of CooperSurgical, Inc. Ultem® and Noryl® are registered trademarks of SABIC Innovation Plastics. Kronner Manipjector®, HUI® and HUMI® are registered trademarks of Richard Kronner. SILASTIC® is a registered trademark of Dow Corning Corporation. VELCRO® is a registered trademark of Velcro Industries B.V.

ALLY Uterine Positioning System®

Advincula Delineator™  
with Koh-Efficient® Technology

Koh-Efficient®

RUMI® II System

Advincula Arch®

Kronner Manipujector®

INSORB® Subcuticular Skin Stapler

Carter-Thomason CloseSure System®

SeeClear® Laparoscopic  
Smoke Evacuation System

Lone Star® Retractor System

OR/Mobius®



CooperSurgical®



800.243.2974

[www.coopersurgical.com](http://www.coopersurgical.com)

95 Corporate Drive, Trumbull, CT 06611