

# Fundus Sound™ XM

## Disposable Uterine Sound

### Key Features and Benefits

- Easy-to-read imprinted and notched centimeter markings for accurate uterine depth determination
- Small diameter for procedural comfort
- Sterile, single-use
- Malleable to adapt to anatomical variations
- Latex-free
- Class IV Medical-Grade Polymer



# Fundus Sound™ XM

Disposable Uterine Sound

## Ergonomic design for easy handling

Designed for effortless hand maneuverability, the Fundus Sound™ XM can simplify procedures.

## Easy to read for fast and precise measurement

Clear depth markings and centimeter notches present a quick, easy-to-read display.

## Balanced flexibility

The advanced polymer material is malleable and offers the correct amount of pliancy to follow the curve of the cervix and uterus for quick and efficient procedures.

## Small diameter

Measuring only 2.8 mm at the very tip and widening to 3.2 mm at the 10 cm mark, the Fundus Sound™ XM is comfortable for the patient.

## Simplify post-procedure cleanup

The Fundus Sound™ XM is designed for single patient use. Simply dispose per local approved Biomedical Waste procedures.

**Made from an advanced medical-grade polymer material, the single-use and sterile Fundus Sound™ XM offers a number of distinct advantages for your procedures.**

The Fundus Sound™ XM can be used for the following procedures:

- To dilate and check patency of the cervix
- To measure the depth of the uterine cavity
- For preparation of proper IUD insertion
- For routine uterine examinations when there is abnormal or dysfunctional uterine bleeding

## Ordering Information

Product number – 1196

Packaging – 25 sterile Fundus Sound™ XM per box



CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products, which aid clinicians in the management and treatment of commonly seen conditions.

800.243.2974 | 203.601.5200 | [coopersurgical.com](http://coopersurgical.com) | 95 Corporate Drive | Trumbull, CT 06611  
Also on YouTube.com/coopersurgical