

Milex® Wide-Seal Diaphragms

Latex-free brand of choice

Features and Benefits

- Hypoallergenic, silicone diaphragms with a unique wide seal
- Autoclavable fitting kit
- Eight sizes for true fit
- Easy-to-use Omniflex style or Arcing style
- Distortion-free and Cadmium-free tension coil



*Must be used with spermicidal gel.

 CooperSurgical

Milex® Wide-Seal Diaphragms*

Latex-free



Proven Performance For Over 40 Years

- Full line of latex-free, hypoallergenic silicone vaginal diaphragms
- Wide-Seal provides increased suction for added protection
- Durable, autoclavable fitting kits and diaphragms

"True Fit" Technology

- Available in both Arcing Style and OmniFlex Style, assuring the flexibility needed for both ease-of-use and comfortable fit
- Both styles available in 8 sizes (60-95 mm in 5 mm increments)
- Combined with the unique Wide-Seal, Milex Diaphragms offer your patients a true fit

*Must be used with spermicidal gel.

Milex® Wide-Seal Diaphragms Ordering Information

Code Product Description

Milex® Wide-Seal Arcing-Style Diaphragm

Tension adjusted spring curves in one place

MXWSFS	Arcing -Style Diaphragm Fitting Set	<i>*includes these 5 sizes</i>
MXWS60	Arcing-Style Diaphragm 60 mm	
MXWS65	Arcing-Style Diaphragm 65 mm*	
MXWS70	Arcing-Style Diaphragm 70 mm*	
MXWS75	Arcing-Style Diaphragm 75 mm*	
MXWS80	Arcing-Style Diaphragm 80 mm*	
MXWS85	Arcing-Style Diaphragm 85 mm*	
MXWS90	Arcing-Style Diaphragm 90 mm	
MXWS95	Arcing-Style Diaphragm 95 mm	

Milex® Wide-Seal OmniFlex-Style Diaphragm

Distortion-free spring provides arc no matter where the rim is compressed

MXWFFS	OmniFlex-Style Diaphragm Fitting Set	<i>*includes these 5 sizes</i>
MXWF60	OmniFlex-Style Diaphragm 60 mm	
MXWF65	OmniFlex-Style Diaphragm 65 mm*	
MXWF70	OmniFlex-Style Diaphragm 70 mm*	
MXWF75	OmniFlex-Style Diaphragm 75 mm*	
MXWF80	OmniFlex-Style Diaphragm 80 mm*	
MXWF85	OmniFlex-Style Diaphragm 85 mm*	
MXWF90	OmniFlex-Style Diaphragm 90 mm	
MXWF95	OmniFlex-Style Diaphragm 95 mm	

CooperSurgical

CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products that aid clinicians in the management and treatment of commonly seen conditions.

Find out how Milex® Wide-Seal Diaphragms can benefit your practice... and your patients. Call 800.243.2974 or visit www.coopersurgical.com.

800.243.2974 | 203.601.5200 | coopersurgical.com | 95 Corporate Drive | Trumbull, CT 06611
Also on [YouTube.com/coopersurgical](https://www.youtube.com/coopersurgical)