

Wallach®

In-Office Cryosurgical Systems Offer Fast Freeze and Precise Control With a Large Selection of Cryotips

Wallach LL100™ System (N₂O) Wallach LLCO₂™ System (CO₂)

Wallach's exclusive two-trigger cryosurgical system performs freeze and defrost in less time than passive defrost systems. The physician has ultimate control with the LL100's two-trigger design because it allows for freeze and defrost only when triggered. A wide selection of highly conductive, thermal and reusable cryotips can be changed at any time without venting.

- Lightweight, well-balanced and maneuverable in-office system
- Provides full view of the target tissue for cryoablation
- Choose from 23 autoclavable cryotips with disposable plastic shields

Single-Trigger Design is Optimized for CO₂ Gas

For those physicians who prefer a CO₂ system, Wallach offers the LLCO₂™ single-trigger cryosurgical system. The system maximizes temperature for ablation while minimizing the clogging and sputtering commonly associated with CO₂ gas.



Wallach®

Cryosurgical Tips for Use with LL100™ and LLCO₂™
Highly conductive, thermal, reusable, sterilizable cryotips with disposable plastic shields.



Item #	Description
900001	LL100 Cryosurgical System – N ₂ O (Freezer Only)
900160	LLCO ₂ Cryosurgical System – CO ₂ (Freezer Only)
901060	N ₂ O Cylinder (unfilled)
901061	CO ₂ Cylinder (unfilled)
901011	Cart for 20 lb. Cylinder

Item #	Description
900506-1-N2O	WA1000B Cryosurgical Console – Blank N ₂ O LL100
900506-1-CO2	WA1000B Cryosurgical Console – Blank CO ₂ LL100
900506-2-N2O	WA1000B Cryosurgical Console – Digital Temp. (N ₂ O) Indicator, Timer, Base, and LL100
900506-2-CO2	WA1000B Cryosurgical Console – Digital Temp. (CO ₂) Indicator, Timer, Base, and LLCO ₂



WA1000B™ Cryo Console System

Uses LL100™ two-trigger or LLCO₂™ single-trigger freezers.

- High-performance features and easy-to-use controls
- Temperature Indicator with integrated timer for precise documentation
- Cart carries a single 20 lb. cylinder
- Use with Wallach's large selection of cryotips and shields



800.243.2974 | 203.601.5200
www.coopersurgical.com/distribution
95 Corporate Drive, Trumbull, CT 06611