



INSTRUMENT DIMENSIONS

| Catalog No. | Size | A | | B | | Wire Length | |
|-------------|---------------------|----|------|----|------|-------------|-------|
| | | mm | in. | mm | in. | mm | in. |
| 900-150 | Small | 9 | .354 | 12 | .472 | 15 | .590 |
| 900-151 | Medium | 12 | .472 | 15 | .590 | 19 | .756 |
| 900-152 | Large | 15 | .590 | 18 | .708 | 23 | .921 |
| 900-154 | Medium Extended | 12 | .472 | 20 | .787 | 23 | .918 |
| 900-155 | Large Shallow | 15 | .590 | 8 | .315 | 17 | .669 |
| 900-157 | Small Wide | 9 | .354 | 12 | .472 | 17.9 | .705 |
| 900-158 | Small Wide Extended | 9 | .354 | 20 | .787 | 25.8 | 1.020 |




BIOPSY SIZES


900-150
Small




900-151
Medium



900-152
Large



900-154
Medium Extended



900-155
Large Shallow

