

Lone Star[®]



RETRACTOR SYSTEM

Exceptional
Visibility
and Access


Self-retaining retractor system
offers unobstructed visualization,
and access for a wide range
of surgical procedures.




3304G – Sterile plastic, 32.5 cm x 18.3 cm



3307G – Sterile plastic, 14.1 cm x 14.1 cm




3308G – Sterile plastic, 28.3 cm x 18.3 cm




3309G – Sterile plastic, 18.6 cm x 8.9 cm



3310G – Sterile plastic, 25.0 cm x 25.0 cm



3715 – Sterile plastic, 16.6 cm x 16.6 cm



3301G – Sterile plastic, 25.3 cm x 11.3 cm



3302G – Sterile plastic, 22.0 cm x 11.3 cm



SHARP HOOKS



3311-8G*
5 mm sharp hook
(8/pack)



3314-8G*
3 mm sharp hook
(8/pack)



3325-4G
5 mm two-finger
(4/pack)



3327-4G
7 mm two-finger
(4/pack)



3329-4G
9 mm two-finger
(4/pack)


RAKE AND BLADE HOOKS



3330-4G
13 mm x 16 mm solid
blade (4/pack)



3331-4G
13 mm x 20 mm Spira
blunt hook (4/pack)




3332-4G
14 mm x 16 mm
two-finger (4/pack)



3333-4G
16 mm x 20 mm Spira
two-finger (4/pack)



3334-4G
20 mm x 16 mm
four-finger (4/pack)



3335-4G
20 mm x 20 mm Spira
three-finger (4/pack)



3338-4G
6.5 mm x 16 mm
solid blade (4/pack)

BLUNT HOOKS



3350-8G*
12 mm blunt hook
(8/pack)



3350L-4G
12 mm blunt hook,
elongated lead (4/pack)



3384-4
12 mm two-finger
blunt hook (4/pack)



3550-1G
12 mm blunt hook,
dual lead (50/case)



3316-8G*
5 mm blunt hook
(8/pack)



3205-8G
5 mm semi-blunt
PRO-Stay (8/pack)



Lone Star®

RETRACTOR SYSTEM

Self-Retaining Retractor System That Provides Exceptional Visibility and Access

Disposable Retractor Ring Specifications

NUMBER	DESCRIPTION
3301G	Sterile plastic, 25.3 cm x 11.3 cm
3302G	Sterile plastic, 22.0 cm x 11.3 cm
3304G*	Sterile plastic, 32.5 cm x 18.3 cm
3307G	Sterile plastic, 14.1 cm x 14.1 cm
3308G*	Sterile plastic, 28.3 cm x 18.3 cm
3309G	Sterile plastic, 18.6 cm x 8.9 cm
3310G	Sterile plastic, 25.0 cm x 25.0 cm
3715*	Sterile plastic, 16.6 cm x 16.6 cm

Made of Noryl® resin, the disposable retractor rings are pre-sterilized and packaged singly in double-peel pouches for easy delivery to the sterile surgical field.

*Packaged singly with 2 catheter clips

Lone Star Retractor Kits

NUMBER	DESCRIPTION
3704	Gynecology Kit - One large snowman disposable retractor ring, 32.5 cm x 18.3 cm – Model 3304G - One eight-pack of 5 mm sharp stay hooks – Model #3311-8G (5 kits/box)
3707	Colorectal Kit - One square disposable retractor ring, 14.1 cm x 14.1 cm – Model #3307G - One eight-pack of 5 mm sharp stay hooks – Model #3311-8G (5 kits/box)

US Patent Nos. 5,785,649; 5,964,697

Lone Star® Retractor System is a registered trademark of CooperSurgical, Inc.
Noryl® is a registered trademark of SABIC Innovation Plastics.

Stay Hook Specifications

NUMBER	DESCRIPTION
3205-8G	5 mm semi-blunt PRO-Stay (8/pack)
3311-1G	5 mm sharp hook (50/case)
3311-8G	5 mm sharp hook (8/pack)
3314-1G	3 mm sharp hook (50/case)
3314-8G	3 mm sharp hook (8/pack)
3316-1G	5 mm blunt hook (50/case)
3316-8G	5 mm blunt hook (8/pack)
3325-4G	5 mm two-finger (4/pack)
3327-4G	7 mm two-finger (4/pack)
3329-4G	9 mm two-finger (4/pack)
3330-4G	13 mm x 16 mm solid blade (4/pack)
3331-4G	13 mm x 20 mm Spira blunt hook (4/pack)
3332-4G	14 mm x 16 mm two-finger (4/pack)
3333-4G	16 mm x 20 mm Spira two-finger (4/pack)
3334-4G	20 mm x 16 mm four-finger (4/pack)
3335-4G	20 mm x 20 mm Spira three-finger (4/pack)
3338-4G	6.5 mm x 16 mm solid blade (4/pack)
3350-1G	12 mm blunt hook (50/case)
3350-8G	12 mm blunt hook (8/pack)
3350L-4G	12 mm blunt hook, elongated lead (4/pack)
3384-4	12 mm two-finger blunt hook (4/pack)
3550-1G	12 mm blunt hook, dual lead (50/case)

50/case stays are individually wrapped. All stay hooks are single-use.

800.243.2974 | 203.601.5200
www.coopersurgical.com
95 Corporate Drive, Trumbull, CT 06611

