

# SureGrasp™


## SureGrasp™ Suture Passer Offers Precise Control and Safety for Port Site Closure

Our advanced SureGrasp Suture Passer features long, retractable grasping fingers that easily and safely capture sutures during laparoscopic and other surgical procedures. When the grasping fingers are closed, the SureGrasp tip seals to quickly and safely pass suture. Both of our Carter-Thomason Systems feature the new SureGrasp technology, which helps ensure:

- Complete closure of fascia and peritoneum
- Pneumoperitoneum is maintained
- System works under direct visualization
- Reproducibility through easy-to-learn techniques
- Precise closure for a range of surgical procedures, including hernia repair




### SureGrasp Fingers Opened



Extendable and blunt grasping fingers facilitate suture capture

### SureGrasp Fingers Closed




Retractable grasping fingers seal the needle tip to prevent tissue coring

For more information about the Carter-Thomason SG, call 800.243.2974 or 203.601.5200 or visit [coopersurgical.com](http://coopersurgical.com)

**CooperSurgical**

# Choose the System That's Right for Your Surgical Procedure

## Carter-Thomason SG



### Pilot Guide

Proven in numerous laparoscopic procedures, the Carter-Thomason® SG Pilot Guide closes port sites efficiently and reliably.

## Carter-Thomason II Port Closure System



### Improved Suture Guide

The advanced Carter-Thomason® II Suture Guide provides precise, consistent closure on any body type.

## Ordering Information

Carter-Thomason SG with SureGrasp Technology	
Product / Order #	Contents
Carter-Thomason SG CloseSure System CTI-1015SG	<ul style="list-style-type: none"> <li>• Carter-Thomason CloseSure SureGrasp Suture Passer</li> <li>• Pilot Guides: (1) 10/12 mm, (1) 15 mm</li> <li>• Sterile/Single-use</li> <li>• Qty: 5 systems per box</li> </ul>
Carter-Thomason SG CloseSure System XL® CTXL-SG	<ul style="list-style-type: none"> <li>• Carter-Thomason CloseSure SureGrasp Suture Passer XL</li> <li>• Pilot Guides: (1) 10/12 mm XL, (1) 15 mm XL</li> <li>• Sterile/Single-use</li> <li>• Qty: 3 systems per box</li> </ul>

## Ordering Information

Carter-Thomason II Port Closure System	
Product / Order #	Contents
Carter-Thomason II Port Closure System CTI-1012P	<ul style="list-style-type: none"> <li>• (1) Carter-Thomason II Port Closure System Suture Passer</li> <li>• (1) 10/12 mm Suture Guide</li> <li>• Qty: 5 systems per box</li> </ul>
Carter-Thomason II Port Closure System CTI-15P	<ul style="list-style-type: none"> <li>• (1) Carter-Thomason II Port Closure System Suture Passer</li> <li>• (1) 15 mm Suture Guide</li> <li>• Qty: 5 systems per box</li> </ul>
Carter-Thomason II Port Closure System CTI-1015P	<ul style="list-style-type: none"> <li>• (1) Carter-Thomason II Port Closure System Suture Passer</li> <li>• (1) 10/12 mm Suture Guide</li> <li>• (1) 15 mm Suture Guide</li> <li>• Qty: 5 systems per box</li> </ul>
Carter-Thomason II Port Closure System XL CTXL-P	<ul style="list-style-type: none"> <li>• (1) Carter-Thomason II Port Closure System XL Suture Passer</li> <li>• (1) 10/12 mm XL Suture Guide</li> <li>• (1) 15 mm XL Suture Guide</li> <li>• Qty: 3 systems per box</li> </ul>

Find out how the Carter-Thomason® family of port closure systems can benefit your practice... and your patients. Call 800.243.2974 or 203.601.5200 or visit [www.coopersurgical.com](http://www.coopersurgical.com). Also, see us on YouTube at [YouTube.com/CooperSurgical](http://YouTube.com/CooperSurgical).

**CooperSurgical**